

Chemical Compatibility of ChemoClave™ and ChemoLock™

Thank you for your interest in the ICU Medical ChemoClave and ChemoLock closed system transfer devices (CSTDs). ICU has been innovating CSTD products since 2007 and has extensive experience with medical-grade consumable devices that are exposed to a wide variety of drug formulations.

ICU Medical manufactures all ChemoClave and ChemoLock products exclusively without the use of ortho-phthalates, a class of organic plasticizers frequently used in PVC. These plasticizers can be added in concentrations up to 40% to render the rigid PVC more pliable and suitable for extrusion into components used for IV infusions.^{1,2} These ortho-phthalates, particularly the most commonly used form di(2-ethylhexyl) phthalate (DEHP), are linked to developmental and reproductive toxicity, endocrine disruption, and are listed as probable carcinogens.^{3,4} These plasticizers do not chemically bind to the PVC backbone and can therefore migrate (leach) out when in contact with certain drug formulations or under certain environmental conditions and travel down into the patient's bloodstream.⁵ For instance, high amounts of DEHP were extracted from PVC bags containing paclitaxel, cyclosporine, teniposide, 5-fluorouracil, and quinine.⁶ With DEHP being fat soluble, the leaching is particularly enhanced from PVC bags and tubing in the presence of surfactants like Kolliphor EL and Polysorbate, which are commonly used in chemotherapeutic formulations.^{7,8} The absence of ortho-phthalates like DEHP in ICU Medical CSTD products protects patients from potential exposure to these chemical plasticizers.

Microbial ingress and sterility integrity have been verified independently for a 7-day period, where under a simulated use protocol, the ChemoClave system was able to maintain complete sterility of the drug.⁹⁻¹¹

Many of the antineoplastic drugs, such as etoposide and paclitaxel, contain chemical solvents, surfactants, or fat emulsions, which have the potential to damage certain medical plastics. ICU Medical developed the ChemoClave and ChemoLock using medical-grade plastics that are compatible with such drug formulations—since the inception of the product line in 2007. ICU Medical has extensively tested ChemoClave and ChemoLock products with a number of representative drug formulations that demonstrate compatibility for several therapeutic classes.¹²⁻¹⁶ These test results are applicable to the following products within the ChemoClave and ChemoLock product lines that use identical materials of construction:

- › Spiros™ closed male luer
- › Clave™ and MicroClave™ and Neutron™ needlefree connectors
- › ChemoClave bag spikes
- › ChemoClave vial spikes
- › ChemoClave administration sets
- › ChemoLock injector
- › ChemoLock vial and bag spikes
- › ChemoLock administration sets
- › Diana™ cassettes and diluent sets

Drug formulations were selected for chemical compatibility testing based on their widespread use, chemical composition, and known potential to interact with or damage plastics. These studies were designed to include both single-dose and multi-dose drugs that allow for extended storage periods, such as Cisplatin (up to 28 days), thus offering the longest potential opportunity for interaction with the devices.

The ChemoClave system, as well as the ChemoLock system that is constructed of materials identical to ChemoClave, were exposed to over 13 antineoplastic drugs in their undiluted form.¹²⁻¹⁶ Administration sets were similarly exposed to drug formulations at three times their therapeutic dose in order to simulate the worst-case exposure for administration equipment. Three separate tests were completed including a functional integrity (leak) test, drug stability test, and plastic migration test. The results demonstrate that the ChemoClave and ChemoLock products are compatible with common therapeutic classes for a minimum of 24-hours and a maximum of 28-days for multi-use drugs.

The following table includes a list of drugs that are known to fall into these common therapeutic classes and may therefore be considered compatible with the ChemoClave and ChemoLock systems.

5-Fluorouracil	Dacarbazine	Liposomal Daunorubicin &	Porfimer sodium
Aflibercept	Dactinomycin	Cytarabine	Pralatrexate
Aldesleukin	Daunorubicin	Liposomal Doxorubicin	Raltitrexed
Alemtuzumab	Decitabine	Liposomal Vincristine	Rituximab
Amsacrine	Deferoxamine	Melphalan	Romidepsin
Arsenic trioxide	Denosumab	Mesna	Streptozocin
Atezolizumab	Desferal	Mechlorethamine HCl	Tacrolimus
Azacitidine	Dexrazoxane	Methotrexate	Temozolomide
BCG (Bacillus Calmette-Guérin)	Docetaxel	Mitomycin	Temsirolimus
Belimumab	Doxorubicin	Mitoxantrone	Teniposide
Bendamustine HCl	Epirubicin	Mustine	Testosterone cypionate (DEPO-TESTOSTERONE)
Bevacizumab (Avastin)	Eribulin	Mycophenolate mofetil	Thiotepa
Bleomycin	Etoposide	Natalizumab	Topotecan
Blinatumomab	Fludarabine	Nelarabine	Tocilizumab
Brentuximab vedotin	Gemcitabine	Ofatumumab	Trastuzumab
Busulfan	Gemtuzumab ozogamicin	Omalizumab	Trastuzumab emtansine
Cabazitaxel	Hydrocortisone	Oxaliplatin	Treosulfan
Carboplatin	Idarubicin	Paclitaxel	Triptorelin
Carfilzomib	Ifosfamide	Paclitaxel (Protein-Bound)	Ustekinumab
Carmustine	Infliximab	Panitumumab	Vedolizumab
Cetuximab	Interferon alfa-2b	Pegaspargase	Vinblastine
Cidofovir	Intron A	Peginterferon alfa-2a	Vincristine
Cisplatin	Ipilimumab	Pembrolizumab	Vindesine
Cladribine	Irinotecan	Pemetrexed	Vinflunine
Clofarabine	Ixabepilone	Pentostatin	Vinorelbine
Cyclophosphamide	L-Asparaginase	Pertuzumab	
Cyclosporine	Leuprolide acetate	Plerixafor	
Cytarabine	Lipiodol	Prexasertib (LY2606368)	

Busulfan (Busulfex®) and Amsacrine (AMSA PD®) Specifics: All ChemoClave and ChemoLock products are compatible with the diluted solution of Busulfan at a 0.536 mg/mL concentration, which corresponds to a 0.7 mg/mL concentration in 500 mL of 0.9% sodium chloride, as instructed in the prescribing information.^{17,18} Similarly, all ChemoClave and ChemoLock products are compatible with the diluted solution of Amsacrine up to a concentration of about 0.146 mg/mL in 5% dextrose, as instructed in the prescribing information.¹⁹ These concentrations represent the typical administration scenarios used for delivery of Busulfan and Amsacrine to patients. Test results show that ChemoClave and ChemoLock administration sets meet functional requirements following exposure to 0.536 mg/mL Busulfan and 0.146 mg/mL Amsacrine solutions over a 24-hour period and are therefore deemed compatible for administration within that time period to the patient when in the diluted form.

The ChemoClave 13 mm vial spike CH-82, CH-62, and vented vial spike CH-72 are compatible with the undiluted solutions of Busulfan and Amsacrine. No other ChemoClave products are compatible with undiluted Busulfan and Amsacrine solutions. This means that before infusing Busulfan or Amsacrine using a ChemoClave bag spike, the drugs must first be diluted following the recommended methods:

Busulfan: After drawing the drug out of a vial through the CH-82, CH-62, or CH-72 vial spike, connect the drug syringe to the Clave connector on the bag spike, aspirate back a minimum of 6.5 mL of solution from the IV bag into the drug syringe, and re-inject into the IV bag through the Clave bag spike.

Amsacrine: Amsacrine STOCK SOLUTION should be prepared as per prescribing information.¹⁹ After drawing the STOCK SOLUTION out of the vial through the CH-82, CH-62, or CH-72 vial spike, connect the drug syringe to the Clave connector on the bag spike, aspirate back a minimum of 6.5 mL of solution into the drug syringe, and re-inject into the IV bag through the Clave bag spike.

The ChemoLock 13 mm vial spike CL-82, CL-62, and vented CL-72 and CL2000S syringe injector are compatible with undiluted Busulfan and Amsacrine for a maximum duration of 1 hour. No other ChemoLock products are compatible with undiluted Busulfan and Amsacrine. This means that before infusing Busulfan or Amsacrine using a ChemoLock bag spike, the drugs must first be diluted using the following recommended methods:

Busulfan: After drawing the drug out of a vial through a CL-82, CL-62, or CL-72 vial spike and CL2000S syringe injector, connect the CL2000S to the port connector on the bag spike, aspirate back a minimum of 6.5 mL of solution from the IV bag into the drug syringe, and re-inject into the IV bag through the port bag spike.

Amsacrine: The Amsacrine STOCK SOLUTION should be prepared as per prescribing information.¹⁹ After drawing the STOCK SOLUTION out of the vial through the CL-82, CL-62, or CL-72 vial spike and CL2000S syringe injector, connect the CL2000S to the port connector on the bag spike, aspirate back a minimum of 6.5 mL of solution into the drug syringe, and re-inject into the IV bag through the Clave bag spike.

Thank you again for your continued support of ICU Medical CSTDs. For further questions or concerns, please visit our website at www.icumed.com or contact your local sales representative.

1. Marie C, Hamlaoui S, Bernard L, Bourdeaux D, Sautou V, Lémyery D, Vendittelli F, Sauvart-Rochat M-P. Exposure of hospitalised pregnant women to plasticizers contained in medical devices. *BMC Women's Health*. 2017;17:45. DOI: <https://dx.doi.org/10.1186%2F12905-017-0398-7>.
2. Treleano A, Wolz G, Brandsch R, Welle F. Investigation into the sorption of nitroglycerin and diazepam into PVC tubes and alternative tube materials during application. *International journal of pharmaceutics*. 2009;369(1):30-7. DOI: <https://doi.org/10.1016/j.ijpharm.2008.10.024>.
3. Scientific Committee on Emerging and Newly-Identified Health Risks (SCENIHR). The safety of medical devices containing DEHP plasticized PVC or other plasticizers on neonates and other groups possibly at risk (2015 update) Opinion. Luxembourg: European Commission, 2016.
4. Center for Devices and Radiological Health. Safety Assessment of Di(2-ethylhexyl)phthalate (DEHP) Released from PVC Medical Devices. Rockville, MD: U.S. Food and Drug Administration.
5. Loff S, Kabs F, Subotic U, Schaible T, Reinecke F, Langbein M. Kinetics of diethylhexyl phthalate extraction from polyvinylchloride infusion lines. *Journal of Parenteral and Enteral Nutrition*. 2002;26(5):305-9. DOI: [10.1177/0148607102026005305](https://doi.org/10.1177/0148607102026005305).
6. Kambia NK, Dine T, Dupin-Spriet T, Gressier B, Luyckx M, Goudaliez F, Brunet C. Compatibility of nitroglycerin, diazepam and chlorpromazine with a new multilayer material for infusion containers. *Journal of Pharmaceutical and Biomedical Analysis*. 2005;37(2):259-64. DOI: <https://doi.org/10.1016/j.jpba.2004.10.020>.
7. Faouzi Me-A, Dine T, Luyckx M, Brunet C, Mallevais ML, Goudaliez F, Gressier B, Cazin M, Kablan J, Cazin JC. Stability, compatibility and plasticizer extraction of miconazole injection added to infusion solutions and stored in PVC containers. *Journal of Pharmaceutical and Biomedical Analysis*. 1995;13(11):1363-72. DOI: [https://doi.org/10.1016/S0731-7085\(99\)00204-6](https://doi.org/10.1016/S0731-7085(99)00204-6).
8. Bagel-Boithias S, Sautou-Miranda V, Bourdeaux D, Tramier V, Boyer A, Chopineau J. Leaching of diethylhexyl phthalate from multilayer tubing into etoposide infusion solutions. *American Journal of Health-System Pharmacy*. 2005;62(2):182-8.
9. ICU Medical. Evaluation of the Clave® technology and resistance to microbial ingress. 2012. Available from: http://www.icumed.com/media/73891/m1-1212_clave_microbial_ingress_white_paper_rev03.pdf.
10. ICU Medical. Microbial Ingress Study for ChemoClave® Devices. 2014. Available from: http://www.icumed.com/media/485714/M1-1474-ChemoClave_Microbial_Ingress-Rev01-Web.pdf.
11. ICU Medical. Microbial Ingress Study for ChemoLock® Devices. 2016. Available from: http://www.icumed.com/media/517659/m1-1452-chemolock7daymicrobialingress-rev02_web.pdf.
12. ICU Medical. Antineoplastic Drug Compatibility with the ChemoClave™ System. 2012. Available from: <http://www.icumed.com/media/122025/m1-1185%20antineoplastic%20compatibility%20rev%2005-web.pdf>.
13. ICU Medical. Extended Chemical Compatibility of CLAVE Based Products and Spiros. 2007. ENG-33.
14. ICU Medical. Universal Vial Spike Chemotherapy Drugs Compatibility Tests. 2008. ENG-219.
15. ICU Medical. CH-10, CH2000, CH-77 Compatibility Tests with HERCEPTIN/TRASTUZUMAB. 2010. ENG-317.
16. ICU Medical. CH-12, CH2000S, CH-77 Compatibility Tests with DOXORUBICIN, MITOMYCIN, GEMCITABINE. 2010. ENG-318.
17. Accord Healthcare Ltd. Busulfan 6 mg/ml concentrate for solution for infusion. 2018. Available from: <https://www.medicines.org.uk/emc/product/2160/smpc#companyDetails>.
18. Actavis Pharma Inc. Full Prescribing Information - BUSULFAN- busulfan injection, solution, concentrate. 2018. Available from: <https://dailymed.nlm.nih.gov/dailymed/fda/fdaDrugXsl.cfm?setid=45650bfe-5f73-4591-812a-14d0feeecaad&type=display>.
19. Erf Canada 2012 Inc. Full Prescribing Information - AMSA PD™ Amsacrine Injection 75 mg/1.5 ml Ampoule (50 mg/ml). 2005. Available from: <http://eci2012.net/wp-content/uploads/2015/03/Amsa-PM-En-100449.02-19Feb2015.pdf>.