



Product Reference Sheet

NanoClave™

Low Profile Neutral Displacement
Needlefree Connector



Disinfect to Protect

- › When placing a new NanoClave on a catheter, disinfect the catheter hub and prime the NanoClave.
- › Before accessing a NanoClave, always disinfect the injection site with the approved antiseptic per facility protocol.
- › Scrub the injection site in accordance with facility protocol for appropriate scrubbing and dry times.



Administer or Aspirate

- › Attach IV tubing, syringe or blood tube holder to NanoClave by inserting the luer and twisting 1/4 turn, or until a friction fit is achieved.
- › Do not over-tighten a luer beyond the friction fit as this may damage both the luer and the NanoClave.
- › To disconnect, grasp NanoClave and then twist mating luer away from NanoClave until loose. Do not hold catheter hub during disconnect as this may cause accidental removal of NanoClave from hub.



Flush After Each Use

- › Flush the NanoClave with normal saline or in accordance with facility protocol. After blood use, the NanoClave can be flushed clean and does not require change-out.
- › Use routine flushing in accordance with facility protocol in order to maintain catheter patency.
- › Change NanoClave in accordance with facility protocol and CDC Guidelines.

Functional Attributes

- › Lipid and Blood Compatible
- › Radiographic Imaging Compatible (CT Compatible)
- › Contains No Latex, Phthalate's (DEHP) or Metal Components, MRI Compatible



NanoClave's saline flush option can eliminate the risk of Heparin Induced Thrombocytopenia (HIT).

icumedical
human connections