



Product Reference Guide

Tego™

Needlefree Haemodialysis Connector

Attach Tego to Patient's Catheter



- › Use aseptic technique to remove Tego from package.
- › Prime Tego in accordance with facility procedure/protocol.
- › Remove cap and attach Tego to desired vascular access device.

Swab Tego Before Each Access



- › Swab silicone seal of Tego before each use in accordance with facility procedure/protocol.
- › Do not use needles, blunt cannulas, or luer-type caps to access Tego.

Access Tego



- › Attach administration devices or syringes by pushing straight into Tego access devices for infusions.
- › When using Luer slip connection, only insert halfway while turning slightly clockwise. When removing, apply same clockwise turning motion. Do not overtighten.

Flush Tego



- › Clamp line before disconnecting from Tego and between dialysis sessions per established clinical and catheter manufacturer's practices.
- › Dispose of Tego in accordance with facility procedure/protocol.
- › Precautions: Device should be changed in accordance with current, recognised guidelines of IV therapy. Tego has been validated for functional and microbiological efficacy for up to seven days.

Functional Attributes

- › Straight fluid path accommodates flow rates of greater than 600 mL/min.
- › Remains in place during the entire haemodialysis treatment period.
- › Validated for continuous patient use for up to seven days.
- › Silicone seal remains closed when not activated, closing the fluid path and eliminating blood exposure.



Tego's saline flush option is designed to help reduce the risk of heparin-induced thrombocytopenia.

icumedical
human connections