

Gripper™

Port Access Needle



Non-Safety Huber Needle

- › A non-safety port access needle to administer or withdraw fluids from implanted ports
- › Available in different sizes and lengths
- › Small footprint base with thick cushioned needle platform sits against the skin when in-situ
- › Removable contoured grip for easy insertion
- › Non-siliconized, non-coring needle
- › The 20 cm-long extension set is available with or without Y-site with removable injection cap on luer lock
- › Extension set with color-coded clamp(s)
- › Priming volume: 0.3 mL or 0.5 mL, depending on the extension set is with or without the Y-site
- › Pressure-rated up to 45 psi

Ordering Information

Y-site — Removable injection cap on luer-lock Y site

Product Code	Gauge	Length	Units
21-2938-24	● 22	5/8 in	12/box
21-2939-24		3/4 in	
21-2940-24		1 in	
21-2941-24		1 1/4 in	
21-2946-24	● 20	5/8 in	12/box
21-2947-24		3/4 in	
21-2948-24		1 in	
21-2949-24		1 1/4 in	
21-2955-24	● 19	3/4 in	12/box
21-2956-24		1 in	
21-2957-24		1 1/4 in	

Without Y-site

Product Code	Gauge	Length	Units
21-2737-24	● 22	5/8 in	12/box
21-2733-24		3/4 in	
21-2714-24		1 in	
21-2715-24		1 1/4 in	
21-2736-24	● 20	5/8 in	12/box
21-2734-24		3/4 in	
21-2717-24		1 in	
21-2718-24		1 1/4 in	
21-2735-24	● 19	3/4 in	12/box
21-2720-24		1 in	
21-2721-24		1 1/4 in	

