



Rio™

Drug Reconstitution Transfer Device

Facilitates preparation of IV medication admixtures when and where you need them

icumedical
human connections

Meeting your admixture needs in the pharmacy or at the point of care

- For reconstitution of liquid or lyophilized 20 mm standard drug vials and aseptic transfer into ICU Medical's multi-port LifeCare™ IV container systems (up to 250 mL)
- Improve operational flexibility with immediate or future activation
- After reconstitution, the stopcock allows the fluid pathway to be closed, preventing reconstituted drug from re-entering the vial
- Needlefree design helps reduce needlestick related injuries and sharps disposal costs

Latex-free | Non-PVC/Non-DEHP
Single-use only | Individual sterile packaging

Product Information	Item #	Description	Case Qty
		VB-20	Rio - Drug Reconstitution Transfer Device



In-Service Video

Your guide to the safe and effective use of the Rio Drug Reconstitution Transfer Device



To learn how to put the Rio drug reconstitution transfer device to work for you, please call 866.829.9025 or visit www.icumed.com

icumedical
human connections