

500 mL

RX ONLY

List No. 7951-13

# EMPTY CONTAINER

*SINGLE-USE I.V. CONTAINER  
FOR ADMIXING  
COMPATIBLE MEDICATIONS.*

FLUID PATH AND AREAS BENEATH UNDISTURBED PROTECTIVE COVERS ARE STERILE AND NONPYROGENIC IN INTACT UNIT PACKAGE. USE ASEPTIC TECHNIQUE. MIX CONTAINER CONTENTS THOROUGHLY. ALL GRADUATIONS ARE APPROXIMATE. LABEL APPROPRIATELY. BEFORE ADMINISTRATION USE, CHECK FOR MINUTE LEAKS BY SQUEEZING CONTAINER FIRMLY. DO NOT USE IF CONTAINER IS DAMAGED. ADHERE TO STORAGE REQUIREMENTS OF ADDED MEDICATIONS. MUST NOT BE USED IN SERIES CONNECTIONS. MADE IN COSTA RICA. PRODUCT INQUIRIES SHOULD BE DIRECTED TO ICU MEDICAL, INC., LAKE FOREST, ILLINOIS, 60045, USA



Not made with  
natural rubber latex



CONTAINS DEHP  
Recyclable as Type-3 Plastic

IM-4364 (06/2018)

**icumedical**