




ICU Medical, Inc.

♥ CRITICAL CARE SYSTEMS

Because positive
outcomes depend on
reliable monitoring.



TRANSPAC® IV

DISPOSABLE PRESSURE TRANSDUCERS

Transpac® IV gives you the reliable hemodynamic and
cardiovascular data you need to manage your critical patients

Advancing reliability and flexibility in patient monitoring.



Clear fluid pathway
facilitates setup
and priming



Two proven continuous flush devices to choose from
Intraflo® flush



Squeeze flush



**Patient or
pole mounting**
for clinical flexibility



Easy-to-use
"telephone style"
cable connection

Multiple Standard Configurations

- Adult and neonatal/pediatric
- Intracranial and intercompartmental
- Various single-, double-, and triple-line options
- Large made-to-order capability

Available Safety Enhancements

- SAFESSET™ integral closed blood sampling and blood conservation system
- Needle-free systems

Transpac® IV Product Specifications

Operating Pressure: -50 to 300 mm Hg
Overpressure Limits: -400 to 5000 mm Hg
Sensitivity: 5 microvolt/Volt/mm Hg \pm 1%
@ 6.0 VDC and 22°C

Zero Offset: <25 mm Hg
Zero Drift: <2 mm Hg in 8 hours
Input Impedance: 300-350 ohms
Output Impedance: 300 ohms Nominal

Exceeds AAMI Requirements for Accuracy

TRANSPAC® IV

DISPOSABLE PRESSURE TRANSDUCERS



ICU Medical, Inc.

951 Calle Amanecer, San Clemente, CA 92673, U.S.A.
Tel: +1 949 366 2183 • www.icumed.com

M1-1224 Rev. 01