



**ICU Medical, Inc.**

♥ CRITICAL CARE SYSTEMS

Because positive outcomes  
depend on the safety  
of your interventions.




**SAFESET™**  
BLOOD SAMPLING SYSTEM

**SAFESET™ is a closed, needle-free system that protects both you and your patients**

- Protects against exposure to bloodborne pathogens and IV line contamination
- Reduces blood waste and therefore may help reduce the risk of iatrogenic anemia
- May increase the efficiency and safety of the blood-sampling process
- Available integrated with the Transpac® IV Pressure Monitoring Kit or as a stand-alone system

# Advancing safety for patients and providers.



## In-line Integration

- Safe, convenient blood withdrawal from pressure-monitoring lines
  - Engineered to optimize dynamic response when used with Transpac® IV pressure-monitoring kit
- May reduce risk and expense associated with blood waste
  - Discard or clearing volume can be reinfused, which may reduce blood transfusions due to phlebotomy



## Multiple Standard Configurations

- Adult and neonatal/pediatric
- Various single-, double-, and triple-line options
- Large made-to-order capability


## SAFESET™ Product Specifications

Pressure tubing: Internal volume, .055 mL/in  
 Reservoir: Adult, 12 mL (graduated to 10 mL); Neonatal/pediatric, 3 mL (graduated)  
 Latex free: Yes


**SAFESET™**  
 BLOOD SAMPLING SYSTEM

## SAFESET™ Accessories

42303-02  
 Shielded Blunt Cannula



42327-02  
 Shielded Blunt Cannula with Tube Holder



To place an order, please contact your ICU Medical Product Specialist at 1-800-824-7890.

**CAUTION:** Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician or other licensed practitioner.