



In Living Color

The ICU Medical range of polycarbonate syringes answers the demand for top quality latex free syringes with custom color plungers and print. We use crystal clear polycarbonate for excellent visibility. The construction is smooth and sturdy for a "glass-like" syringe feel.

Improving Safety



Our syringes feature pre-printed and blank barrels, seven different color coded plungers, and six different sizes. Along with our other design benefits, these features add up to improved safety by minimizing medication errors. The blank syringes are available in stand-alone configurations as well as ready to be put into a custom kit. Pre-printed and custom labeled syringes are available to fulfill your custom kit needs.

Catalog #	Size	Color
42900-01	1 ml	White
42900-02	1 ml	Blue
42900-03	1 ml	Yellow
42900-05	1 ml	Red
42901-01	3 ml	White
42901-02	3 ml	Blue
42901-03	3 ml	Yellow
42901-05	3 ml	Red
42902-01	6 ml	White
42902-02	6 ml	Blue
42902-05	6 ml	Red
42903-01	10 ml	White
42903-02	10 ml	Blue
42903-03	10 ml	Yellow
42903-05	10 ml	Red
42903-06	10 ml	Light Green
42903-07	10 ml	Purple
42903-08	10 ml	Dark Green
42904-01	20 ml	White
42904-02	20 ml	Blue
42904-03	20 ml	Yellow
42904-05	20 ml	Red
42905-01	30 ml	White
42905-02	30 ml	Blue
42905-07	30 ml	Purple

Sold 25 per box