



RE: Use of the Antimicrobial CLAVE® and MicroCLAVE® Connector with MRI

Thank you for your interest in the Antimicrobial CLAVE devices. The antimicrobial agent that is used in these products is a silver ion, not a silver metal. The silver ion has no ferromagnetic activity and therefore will not have either a magnetic or heat related reaction within the MRI chamber or procedure. Furthermore, studies on products like wound care dressings which use the same type of silver ions, have confirmed that these types of silver ions do not cause a disruption in the image itself.

Thank you again for your interest, if there are further questions or concerns please visit our website at www.icumed.com, or contact the corporate offices at 949-366-2183.

Technical Services
ICU Medical, Inc.