



RE: Stopcock Compatibility with CLAVE® Connector

Thank you for your interest in the CLAVE Connector by ICU Medical. As you are aware some stopcocks do not operate favorably with the CLAVE products and ICU makes best efforts to keep an updated list of those which are fully compatible. At times manufactures may make changes in the design to stopcocks which could affect compatibility and we appreciate your bringing this to our attention. In regards to Medex manufactured stopcocks ICU Medical verifies that the following part numbers are compatible with the CLAVE and MicroCLAVE.

Small Body Stopcocks:

- MX531-1L Three-way with swivel luer
- MX543-1LN Four-way with swivel luer and end cap
- MX531 FL Three-way with fixed luer lock
- MX534 FL Four-way with fixed luer lock

Large Body Stopcocks:

- MX411 Three-way with luer slip
- MX421 Four-way with luer slip
- MX431-1L Three-way with fixed luer lock
- MX434-1L Four-way with fixed luer lock

Thank you again for your interest, if there are further questions or concerns please visit our website at www.icumed.com, or contact the corporate offices at 949-366-2183 or 800-824-7890.

Technical Services
ICU Medical, Inc.