

# Manage hazardous drug vapor your way

Standardise on a preferred vial adapter technology that best fits your facility's needs when considering cost, workflow and ease of use.

Air-capture  
barrier type CSTD

Air-cleaning  
filter type CSTD

Access a broad range of vial sizes and stopper types.



20 mm



28 mm



20 mm or 28 mm



Skirted vial access devices securely lock on to 20 mm and 28 mm vials.

Universal vial adapter allows access to 20 mm or 28 mm vials with a single product, helping you reduce SKUs.

Longer fluid path on spike supports more thorough draining and reduced drug waste.

## Air-capture barrier type CSTD

Product Code	Description
011-CL-80	ChemoLock Closed Vial Spike
011-CL-80S	ChemoLock Vial Spike, 20 mm
011-CL-80SL	ChemoLock Vial Spike 28 mm

## Air-cleaning filter type CSTD

Product Code	Description
011-CL-70	ChemoLock™ Universal Vented Vial Spike
011-CL-70S	ChemoLock Vented Vial Spike 20 mm
011-CL-72	ChemoLock Vented Vial Spike, 13 mm

To learn more about ICU Medical's full line of closed system transfer devices (CSTDs), visit [www.icumed.com/CSTDs](http://www.icumed.com/CSTDs)

**icumedical**  
human connections