



Implantable Ports and Port Access Needles

Product Catalogue

icumedical
human connections

Contents

Implantable Access Systems

Port-A-Cath II Power P.A.C.™ Polysulfone/Titanium Implantable Venous Access Systems	3
Port-A-Cath II Power P.A.C.™ Systems with Valved Peel-Away Introducers	3
Port-A-Cath Power P.A.C.™ Titanium Implantable Venous Access Systems	4
P.A.S. Port™ T2 Power P.A.C. Titanium Implantable Peripheral Venous Access Systems	4
Port-A-Cath II™ Systems with Valved Peel-Away Introducers	5
Port-A-Cath II™ Polysulfone/Titanium Implantable Venous Access Systems	5
Port-A-Cath™ Titanium Implantable Venous Access Systems	6
P.A.S. Port™ T2 Titanium Implantable Peripheral Venous Access Systems	6
P.A.S. Port™ Elite Plastic/Titanium Implantable Venous Access Systems	6
Port-A-Cath™ System Accessories	7

Port Access Needles

Gripper Plus™ Power P.A.C. Safety Huber Needles	8
Gripper Micro™ Blunt Cannula, Non-Coring Safety Needles	9
Gripper Plus™ Safety Huber Needles	10
Gripper™ Huber Needles	11
Port-A-Cath™ Plastic Hub Huber Needles	12

Port-A-Cath II Power P.A.C. Polysulfone/Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Polysulfone/Titanium Portals with PolyFlow™ Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4452-24 Kit	1 0610586 03267 7	each
						21-4453-24 Tray	1 0610586 03268 4	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4454-24 Kit	1 0610586 03269 1	each
						21-4455-24 Tray	1 0610586 03270 7	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4482-24 Kit	10610586 03276 9	each
						21-4483-24 Tray	1 0610586 03277 6	each
Single-lumen Low Profile Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4470-24 Kit	1 0610586 03273 8	each
						21-4471-24 Tray	1 0610586 03274 5	each
Dual-lumen Low Profile Portal	7 Fr	Polyurethane	2.2	1.0/1.0	76	21-8466-24 Tray	1 0610586 03283 7	each
Dual-lumen Low Profile Portal	10 Fr	Polyurethane	3.2	1.4/1.4	76	21-8468-24 Tray	1 0610586 03285 1	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, and Point-Lok™ device.

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus Power P.A.C. safety needle, 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

Port-A-Cath II Power P.A.C. Systems with Valved Peel-Away Introducers

Introducer sheath length – 13 cm – Dilator sheath length – 18 cm



Port-A-Cath Power P.A.C. Systems with Valved Peel-Away Introducers									
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Connector	Reorder No.	GTIN	Unit
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	Ultra-Lock™	21-4474-24 Tray	1 0610586 03996 6	each
Single-lumen Low Profile Portal	8.5 Fr	Polyurethane	2.6	1.6	76	Ultra-Lock	21-4475-24 Tray	1 0610586 03997 3	each
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	Ultra-Lock	21-4476-24 Tray	1 0610586 03998 0	each
Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	1.0	76	Ultra-Lock	21-4477-24 Tray	1 0610586 03999 7	each
Dual-lumen Low Profile Portal	7 Fr	Polyurethane	2.2	1.0/1.0	76	Wing-Lock™	21-8469-24 Tray	1 0610586 04000 9	each
Dual-lumen Low Profile Portal	10 Fr	Polyurethane	3.2	1.4/1.4	76	Wing-Lock	21-8470-24 Tray	1 0610586 04001 6	each

Valved peel-away introducers are available in ICU Medical's Port-A-Cath™ system trays and also as separate sets.

Port-A-Cath Power P.A.C. Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Titanium Contoured Portals with PolyFlow Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4424-24 Kit	10610586032639	each
						21-4425-24 Tray	10610586032646	each
Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4436-24 Kit	10610586032653	each
						21-4437-24 Tray	10610586032653	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4422-24 Kit	10610586032615	each
						21-4423-24 Tray	10610586032622	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, and Point-Lok device.

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus Power P.A.C. safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

P.A.S. Port T2 Power P.A.C. Titanium Implantable Peripheral Venous Access Systems

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Titanium Contoured Portals with PolyFlow Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Peripheral Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4872-24 Kit	10610586032806	each
						21-4873-24 Tray	10610586032813	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, and Point-Lok device.

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus Power P.A.C. safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

Port-A-Cath II Systems with Valved Peel-Away Introducers



Port-A-Cath II Systems with Valved Peel-Away Introducers Introducer sheath length – 13 cm – Dilator sheath length – 18 cm									
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Connector	Reorder No.	GTIN	Unit
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	Ultra-Lock	21-4066-24 Tray	1 0610586 03990 4	each
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	Ultra-Lock	21-4067-24 Tray	1 0610586 03991 1	each
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	1.6	76	Ultra-Lock	21-4068-24 Tray	1 0610586 03992 8	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	Ultra-Lock	21-4069-24 Tray	1 0610586 03993 5	each

Valved peel-away introducers are available in ICU Medical's Port-A-Cath system trays and also as separate sets.

Port-A-Cath II Polysulfone/Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla



Polysulfone/Titanium Portals with PolyFlow Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4052-24 Kit	1 0610586 02352 1	each
						21-4053-24 Tray	1 0610586 02353 8	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4054-24 Kit	1 0610586 02354 5	each
						21-4055-24 Tray	1 0610586 02355 2	each
Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4082-24 Kit	1 0610586 02385 9	each
						21-4083-24 Tray	1 0610586 02361 3	each
Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4070-24 Kit	1 0610586 01782 7	each
						21-4071-24 Tray	1 0610586 01784 1	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick and Point-Lok device(s).

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle, 2 syringes, tunneling tool and Point-Lok device.

All systems are unassembled.

Polysulfone/Titanium Portals with Silicone Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	9 Fr	Silicone	2.8	1.0	76	21-4050-24 Kit	1 0610586 02350 7	each
						21-4051-24 Tray	1 0610586 02351 4	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick and Point-Lok device(s).

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle, 2 syringes, tunneling tool and Point-Lok device.

All systems are unassembled.

Port-A-Cath Titanium Implantable Venous Access Systems

MR Conditional up to 3-Tesla



Titanium Contoured Portals with PolyFlow Polyurethane Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4024-24 Kit	1 0610586 02379 8	each
						21-4025-24 Tray	1 0610586 02380 4	each
Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4022-24 Kit	1 0610586 02377 4	each
						21-4023-24 Tray	1 0610586 02378 1	each
Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4036-24 Kit	1 0610586 02383 5	each
						21-4037-24 Tray	1 0610586 02384 2	each

Titanium Contoured Portals with Silicone Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	9 Fr	Silicone	2.8	1.0	76	21-4008-24 Kit	1 0610586 02374 3	each
						21-4009-24 Tray	1 0610586 02375 0	each
Single-lumen Low Profile Portal	9 Fr	Silicone	2.8	1.0	76	21-4034-24 Kit	1 0610586 02381 1	each
						21-4035-24 Tray	1 0610586 02382 8	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, and Point-Lok device(s).

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

P.A.S. Port T2 Titanium Implantable Peripheral Venous Access Systems

MR Conditional up to 3-Tesla



Titanium Portal with Polyurethane Catheter								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4572-24 Kit	1 0610586 02516 7	each
						21-4573-24 Tray	1 0610586 02517 4	each

P.A.S. Port Elite – Plastic/Titanium Implantable Peripheral Venous Access System



Plastic/ Titanium Portal with Polyurethane Catheter								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4591-24 Kit	1 0610586 02154 1	each
						21-4591-24 Tray	1 0610586 02156 5	each

Port-A-Cath System Accessories



Non-Valved Peel-Away Introducer Sets					
	Sheath Length (cm)	Dilator Length (cm)	Reorder No.	GTIN	Unit
6 French Includes two guidewires	16	19	21-2320-24	3 0610586 02511 6	5/box
7 French	16	19	21-2471-24*	3 0610586 04015 7	5/box
8.5 French	16	19	21-2470-24	3 0610586 02514 7	5/box
9 French	16	19	21-2340-24	3 0610586 02512 3	5/box
10 French	16	19	21-2472-24*	3 0610586 04016 4	5/box
11 French	16	19	21-2386-24	3 0610586 02513 0	5/box

Non-valved peel-away introducers are available in ICU Medical's Port-A-Cath™ system trays and also as separate sets.

*Subject to availability



Valved Peel-Away Introducer Sets					
	Sheath Length (cm)	Dilator Length (cm)	Reorder No.	GTIN	Unit
6 French	13	18	21-2476-24*	3 0610586 04020 1	5/box
7 French	13	18	21-2477-24*	3 0610586 04021 8	5/box
8.5 French	13	18	21-2478-24*	3 0610586 04022 5	5/box
9 French	13	18	21-2479-24*	3 0610586 04023 2	5/box
10 French	13	18	21-2480-24*	3 0610586 04024 9	5/box
11 French	13	18	21-2481-24*	3 0610586 04025 6	5/box

*Subject to availability



Subcutaneous Tunneling Device			
	Reorder No.	GTIN	Unit
3.2 mm OD x 24 cm length, solid malleable stainless steel, disposable	21-0125-24	3 0610586 01319 9	5/box

Introducer sets include: 18G extra thin wall introducer needle, "J" guidewire, vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle (except 21-2320), 2 syringes, and Point-Lok device.

Gripper Plus Power P.A.C. Safety Huber Needles

MR Conditional up to 3-Tesla, power injection up to 5 mL/sec and 300 psi



Gripper Plus Power P.A.C. Safety Huber Needles Not made with natural rubber latex, non-coring, non-siliconised, removable contoured grip, cushioned needle platform, 20 cm extension tubing* with or without Y-site				
Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	GTIN	Unit
● 20 gauge	16 mm	21-3464-24	3 0610586 03363 6	12/box
	19 mm	21-3465-24	3 0610586 03304 3	12/box
	25 mm	21-3466-24	3 0610586 03305 0	12/box
	32 mm	21-3467-24	3 0610586 03306 7	12/box
● 19 gauge	19 mm	21-3468-24	3 0610586 03307 4	12/box
	25 mm	21-3469-24	3 0610586 03308 1	12/box
	32 mm	21-3470-24	3 0610586 03309 8	12/box
Without Y-site	Length	Reorder No.	GTIN	Unit
● 20 gauge	16 mm	21-3369-24	3 0610586 03302 9	12/box
	19 mm	21-3367-24	3 0610586 03300 5	12/box
	25 mm	21-3362-24	3 0610586 03296 1	12/box
	32 mm	21-3363-24	3 0610586 03297 8	12/box
● 19 gauge	19 mm	21-3368-24	3 0610586 03301 2	12/box
	25 mm	21-3364-24	3 0610586 03298 5	12/box
	32 mm	21-3365-24	3 0610586 03299 2	12/box

* Manufactured with Non-DEHP plasticised tubing

Gripper Micro Blunt Cannula, Non-Coring Safety Needles

MR Conditional up to 3-Tesla



Gripper Micro Blunt Cannula, Non-Coring Safety Needles Not made with natural rubber latex, non-coring, non-siliconised, cushioned infusion site, 20 cm extension tubing* with or without Y-site				
Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	GTIN	Unit
● 22 gauge	19 mm	21-3271-24	3 0610586 02988 6	12/box
	25 mm	21-3272-24	3 0610586 02989 3	12/box
	32 mm	21-3273-24	3 0610586 02990 9	12/box
● 20 gauge	19 mm	21-3276-24	3 0610586 02993 0	12/box
	25 mm	21-3277-24	3 0610586 02994 7	12/box
	32 mm	21-3278-24	3 0610586 02995 4	12/box
● 19 gauge	19 mm	21-3280-24	3 0610586 02997 8	12/box
	25 mm	21-3281-24	3 0610586 02998 5	12/box
	32 mm	21-3282-24	3 0610586 02999 2	12/box
Without Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	19 mm	21-3251-24	3 0610586 02974 9	12/box
	25 mm	21-3252-24	3 0610586 02975 6	12/box
	32 mm	21-3253-24	3 0610586 02976 3	12/box
● 20 gauge	19 mm	21-3256-24	3 0610586 02979 4	12/box
	25 mm	21-3257-24	3 0610586 02980 0	12/box
	32 mm	21-3258-24	3 0610586 02981 7	12/box
● 19 gauge	19 mm	21-3260-24	3 0610586 02983 1	12/box
	25 mm	21-3261-24	3 0610586 02984 8	12/box
	32 mm	21-3262-24	3 0610586 02985 5	12/box

* Manufactured with Non-DEHP plasticized tubing

Gripper Plus Safety Huber Needles

MR Conditional up to 3-Tesla



Gripper Plus Safety Huber Needles Not made with natural rubber latex, non-coring, non-siliconised, removable contoured grip, cushioned needle platform, 20 cm extension tubing* with or without Y-site				
Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2960-24	3 0610586 02768 4	12/box
	19 mm	21-2961-24	3 0610586 02769 1	12/box
	25 mm	21-2962-24	3 0610586 02770 7	12/box
	32 mm	21-2963-24	3 0610586 02771 4	12/box
● 20 gauge	16 mm	21-2964-24	3 0610586 02772 1	12/box
	19 mm	21-2965-24	3 0610586 02773 8	12/box
	25 mm	21-2966-24	3 0610586 02774 5	12/box
	32 mm	21-2967-24	3 0610586 02775 2	12/box
● 19 gauge	19mm	21-2968-24	3 0610586 02776 9	12/box
	25 mm	21-2969-24	3 0610586 02777 6	12/box
	32 mm	21-2970-24	3 0610586 02778 3	12/box
Y-site Removable split septum injection cap on Luer-lock Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2860-24	3 0610586 02564 2	12/box
	19 mm	21-2861-24	3 0610586 02565 9	12/box
	25 mm	21-2862-24	3 0610586 02566 6	12/box
	32 mm	21-2863-24	3 0610586 02567 3	12/box
● 20 gauge	16 mm	21-2864-24	3 0610586 02568 0	12/box
	19 mm	21-2865-24	3 0610586 02569 7	12/box
	25 mm	21-2866-24	3 0610586 02570 3	12/box
	32 mm	21-2867-24	3 0610586 02571 0	12/box
● 19 gauge	19 mm	21-2868-24	3 0610586 02572 7	12/box
	25 mm	21-2869-24	3 0610586 02573 4	12/box
	32 mm	21-2870-24	3 0610586 02574 1	12/box
Without Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2770-24	3 0610586 02563 5	12/box
	19 mm	21-2766-24	3 0610586 02559 8	12/box
	25 mm	21-2760-24	3 0610586 02553 6	12/box
	32 mm	21-2761-24	3 0610586 02554 3	12/box
● 20 gauge	16 mm	21-2769-24	3 0610586 02562 8	12/box
	19 mm	21-2767-24	3 0610586 02560 4	12/box
	25 mm	21-2762-24	3 0610586 02555 0	12/box
	32 mm	21-2763-24	3 0610586 02556 7	12/box
● 19 gauge	19 mm	21-2768-24	3 0610586 02561 1	12/box
	25 mm	21-2764-24	3 0610586 02557 4	12/box
	32 mm	21-2765-24	3 0610586 02558 1	12/box

* Manufactured with Non-DEHP plasticised tubing

Gripper Huber Needles

MR Conditional up to 3-Tesla



Gripper Huber Needles

Not made with natural rubber latex, non-coring, non-siliconised, removable contoured grip, cushioned needle platform, 20 cm extension tubing* with or without Y-site

Y-site Removable injection cap on Luer-lock Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2938-24	3 0610586 02321 1	12/box
	19 mm	21-2939-24	3 0610586 01674 9	12/box
	25 mm	21-2940-24	3 0610586 01675 6	12/box
	32 mm	21-2941-24	3 0610586 01676 3	12/box
● 20 gauge	16 mm	21-2946-24	3 0610586 02323 5	12/box
	19 mm	21-2947-24	3 0610586 01677 0	12/box
	25 mm	21-2948-24	3 0610586 01678 7	12/box
	32 mm	21-2949-24	3 0610586 01679 4	12/box
● 19 gauge	19 mm	21-2955-24	3 0610586 01680 0	12/box
	25 mm	21-2956-24	3 0610586 01681 7	12/box
	32 mm	21-2957-24	3 0610586 01682 4	12/box
Without Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2737-24	3 0610586 02320 4	12/box
	19 mm	21-2733-24	3 0610586 01584 1	12/box
	25 mm	21-2714-24	3 0610586 01578 0	12/box
	32 mm	21-2715-24	3 0610586 01579 7	12/box
● 20 gauge	16 mm	21-2736-24	3 0610586 02322 8	12/box
	19 mm	21-2734-24	3 0610586 01585 8	12/box
	25 mm	21-2717-24	3 0610586 01580 3	12/box
	32 mm	21-2718-24	3 0610586 01581 0	12/box
● 19 gauge	19 mm	21-2735-24	3 0610586 01586 5	12/box
	25 mm	21-2720-24	3 0610586 01582 7	12/box
	32 mm	21-2721-24	3 0610586 01583 4	12/box

* Manufactured with Non-DEHP plasticised tubing

Port-A-Cath Plastic Hub Huber Needles

MR Conditional up to 3-Tesla



Port-A-Cath Plastic Hub Huber Needles Non-coring, non-siliconised				
Straight	Length	Reorder No.	GTIN	Unit
● 22 gauge	25 mm	21-2005-24	3 0610586 02399 0	12/box
	38 mm	21-2007-24	3 0610586 02401 0	12/box
● 20 gauge	25 mm	21-2013-24	3 0610586 02403 4	12/box
	38 mm	21-2015-24	3 0610586 02405 8	12/box
● 19 gauge	25 mm	21-2021-24	3 0610586 02407 2	12/box
	51 mm	21-2027-24	3 0610586 02846 9	12/box
90° Bent	Length	Reorder No.	GTIN	Unit
● 22 gauge	13 mm	21-2316-24	3 0610586 02411 9	12/box
	25 mm	21-2006-24	3 0610586 02400 3	12/box
	38 mm	21-2008-24	3 0610586 02402 7	12/box
● 20 gauge	13 mm	21-2317-24	3 0610586 02412 6	12/box
	19 mm	21-2287-24	3 0610586 02408 9	12/box
	25 mm	21-2014-24	3 0610586 02404 1	12/box
	38 mm	21-2016-24	3 0610586 02406 5	12/box
● 19 gauge	25 mm	21-2022-24	3 0610586 02805 6	12/box
	38 mm	21-2026-24	3 0610586 02806 3	12/box

The product complies with current legislation and has the corresponding CE marking.
For additional information, warnings and / or safety precautions, refer to the manufacturer’s user manual.