

NanoClave® for Use with Magnetic Resonance Imaging (MRI)

The NanoClave is founded in original MicroClave technology, which is the most widely used needlefree connector worldwide. The NanoClave, including the NanoClave Manifold and NanoClave Stopcock, is comprised of medical grade plastics and silicone, which are fully compatible with the MRI procedure, such that they will in no way disrupt the image. Tens of thousands of technologically similar MicroClaves have been in use during MRI procedures in the past decade, with no reported complaints of image disruption.

If there are further questions or concerns, please visit our website at www.icumed.com or contact the corporate offices at 949-366-2183 or 800-824-7890.

Technical Services
ICU Medical Inc.