

## Use of the Tego<sup>®</sup> Connector with MRI

The Tego is a Closed Needlefree Access device specifically designed for use with Hemodialysis Central Venous Catheters and related applications. The Tego is comprised of medical grade plastics and silicone, which are fully compatible with the MRI procedure, such that they will in no way disrupt the image. Tens of thousands of technologically similar Tegos have been in use during MRI procedures in the past decade, with no reported complaints of image disruption. Since there are no metal components, the Tego will not have any magnetic effect within the MRI chamber or procedure.

If there are further questions or concerns, please visit our website at [www.icumed.com](http://www.icumed.com) or contact the corporate offices at 949-366-2183 or 800-824-7890.

Technical Services  
ICU Medical Inc.