

ICU Medical MedNet™ and Connexall® Seamless Integration

Combining the power of ICU Medical MedNet safety software and Connexall alarm management platform in your infusion therapy practice provides intelligent alert routing and enables clinicians to prioritize alarms, contributing to a reduction of alarm fatigue and helping you improve clinician awareness and response times.



Plum 360



LifeCare PCA



Connexall®



This integration is a powerful tool for allowing a network-based alarm message from the Plum 360™ or LifeCare PCA™ infusion pumps via ICU Medical MedNet safety software to be forwarded to Connexall for intelligent routing-based customized workflows.

Primary features of this integration include:

- > Clinician assignment
- > Rules-based secondary notification to caregivers
- > Intelligent escalation of alerts
- > Reporting platform

To learn more about ICU Medical IV smart infusion pumps with MedNet safety software, please call 800.824.7890 or visit www.icumed.com

To learn more about the Connexall alarm management and event notification platform, please call 866.556.3377 or visit www.connexall.com

Connexall®

icumedical
human connections