



Bivona™ Silicone Tracheostomy Tubes

Product Catalog

Contents

Introduction	4
Bivona™ TTS™	5
Bivona TTS Adult Tracheostomy Tubes (Cuffed)	5
Bivona TTS Neonatal Tracheostomy Tubes (Cuffed "V" Neck Flange)	5
Bivona TTS Neonatal Tracheostomy Tubes (Cuffed Straight Neck Flange)	5
Bivona TTS Pediatric Tracheostomy Tubes (Cuffed "V" Neck Flange)	5
Bivona TTS Pediatric Tracheostomy Tubes (Cuffed Straight Neck Flange)	6
Bivona™ TTS™ FlexTend™	7
Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuffed "V" Neck Flange)	7
Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuffed Straight Neck Flange)	7
Bivona TTS Pediatric FlexTend Tracheostomy Tubes (Cuffed "V" Neck Flange)	7
Bivona TTS Pediatric FlexTend Tracheostomy Tubes (Cuffed Straight Neck Flange)	7
Bivona™ TTS™ FlexTend Plus™	8
Bivona TTS Pediatric FlexTend Plus Tracheostomy Tubes (Cuffed "V" Neck Flange)	8
Bivona TTS Pediatric FlexTend Plus Tracheostomy Tubes (Cuffed Straight Neck Flange)	8
Bivona™ FlexTend Plus™	9
Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)	9
Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes (Extra Length, Straight Neck Flange)	9
Bivona Cuffless Pediatric FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)	9
Bivona Cuffless Pediatric FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)	9
Bivona Cuffless FlexTend Pediatric Tracheostomy Tubes (Straight Neck Flange, Extra Length)	10
Bivona Cuffless FlexTend Pediatric Tracheostomy Tubes (Straight Neck Flange, Extra Length)	10
Bivona™ Aire-Cuf™	11
Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes	11
Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes with Talk Attachment	11
Bivona Aire-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)	11
Bivona Aire-Cuf Neonatal Tracheostomy Tubes (Straight Neck Flange)	12
Bivona Aire-Cuf Pediatric Tracheostomy Tubes ("V" Neck Flange)	12
Bivona Aire-Cuf Pediatric Tracheostomy Tubes (Straight Neck Flange)	12
Bivona™ Fome-Cuf™	13
Bivona Fome-Cuf Adult Tracheostomy Tubes	13
Bivona Fome-Cuf Adult Tracheostomy Tubes (with Talk Attachment)	13
Bivona Fome-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)	13
Bivona Fome-Cuf Pediatric Tracheostomy Tubes ("V" Neck Flange)	13

Bivona™ Adjustable Hyperflex™	14
Bivona TTS Adjustable Neck Flange Hyperflex Tracheostomy Tubes.....	14
Bivona Mid-Range Aire-Cuf Adjustable Neck Flange Hyperflex Tracheostomy Tubes.....	14
Bivona Cuffless Tube Adjustable Neck Flange Hyperflex Pediatric “V” Neck Flange Tracheostomy Tubes.....	14
Bivona™ Fixed Flange Hyperflex™	15
Bivona Cuffless Extra Length Fixed Hyperflex Tracheostomy Tubes	15
Bivona TTS Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes.....	15
Bivona Cuffed Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes.....	15
Bivona™ Cuffless.....	16
Bivona Cuffless Adult Tracheostomy Tubes	16
Bivona Cuffless Neonatal Tracheostomy Tubes (“V” Neck Flange).....	16
Bivona Cuffless Neonatal Tracheostomy Tubes (Straight Neck Flange)	16
Bivona Cuffless Pediatric Tracheostomy Tubes (“V” Neck Flange)	16
Bivona Cuffless Pediatric Tracheostomy Tubes (Straight Neck Flange)	17
Bivona™ Customized Silicone Tracheostomy Tubes.....	18
Accessories.....	18

Each tube comes individually packaged with obturator, disconnection wedge, and twill tape.

Bivona™ Silicone Tracheostomy Tubes

Silicone offers important advantages over conventional tube materials:

Softness and flexibility

For comfort, ease of use, and greater range of mobility.

Biocompatibility

Hydrophobic and hypoallergenic, tissue-friendly, and easy to clean.

Durability

Robust for reliable reprocessing, withstands repeated sanitization up to 10 times (adult) and 5 times (for neonatal and pediatrics).



Bivona silicone tracheostomy tubes are available in a range of tube configuration options to suit each, age, size, and application requirements. Customized Bivona silicone tubes are available.

For comfort

A variety of tubes designed to aid in patient comfort:

- › Allows 15 mm connector to be distal from chin.
- › Reduces the risk of device-related pressure ulcers.
- › Option to utilize self-expanding cuff to mitigate risks of over inflation.

For design

Integrating unique features to aid patient comfort, ease of use and avoidance of contamination.

For usability and safety

Hydrophobic material for easy cleaning, free from DEHP, latex, and plasticizers.

For cost-efficiency

Bivona tubes can be reprocessed for single patient use up to 10 times for adult-sized tubes and up to 5 times for neonatal- and pediatric-sized tubes.



For choice

Combining the softness and flexibility of silicone, with an extensive range of size and shape options for all ages.

For safe and proper use, refer to the instructions for Use.

Bivona™ Tight-to-Shaft™ (TTS)

The TTS cuff design offers the benefits of a cuff with the profile of a cuffless tube. When the silicone cuff is totally deflated, it adds no significant dimension to the outer diameter of the tube's shaft, yet it can be partially or completely inflated for varied airway management needs.



Hospital



Home



MR Conditional



670180

Bivona TTS Adult Tracheostomy Tubes (Cuffed)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
670150	5.0	5.0	7.4	60.0	100°	9.5	1
670160	6.0	6.0	8.8	70.0	100°	11.0	1
670170	7.0	7.0	10.0	80.0	100°	12.0	1
670175	7.5	7.5	10.5	80.0	100°	12.5	1
670180	8.0	8.0	11.0	88.0	100°	13.0	1
670185	8.5	8.5	11.8	88.0	100°	14.0	1
670190	9.0	9.0	12.3	98.0	100°	14.5	1



67N030

Bivona TTS Neonatal Tracheostomy Tubes (Cuffed "V" Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67N025	2.5	2.5	4.0	30.0	120°	4.0	1
67N030	3.0	3.0	4.7	32.0	120°	4.7	1
67N035	3.5	3.5	5.3	34.0	120°	5.3	1
67N040	4.0	4.0	6.0	36.0	120°	6.0	1



67SN030

Bivona TTS Neonatal Tracheostomy Tubes (Cuffed Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67SN025	2.5	2.5	4.0	30.0	120°	4.0	1
67SN030	3.0	3.0	4.7	32.0	120°	4.7	1
67SN035	3.5	3.5	5.3	34.0	120°	5.3	1
67SN040	4.0	4.0	6.0	36.0	120°	6.0	1



67P045

Bivona TTS Pediatric Tracheostomy Tubes (Cuffed "V" Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67P025	2.5	2.5	4.0	38.0	120°	4.0	1
67P030	3.0	3.0	4.7	39.0	120°	4.7	1
67P035	3.5	3.5	5.3	40.0	120°	5.3	1
67P040	4.0	4.0	6.0	41.0	120°	6.0	1
67P045	4.5	4.5	6.7	42.0	120°	6.7	1
67P050	5.0	5.0	7.3	44.0	120°	7.3	1
67P055	5.5	5.5	8.0	46.0	120°	8.0	1

Bivona™ TTS™



Bivona TTS Pediatric Tracheostomy Tubes (Cuffed Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67SP025	2.5	2.5	4.0	38.0	120°	4.0	1
67SP030	3.0	3.0	4.7	39.0	120°	4.7	1
67SP035	3.5	3.5	5.3	40.0	120°	5.3	1
67SP040	4.0	4.0	6.0	41.0	120°	6.0	1
67SP045	4.5	4.5	6.7	42.0	120°	6.7	1
67SP050	5.0	5.0	7.3	44.0	120°	7.3	1
67SP055	5.5	5.5	8.0	46.0	120°	8.0	1

Bivona™ TTS™ FlexTend™

Radiopaque silicone tracheostomy tubes with a wire-reinforced extension, proximal to the neck flange. The silicone cuff lays flat against the shaft when deflated. Proximal extension allows 15 mm connector to be distal from the chin. The tube has a clear inflation balloon with a one-way valve and integral 15 mm connection. It is manufactured with Bivona SuperSlick™ coating.




Hospital



Home



MR
Conditional



67NFP25

Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuffed "V" Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67NFP25	2.5	2.5	4.0	20.0	30.0	120°	4.0	1
67NFP30	3.0	3.0	4.7	20.0	32.0	120°	4.7	1
67NFP35	3.5	3.5	5.3	20.0	34.0	120°	5.3	1
67NFP40	4.0	4.0	6.0	20.0	36.0	120°	6.0	1



67NFPS25

Bivona TTS Neonatal FlexTend Tracheostomy Tubes (Cuffed Straight Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67NFPS25	2.5	2.5	4.0	20.0	30.0	120°	4.0	1
67NFPS30	3.0	3.0	4.7	20.0	32.0	120°	4.7	1
67NFPS35	3.5	3.5	5.3	20.0	34.0	120°	5.3	1
67NFPS40	4.0	4.0	6.0	20.0	36.0	120°	6.0	1



67PFS25

Bivona TTS Pediatric FlexTend Tracheostomy Tubes (Cuffed "V" Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFS25	2.5	2.5	4.0	20.0	38.0	120°	4.0	1
67PFS30	3.0	3.0	4.7	20.0	39.0	120°	4.7	1
67PFS35	3.5	3.5	5.3	20.0	40.0	120°	5.3	1
67PFS40	4.0	4.0	6.0	20.0	41.0	120°	6.0	1
67PFS45	4.5	4.5	6.7	30.0	42.0	120°	6.7	1
67PFS50	5.0	5.0	7.3	30.0	44.0	120°	7.3	1
67PFS55	5.5	5.5	8.0	30.0	46.0	120°	8.0	1




67PFSS25

Bivona TTS Pediatric FlexTend Tracheostomy Tubes (Cuffed Straight Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFSS25	2.5	2.5	4.0	20.0	38.0	120°	4.0	1
67PFSS30	3.0	3.0	4.7	20.0	39.0	120°	4.7	1
67PFSS35	3.5	3.5	5.3	20.0	40.0	120°	5.3	1
67PFSS40	4.0	4.0	6.0	20.0	41.0	120°	6.0	1
67PFSS45	4.5	4.5	6.7	30.0	42.0	120°	6.7	1
67PFSS50	5.0	5.0	7.3	30.0	44.0	120°	7.3	1
67PFSS55	5.5	5.5	8.0	30.0	46.0	120°	8.0	1

Bivona™ TTS™ FlexTend Plus™

Radiopaque silicone tracheostomy tubes with a wire-reinforced extension, proximal to the neck flange. The silicone cuff lays flat against the shaft when deflated. Proximal extension allows 15 mm connector to be distal from the chin. The tube has a clear inflation balloon with a one-way valve and integral 15 mm connection. FlexTend Plus tubes also provide additional distal length for patients requiring additional overall length. It is manufactured with Bivona SuperSlick™ coating.



67PFP40

Bivona TTS Pediatric FlexTend Plus Tracheostomy Tubes (Extra Length, Cuffed "V" Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFP40	4.0	4.0	6.0	30.0	44.0	120°	6.0	1
67PFP45	4.5	4.5	6.7	30.0	48.0	120°	6.7	1
67PFP50	5.0	5.0	7.3	30.0	50.0	120°	7.3	1
67PFP55	5.5	5.5	8.0	30.0	52.0	120°	8.0	1
67PFP60	6.0	6.0	8.7	30.0	56.0	120°	8.7	1



67PFPS40

Bivona TTS Pediatric FlexTend Plus Tracheostomy Tubes (Extra Length, Cuffed Straight Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67PFPS40	4.0	4.0	6.0	30.0	44.0	120°	6.0	1
67PFPS45	4.5	4.5	6.7	30.0	48.0	120°	6.7	1
67PFPS50	5.0	5.0	7.3	30.0	50.0	120°	7.3	1
67PFPS55	5.5	5.5	8.0	30.0	52.0	120°	8.0	1
67PFPS60	6.0	6.0	8.7	30.0	56.0	120°	8.7	1

Bivona™ FlexTend™

Reinforced proximal and distal sections provide a high level of kink and crush resistance. Extended proximal tube allows 15 mm connector to be distal from the chin.



Hospital



Home



MR
Conditional



60NFP30

Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes ("V" Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60NFP25	2.5	2.5	4.0	20.0	30.0	120°	N/A	1
60NFP30	3.0	3.0	4.7	20.0	32.0	120°	N/A	1
60NFP35	3.5	3.5	5.3	20.0	34.0	120°	N/A	1
60NFP40	4.0	4.0	6.0	20.0	36.0	120°	N/A	1



60NFPS35

Bivona Cuffless Neonatal FlexTend Tracheostomy Tubes (Straight Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60NFPS25	2.5	2.5	4.0	20.0	30.0	120°	N/A	1
60NFPS30	3.0	3.0	4.7	20.0	32.0	120°	N/A	1
60NFPS35	3.5	3.5	5.3	20.0	34.0	120°	N/A	1
60NFPS40	4.0	4.0	6.0	20.0	36.0	120°	N/A	1



60PFS40

Bivona Cuffless Pediatric FlexTend Tracheostomy Tubes ("V" Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFS25	2.5	2.5	4.0	20.0	38.0	120°	N/A	1
60PFS30	3.0	3.0	4.7	20.0	39.0	120°	N/A	1
60PFS35	3.5	3.5	5.3	20.0	40.0	120°	N/A	1
60PFS40	4.0	4.0	6.0	20.0	41.0	120°	N/A	1
60PFS45	4.5	4.5	6.7	30.0	42.0	120°	N/A	1
60PFS50	5.0	5.0	7.3	30.0	44.0	120°	N/A	1
60PFS55	5.5	5.5	8.0	30.0	46.0	120°	N/A	1



60PFSS25

Bivona Cuffless FlexTend Pediatric Tracheostomy Tubes (Straight Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFSS25	2.5	2.5	4.0	20.0	38.0	120°	N/A	1
60PFSS30	3.0	3.0	4.7	20.0	39.0	120°	N/A	1
60PFSS35	3.5	3.5	5.3	20.0	40.0	120°	N/A	1
60PFSS40	4.0	4.0	6.0	20.0	41.0	120°	N/A	1
60PFSS45	4.5	4.5	6.7	30.0	42.0	120°	N/A	1
60PFSS50	5.0	5.0	7.3	30.0	44.0	120°	N/A	1
60PFSS55	5.5	5.5	8.0	30.0	46.0	120°	N/A	1

Bivona™ FlexTend Plus™

Reinforced proximal and distal sections provide a high level of kink and crush resistance. Extended proximal tube allows 15 mm connector to be distal from the chin. FlexTend Plus tubes also provide additional distal length for patients requiring additional overall length.



60PFP45

Bivona Cuffless Pediatric FlexTend Tracheostomy Tubes (Extra Length, "V" Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFP40	4.0	4.0	6.0	30.0	44.0	120°	N/A	1
60PFP45	4.5	4.5	6.7	30.0	48.0	120°	N/A	1
60PFP50	5.0	5.0	7.3	30.0	50.0	120°	N/A	1
60PFP55	5.5	5.5	8.0	30.0	52.0	120°	N/A	1
60PFP60	6.0	6.0	8.7	30.0	56.0	120°	N/A	1



60PFPS50

Bivona Cuffless FlexTend Pediatric Tracheostomy Tubes (Extra Length, Straight Neck Flange)								
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Proximal Length (mm)	Distal Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60PFPS40	4.0	4.0	6.0	30.0	44.0	120°	N/A	1
60PFPS45	4.5	4.5	6.7	30.0	48.0	120°	N/A	1
60PFPS50	5.0	5.0	7.3	30.0	50.0	120°	N/A	1
60PFPS55	5.5	5.5	8.0	30.0	52.0	120°	N/A	1
60PFPS60	6.0	6.0	8.7	30.0	56.0	120°	N/A	1

Bivona™ Aire-Cuf™

A flexible, adjustable tube with a soft, secure cuff. Neonatal and pediatric tubes are reinforced to give a high level of kink and crush resistance. Select adult configurations feature a talk attachment allowing vocalization for ventilator-dependent patients with viable vocal cords.




Hospital



Home



MR
Conditional



750160

Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes

Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
750150	5.0	5.0	7.4	60.0	100°	17.0	1
750160	6.0	6.0	8.8	70.0	100°	18.0	1
750170	7.0	7.0	10.0	80.0	100°	21.0	1
750180	8.0	8.0	11.0	88.0	100°	25.0	1
750190	9.0	9.0	12.3	98.0	100°	28.0	1
750195	9.5	9.5	13.3	98.0	100°	28.0	1



755180

Bivona Mid-Range Aire-Cuf Adult Tracheostomy Tubes with Talk Attachment

Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
755150-SM	5.0	5.0	7.3	60.0	100°	17.0	1
755160-SM	6.0	6.0	8.7	70.0	100°	18.0	1
755170-SM	7.0	7.0	10.0	80.0	100°	21.0	1
755180-SM	8.0	8.0	11.0	88.0	100°	25.0	1
755190-SM	9.0	9.0	12.3	98.0	100°	28.0	1
755195-SM	9.5	9.5	13.3	98.0	100°	28.0	1



65N035

Bivona Aire-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)

Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65N025	2.5	2.5	4.0	30.0	120°	12.0	1
65N030	3.0	3.0	4.7	32.0	120°	14.0	1
65N035	3.5	3.5	5.3	34.0	120°	14.0	1
65N040	4.0	4.0	6.0	36.0	120°	18.0	1

Bivona™ Aire-Cuf™



65SN025

Bivona Aire-Cuf Neonatal Tracheostomy Tubes (Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65SN025	2.5	2.5	4.0	30.0	120°	12.0	1
65SN030	3.0	3.0	4.7	32.0	120°	14.0	1
65SN035	3.5	3.5	5.3	34.0	120°	14.0	1
65SN040	4.0	4.0	6.0	36.0	120°	18.0	1



65P030

Bivona Aire-Cuf Pediatric Tracheostomy Tubes ("V" Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65P025	2.5	2.5	4.0	38.0	120°	12.0	1
65P030	3.0	3.0	4.7	39.0	120°	14.0	1
65P035	3.5	3.5	5.3	40.0	120°	14.0	1
65P040	4.0	4.0	6.0	41.0	120°	18.0	1
65P045	4.5	4.5	6.7	42.0	120°	18.0	1
65P050	5.0	5.0	7.3	44.0	120°	20.0	1
65P055	5.5	5.5	8.0	46.0	120°	20.0	1




65SP030

Bivona Aire-Cuf Pediatric Tracheostomy Tubes (Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
65SP025	2.5	2.5	4.0	38.0	120°	12.0	1
65SP030	3.0	3.0	4.7	39.0	120°	14.0	1
65SP035	3.5	3.5	5.3	40.0	120°	14.0	1
65SP040	4.0	4.0	6.0	41.0	120°	18.0	1
65SP045	4.5	4.5	6.7	42.0	120°	18.0	1
65SP050	5.0	5.0	7.3	44.0	120°	20.0	1
65SP055	5.5	5.5	8.0	46.0	120°	20.0	1

Bivona™ Fome-Cuf™

The self-expanding cuff helps reduce the risk of overinflation or underinflation associated with traditional cuffs. Select adult configurations feature a talk attachment allowing vocalization for ventilator-dependent patients with viable vocal cords.



850180

Bivona Fome-Cuf Adult Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
850150	5.0	5.0	7.3	60.0	100°	20.0	1
850160	6.0	6.0	8.7	70.0	100°	30.0	1
850170	7.0	7.0	10.0	80.0	100°	32.0	1
850180	8.0	8.0	11.0	88.0	100°	33.0	1
850190	9.0	9.0	12.3	98.0	100°	33.0	1
850195	9.5	9.5	13.3	98.0	100°	33.0	1



855180

Bivona Fome-Cuf Adult Tracheostomy Tubes (with Talk Attachment)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
855150	5.0	5.0	7.3	60.0	100°	20.0	1
855160	6.0	6.0	8.7	70.0	100°	30.0	1
855170	7.0	7.0	10.0	80.0	100°	32.0	1
855180	8.0	8.0	11.0	88.0	100°	33.0	1
855190	9.0	9.0	12.3	98.0	100°	33.0	1



85N030

Bivona Fome-Cuf Neonatal Tracheostomy Tubes ("V" Neck Flange)*							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
85N030	3.0	3.0	4.7	32.0	N/A	14.0	1
85N035	3.5	3.5	5.3	34.0	N/A	14.0	1
85N040	4.0	4.0	6.0	36.0	N/A	18.0	1

* Products listed above are MR Conditional.



85P035

Bivona Fome-Cuf Pediatric Tracheostomy Tubes ("V" Neck Flange)*							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
85P025	2.5	2.5	4.0	38.0	N/A	12.0	1
85P030	3.0	3.0	4.7	39.0	N/A	14.0	1
85P035	3.5	3.5	5.3	40.0	N/A	14.0	1
85P040	4.0	4.0	6.0	41.0	N/A	18.0	1
85P045	4.5	4.5	6.7	42.0	N/A	18.0	1
85P050	5.0	5.0	7.3	44.0	N/A	20.0	1
85P055	5.5	5.5	8.0	46.0	N/A	20.0	1

* Products listed above are MR Conditional.

Bivona™ Adjustable Hyperflex™


Adjustable flange tubes allow for immediate modification of flange position to meet a variety of patient needs. It can be adjusted for horizontal and vertical shaft drop to accommodate unusual anatomy or pathology. Flexible wire-reinforced silicone shafts provide improved kink and crush resistance. The adjustable flange tubes are single use and cannot be reprocessed.



Hospital



MR
Conditional



67HA60

Bivona TTS Adjustable Neck Flange Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67HA60	6.0	6.0	9.2	110.0	180°	11.0	1
67HA70	7.0	7.0	10.6	120.0	180°	12.0	1
67HA80	8.0	8.0	11.7	130.0	180°	13.0	1
67HA90	9.0	9.0	12.9	140.0	180°	14.5	1



75HA70

Bivona Mid-Range Aire-Cuf™ Adjustable Neck Flange Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
75HA60	6.0	6.0	9.2	110.0	180°	18.0	1
75HA70	7.0	7.0	10.6	120.0	180°	21.0	1
75HA80	8.0	8.0	11.7	130.0	180°	25.0	1
75HA90	9.0	9.0	12.9	140.0	180°	28.0	1



60HA35

Bivona Cuffless Tube Adjustable Neck Flange Hyperflex Pediatric “V” Neck Flange Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Maximum Usable Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60HA25	2.5	2.5	4.0	55.0	180°	N/A	1
60HA30	3.0	3.0	4.7	60.0	180°	N/A	1
60HA35	3.5	3.5	5.3	65.0	180°	N/A	1
60HA40	4.0	4.0	6.0	70.0	180°	N/A	1
60HA45	4.5	4.5	6.7	75.0	180°	N/A	1
60HA50	5.0	5.0	7.3	80.0	180°	N/A	1
60HA55	5.5	5.5	8.0	85.0	180°	N/A	1

Bivona™ Fixed Flange Hyperflex™

The fixed flange tube provides a more permanent solution for patients who require extra vertical and horizontal length. Flexible wire-reinforced silicone shafts provide improved kink and crush resistance.



Hospital



Home



MR
Conditional



60AFHXL70

Bivona Cuffless Extra Length Fixed Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60AFHXL60	6.0	6.0	9.2	100.0	180°	N/A	1
60AFHXL70	7.0	7.0	10.6	110.0	180°	N/A	1
60AFHXL80	8.0	8.0	11.7	120.0	180°	N/A	1
60AFHXL90	9.0	9.0	12.9	130.0	180°	N/A	1



67FHXL70

Bivona TTS™ Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
67FHXL60	6.0	6.0	9.2	100.0	180°	11.0	1
67FHXL70	7.0	7.0	10.6	110.0	180°	12.0	1
67FHXL80	8.0	8.0	11.7	120.0	180°	13.0	1
67FHXL90	9.0	9.0	12.9	130.0	180°	14.5	1



75FHXL60

Bivona Cuffed Extra Length Fixed Neck Flange Hyperflex Tracheostomy Tubes							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
75FHXL60	6.0	6.0	9.2	100.0	180°	18.0	1
75FHXL70	7.0	7.0	10.6	110.0	180°	21.0	1
75FHXL80	8.0	8.0	11.7	120.0	180°	25.0	1
75FHXL90	9.0	9.0	12.9	130.0	180°	28.0	1

Bivona™ Cuffless

When a cuffless tube is indicated, the all-silicone construction of this soft and flexible tube ensures comfort for the patient.



Hospital



Home



MR
Conditional



60A160

Bivona Cuffless Adult Tracheostomy Tubes*							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60A150	5.0	5.0	7.4	60.0	N/A	N/A	1
60A160	6.0	6.0	8.8	70.0	N/A	N/A	1
60A170	7.0	7.0	10.0	80.0	N/A	N/A	1
60A180	8.0	8.0	11.0	88.0	N/A	N/A	1
60A190	9.0	9.0	12.3	98.0	N/A	N/A	1

* Products listed above are MR Safe.



60N030

Bivona Cuffless Neonatal Tracheostomy Tubes ("V" Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60N025	2.5	2.5	4.0	30.0	120°	N/A	1
60N030	3.0	3.0	4.7	32.0	120°	N/A	1
60N035	3.5	3.5	5.3	34.0	120°	N/A	1
60N040	4.0	4.0	6.0	36.0	120°	N/A	1



60SN025

Bivona Cuffless Neonatal Tracheostomy Tubes (Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60SN025	2.5	2.5	4.0	30.0	120°	N/A	1
60SN030	3.0	3.0	4.7	32.0	120°	N/A	1
60SN035	3.5	3.5	5.3	34.0	120°	N/A	1
60SN040	4.0	4.0	6.0	36.0	120°	N/A	1



60P040

Bivona Cuffless Pediatric Tracheostomy Tubes ("V" Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60P025	2.5	2.5	4.0	38.0	120°	N/A	1
60P030	3.0	3.0	4.7	39.0	120°	N/A	1
60P035	3.5	3.5	5.3	40.0	120°	N/A	1
60P040	4.0	4.0	6.0	41.0	120°	N/A	1
60P045	4.5	4.5	6.7	42.0	120°	N/A	1
60P050	5.0	5.0	7.3	44.0	120°	N/A	1
60P055	5.5	5.5	8.0	46.0	120°	N/A	1

Bivona™ Cuffless



Bivona Cuffless Pediatric Tracheostomy Tubes (Straight Neck Flange)							
Product Code	Tube Size (mm)	Tube ID (mm)	Tube OD (mm)	Tube Length (mm)	Tube Angle	Cuff Resting Diameter (mm)	Units Per Case
60SP025	2.5	2.5	4.0	38.0	120°	N/A	1
60SP030	3.0	3.0	4.7	39.0	120°	N/A	1
60SP035	3.5	3.5	5.3	40.0	120°	N/A	1
60SP040	4.0	4.0	6.0	41.0	120°	N/A	1
60SP045	4.5	4.5	6.7	42.0	120°	N/A	1
60SP050	5.0	5.0	7.3	44.0	120°	N/A	1
60SP055	5.5	5.5	8.0	46.0	120°	N/A	1

Bivona™ Customized Silicone Tracheostomy Tubes

Treating a patient with difficult or unusual anatomy can be a challenge; getting that patient a proper-fitting tracheostomy tube shouldn't be. Bivona customization lab technicians can handcraft a tube of soft, durable silicone to conform to the patient's neck, stoma, and trachea.

MyBivona™ Customized Tracheostomy Tube Configurator

With easy-to-use, intuitive menus and powerful, real-time 3D rendering capabilities that let you visualize every detail, it's never been easier to create a uniquely comfortable fit for your patients.



Start customizing today. Scan the QR code or visit <http://www.icumed.com/MyBivona>




Customized Tracheostomy Tubes Online Templates

Designing a Bivona customized tracheostomy tube can be completed using our online templates.

Bivona customized tracheostomy tube delivery times

- › Standard service (sterile): 3-4 weeks
- › Standard service (non-sterile): 7 business days
- › Express service (non-sterile): 3 business days
- › Next-day service (non-sterile): Shipped same day if order is received by 10 a.m. CT.
- › Next-day service (non-sterile): Shipped same day if order is received by 12 p.m. CT (additional fees apply).



Accessories

Accessories		
Product Code	Description	Units Per Case
CON15S	Side Port Connectors NS	50
REDCAP	15 mm Red Cap Plugs	26
DWG300	Disconnect Wedge 15 mm	25
002873	Portex™ HEPA Filter, Straight, 200-1500 mL	50
002815	HME Neonatal, 10 mL	50
525351	Dual-Axis Swivel Adapter (Sterile)	50
570016	Thermovent™ T (Sterile)	50
570017	Thermovent T (Non-Sterile)	50
580011	Thermovent 600 HME	50
100/570/022	Thermovent T2 (Sterile)	50
100/575/010	Oxygen Delivery for HME Thermovent T Thermovent O ₂	10
100/582/000	Thermovent 1200 HME, Adult, Straight	20
625191	Fiber-Optic Bronchoscope Swivel Adapter (Sterile)	25

Visit bivonacustom.com for more information or contact us:

Phone: 1-800-424-8662

Email: bivonacustomorder@icumed.com

