


# Gripper Plus Power P.A.C.<sup>™</sup>

## Port Access Needle



### Safety Huber Needle

- › A safety port access needle to administer or withdraw fluids from implanted ports
- › Available in different sizes and lengths
- › Safety mechanism to be activated when removing the needle — audible click and visible engagement confirm that the safety arm is activated into lock position
- › Thick cushioned needle platform sits against the skin when in-situ
- › Non-siliconized, non-coring needle
- › The 20 cm-long extension set is available with or without Y-site with luer-activated needleless injection site
- › Extension set with blue tubing differentiates from non-power needles
- › Color-coded clamp(s)
- › Priming volume 0.3 mL or 0.8 mL depending on the extension set with or without the Y-site version
- › Power injection up to 5 mL/sec and pressure up to 300 psi
- › Positive ANTT evaluation

### Ordering Information

#### Needleless Y-site — luer-activated needleless injection site

Product Code	Gauge	Length	Units
21-3464-24		5/8 in	
21-3465-24	● 20	3/4 in	12/box
21-3466-24		1 in	
21-3467-24		1 1/4 in	
21-3468-24		3/4 in	
21-3469-24	● 19	1 in	12/box
21-3470-24		1 1/4 in	

#### Without Y-site

Product Code	Gauge	Length	Units
21-3369-24		5/8 in	
21-3367-24	● 20	3/4 in	12/box
21-3362-24		1 in	
21-3363-24		1 1/4 in	
21-3368-24		3/4 in	
21-3364-24	● 19	1 in	12/box
21-3365-24		1 1/4 in	

