



Achieve Efficient Port Access and Enhanced Safety Protocols

Introducing the ICU Medical portfolio of port access needles.




Port-A-Cath™ Plastic Hub

Standard non-safety Huber Needle. Available in three different needle gauge sizes in a straight or bent 90° hub version. MRI conditional up to 3-Tesla.




Gripper™

Non-safety Huber Needle with a small footprint base, and a thick cushioned needle platform to support patient comfort. Removable contoured grip for easy insertion. Available in different needle sizes and lengths. With or without Y-site.




Gripper Plus™

Safety Huber Needle with a thick cushioned needle platform to support patient comfort. Removable contoured grip for easy insertion. With a safety mechanism that activates with an audible click and visible engagement confirmation during needle removal. Gripper Plus is available in different needle sizes and lengths. With or without Y-site.



Gripper Plus™ Power

Gripper Plus Power delivers the features and benefits of Gripper Plus with the addition of a power-rated safety Huber Needle and clear blue tubing for easy differentiation from non-power Huber needles.



Gripper™ Micro

Gripper Micro helps support safety protocols with an audible click, visible engagement confirmation, and a Trocar with Huber cut that leaves the blunt cannula in the port for safe removal and disposal. A tiny footprint and thick cushioned needle platform to support patient comfort. Available in different needle sizes and length. With or without Y-site.

icumedical
human connections