

Gripper Plus™ Power P.A.C.

Port Access Needle

Safety Huber Needle

- › A safety port access needle to administer or withdraw fluids from implanted ports
- › Available in different sizes and lengths
- › Safety mechanism to be activated when removing the needle — audible click and visible engagement confirm that the safety arm is activated into lock position
- › Thick cushioned needle platform sits against the skin when in-situ
- › Non-siliconised, non-coring needle
- › The 20 cm-long extension set is available with or without Y-site with luer-activated needleless injection site
- › Extension set with blue tubing differentiates from non-power needles
- › Colour-coded clamp(s)
- › Priming volume 0.3 mL or 0.8 mL depending on the extension set with or without the Y-site version
- › 19G and 20G power injection up to 5 mL/sec and pressure up to 300 psi
- › 22G power injection up to 2 mL/sec and pressure up to 300 psi
- › Positive ANTT® evaluation

Ordering Information

Needleless Y-site — luer-activated needleless injection site			
Product Code	Gauge	Length	Units
21-3460-24	● 22	5/8 in	12/box
21-3461-24		3/4 in	
21-3462-24		1 in	
21-3463-24		1 1/4 in	
21-3464-24	● 20	5/8 in	12/box
21-3465-24		3/4 in	
21-3466-24		1 in	
21-3467-24		1 1/4 in	
21-3468-24	● 19	3/4 in	12/box
21-3469-24		1 in	
21-3470-24		1 1/4 in	

Without Y-site			
Product Code	Gauge	Length	Units
21-3370-24	● 22	5/8 in	12/box
21-3366-24		3/4 in	
21-3372-24		1 in	
21-3373-24		1 1/4 in	
21-3369-24	● 20	5/8 in	12/box
21-3367-24		3/4 in	
21-3362-24		1 in	
21-3363-24		1 1/4 in	
21-3368-24	● 19	3/4 in	12/box
21-3364-24		1 in	
21-3365-24		1 1/4 in	

Access with Confidence.
De-access with Control.



ANTT® (Aseptic Non Touch Technique) is a trademark of The Association for Safe Aseptic Practice (ASAP)