

Saf-T Wing™ Blood Collection Sets

Promotes end-to-end safety for end-to-end compliance


The Saf-T Wing blood collection set is designed with many innovative features that help promote safety and ease of use.

In-vein activation

Safety mechanism is activated while the needle is in the patient's vein—helping limit the phlebotomist's exposure to a contaminated needle



Advanced features designed to help protect you and your patients



Ordering Information



Saf-T Wing™ Blood Collection and Infusion Sets				
Product Code	Description	Color	Units per Box	Units per Case
972106	21G x 3/4" with 6" Tubing	● Green	50	200
972112	21G x 3/4" with 12" Tubing	● Green	50	200
972306	23G x 3/4" with 6" Tubing	● Blue	50	200
972312	23G x 3/4" with 12" Tubing	● Blue	50	200
972506	25G x 3/4" with 6" Tubing	● Orange	50	200
972512	25G x 3/4" with 12" Tubing	● Orange	50	200



Saf-T Wing™ Blood Collection Sets and Detached Holder				
Product Code	Description	Color	Units per Box	Units per Case
982106D	21G x 3/4" Saf-T Wing with 6" Tubing and Detached Holder	● Green	50	200
982112D	21G x 3/4" Saf-T Wing with 12" Tubing and Detached Holder	● Green	50	200
982306D	23G x 3/4" Saf-T Wing with 6" Tubing and Detached Holder	● Blue	50	200
982312D	23G x 3/4" Saf-T Wing with 12" Tubing and Detached Holder	● Blue	50	200
982506D	25G x 3/4" Saf-T Wing with 6" Tubing and Detached Holder	● Orange	50	200
982512D	25G x 3/4" Saf-T Wing with 12" Tubing and Detached Holder	● Orange	50	200

All Saf-T Holder devices incorporate a multi-sample Luer adapter.



Saf-T Holder™ Devices and Blood Transfer Devices				
Product Code	Description		Units per Box	Units per Case
96000	Saf-T Holder with Multi-Sample Male Luer Adapter		50	200
96000S	Saf-T Holder with Multi-Sample Male Luer Adapter		50	200
96002	Saf-T Holder with Female Luer Adapter		50	200
96002S	Saf-T Holder with Female Luer Adapter		50	200
96004	Saf-T Holder Blood Culture Device with Male Luer Adapter		–	50
96005	Saf-T Holder Blood Culture Device with Female Luer Adapter		–	50
96007	Pediatric Tube Insert		10	100

All Saf-T Holder devices incorporate a multi-sample Luer adapter.