



Rio™

Drug Reconstitution Transfer Device

Facilitates preparation of IV medication admixtures
when and where you need them

icumedical
human connections

Meeting your admixture needs in the pharmacy or at the point of care



- › For reconstitution of liquid or lyophilized 20 mm standard drug vials and aseptic transfer into the multi-port LifeCare™ IV container systems (up to 250 mL)*
- › Improve operational flexibility with immediate or future activation
- › After reconstitution, the stopcock allows the fluid pathway to be closed, preventing reconstituted drug from re-entering the vial
- › Needlefree design helps reduce needlestick related injuries and sharps disposal costs



Latex-free | Non-PVC/Non-DEHP
Single-use only | Individual sterile packaging

Product Information	Item #	Description	Case Qty
	VB-20	Rio - Drug Reconstitution Transfer Device	50



In-Service Video

Your guide to the safe and effective use of the Rio Drug Reconstitution Transfer Device



To learn how to put the Rio drug reconstitution transfer device to work for you, please call 866.829.9025 or visit www.icumed.com

*Rio Drug Reconstitution Transfer Device is manufactured by Otsuka ICU Medical LLC and distributed by ICU Medical Sales, Inc.

*LifeCare IV container systems are manufactured by Otsuka ICU Medical LLC and distributed by ICU Medical Sales, Inc.

© 2025 ICU Medical, Inc. | P25-4873-R1 | **Rx Only** For safe and proper use, refer to the Instructions for Use.

icumedical
human connections

800.824.7890 | www.icumed.com