




Gripper Plus™ Power P.A.C.

19G, 20G, and 22G Safety Port Access Needle

Access with Confidence. De-access with Control.

icumedical
human connections



Safety You Can Feel. Performance They Can Trust.

Every port access carries a moment of unease: Will the needle slip? Is the site sealed? Too much discomfort can turn an already anxious patient into a reluctant one.

By pairing a built-in safety mechanism with a comfort-first design, Gripper™ Plus Power Port Access Needles are designed to help you deliver care with a calm confidence your patients can feel.



Engineered to Deliver Comfort and Confidence

Accidental needlesticks and improper aseptic technique can lead to increased risk of infection, anxiety for staff, and discomfort for patients. Gripper Plus Power P.A.C. is designed to address these challenges every step of the way.

Our safety mechanism activates upon removal with a satisfying click—an audible and visible lock that means you can trust the needle is secured. Meanwhile, a thick cushioned platform gently sits against the skin during infusion, providing both stability and comfort.

Removable Grip Tab

For a comfortable grip during insertion.

Built-In Safety Mechanism

A one-step safety arm activation with a clear click and visible lock to help protect staff and patients.

Versatile Configurations

- › Available in 19G, 20G, and 22G configurations.
- › 20 cm-long extension set with or without Y-site, plus a colour-coded clamp to differentiate from non-power needles.
- › Prime volume of 0.3 mL or 0.8 mL (depending on Y-site) meets a variety of clinical needs.

Non-Coring Design

Sharp, non-siliconised needle prevents coring and preserves the integrity of the port while helping reduce needle retention—because every procedure should be as smooth as possible.

Easy to Identify

Blue tubing helps clearly identify lines that are rated for power injection.

Power Injection Ready up to 300 psi

Engineered for high-pressure injection up to 5 mL/sec, giving you the flexibility to handle contrast studies when needed. 22G needle up to 2 mL/sec.

Comfort Where It Counts

A thick needle cushion rests gently on the skin, helping to reduce skin irritation and maximise patient comfort.

Positive ANTT Evaluation

Recognised compatibility with ANTT guidelines, supporting your efforts to maintain a higher standard of infection control and improve your aseptic technique.



ANTT is a registered trademark of The Association of Safe Aseptic Practice (The-ASAP)
Product availability varies by region. Please contact your local representative for availability in your area.

Take the Next Step in Port Access Safety

*Built to make your procedures safer, simpler, and more efficient.
The difference is in the details—and every detail has been
meticulously crafted to serve you better.*

When used with power injectable implantable venous access ports,
this product is indicated for power injecting contrast media.

With power rating in all sizes, choose Gripper Plus Power P.A.C. for a safe
port access experience—from the very first access click to the very last drop.
Your peace of mind has never been so well-designed.*

Learn more at <http://www.icumed.com.au>

*When used with power injectable implantable venous access ports,
this product is indicated for power injecting contrast media.

icumedical
human connections