




SurgiVet™ VitaLink™ Monitor

Replacement part accessories and optional accessories



Item Code	Description
Monitor	
VM50	SurgiVet VitaLink 15" Monitor, 3 Lead ECG, SpO ₂ , RESP, 2-CH TEMP, NIBP, Printer, Std Battery, Wi-fi, 1/ea
ECG	
VA0724	3 Lead Wire ECG Cable, Lead Wire/Separated, US Type, 1/ea
VA0728	3 Lead ECG Trunk Cable, Separated, US Type, 1/ea
NIBP	
VA0911	3-6 cm Suntech NIBP Cuff, 1/ea
VA0912	4-8 cm Suntech NIBP Cuff, 1/ea
VA0913	6-11 cm Suntech NIBP Cuff, 1/ea
VA0914	7-13 cm Suntech NIBP Cuff, 1/ea
VA0915	8-15 cm Suntech NIBP Cuff, 1/ea
VA0916	17-25 cm, Suntech NIBP Cuff, 1/ea
VA0917	10 ft Suntech NIBP Hose, 1/ea
Other	
VA1057	US Power Cord, 1/ea
VA7813	Operator's Manual
VM6082	1hr Standard Rechargeable Lithium Ion Battery, 1/ea
SPO₂	
VA0266	SpO ₂ Extension Cable, 1/ea
VA0438	SpO ₂ Sensor (WA-103 SpO ₂ Y-SENSOR), 1/ea
VA0452	SpO ₂ Reflectance Sensor, 1/ea
Temperature	
VA0146	Reuseable Temperature Probe, 1/ea
Dual-Gas (Optional Accessory)	
VMPM-NH	SurgiVet VitaLink Dual-Gas Analyzer Sensor Module and Water Trap, 1/ea
VMPM-IH	SurgiVet VitaLink Dual-Gas Interface Module, 1/ea
VA0741	SurgiVet VitaLink Water Trap for ITEC Gas Analyzer, 1/ea
VA0742	SurgiVet VitaLink 8 ft Patient Luer Lock Sample Line for GAS Analyzer, 1/ea
VA0743	SurgiVet VitaLink L-angle Sample Line Adapter for GAS Analyzer, 1/ea
VA0744	SurgiVet VitaLink Straight Sample Line Adapter for GAS Analyzer, 1/ea
ECG (Optional Accessory)	
VA0734	5 Lead ECG Cable, All In One, Grab, US Type, 1/ea
VA0749	SpO ₂ Veterinary SpO ₂ Sensor, 1/ea

Item Code	Description
EtCO₂ Kingst (Optional Accessory)	
VA0703	EtCO ₂ Sidestream Sampling Line, Kingst, 1/ea
VA0704	EtCO ₂ Sidestream Airway Adapter, Kingst, 1/ea
VA0701	EtCO ₂ Sidestream Water Filter, Kingst, 1/ea
VA0706	EtCO ₂ Adult Mainstream Airway Adapter, Kingst, 1/ea
VA0702	EtCO ₂ Sidestream Module Receptacle, Kingst, 1/ea
VA0707	EtCO ₂ Infant Mainstream Airway Adapter, Kingst, 1/ea
VM0083	EtCO ₂ Mainstream Sensor, Kingst, 1/ea
VM0084	EtCO ₂ Sidestream Sensor, Kingst, 1/ea
EtCO₂ Respirationics & Kingst (Optional Accessory)	
VMPM-ICE	SurgiVet VitaLink EtCO ₂ Respirationics and Kingst Interface Module, 1/ea
EtCO₂ Respirationics (Optional Accessory)	
VA0204	Medium/Large Loflo Cannula, 1/ea
VA0229	Module Mounting Bracket, 1/ea
VA0221	Large Airway Adapter, Single Use, 1/ea
VA0222	Extra Small Airway Adapter, Single Use, 1/ea
VA0219	Nafion Straight Sampling Line, 1/ea
VM0019	EtCO ₂ CAPNOSTAT 5 Mainstream Sensor, Respirationics, 1/ea
VM0018	EtCO ₂ Loflo Side Stream Sensor, Respirationics, 1/ea
IBP (Optional Accessory)	
VMPM-TIE	SurgiVet VitaLink IBP Interface Module, 1/ea
VA0196	SurgiVet VitaLink IBP Cable, 1/ea
VA0796	SurgiVet VitaLink IBP Disposable Transducer, 1/ea
Other (Optional Accessory)	
VA0774	ECG Alligator Clip, 1/ea
VA0458	AU Power Cord, 1/ea
VM6083	4hr Premium Rechargeable Lithium Ion Battery, 1/ea
VX6232	47 mm Thermal Printer Paper, 10/cs
VA0935	EU Power Cord, 1/ea
VA0802	SurgiVet VitaLink / VitaLink Mini Monitor Roll Stand, 1/ea
SPO₂ (Optional Accessory)	
VA0771	3 Lead Snap Wire ECG Cable, Lead Snap Wire/ Separated, US Type, 1/ea



SurgiVet™ VitaLink™ Mini Monitor

Replacement part accessories and optional accessories



Item Code	Description
Monitor	
VM40	SurgiVet™ VitaLink™ Mini 10.4" Monitor, 3 Lead ECG, SpO ₂ , RESP, 2-CH TEMP, NIBP, Printer, Std Battery, Wi-fi, 1/ea
ECG	
VA0724	3 Lead Wire ECG Cable, Lead Wire/Separated, US Type, 1/ea
VA0728	3 Lead ECG Trunk Cable, Separated, US Type, 1/ea
NIBP	
VA0911	3-6 cm Suntech NIBP Cuff, 1/ea
VA0912	4-8 cm Suntech NIBP Cuff, 1/ea
VA0913	6-11 cm Suntech NIBP Cuff, 1/ea
VA0914	7-13 cm Suntech NIBP Cuff, 1/ea
VA0915	8-15 cm Suntech NIBP Cuff, 1/ea
VA0916	17-25 cm, Suntech NIBP Cuff, 1/ea
VA0917	10 ft Suntech NIBP Hose, 1/ea
Other	
VA1057	US Power Cord, 1/ea
VA7813	Operator's Manual
VM6082	1hr Standard Rechargeable Lithium Ion Battery, 1/ea
SPO₂	
VA0266	SpO ₂ Extension Cable, 1/ea
VA0438	SpO ₂ Sensor (WA-103 SpO ₂ Y-SENSOR), 1/ea
VA0452	SpO ₂ Reflectance Sensor, 1/ea
Temperature	
VA0146	Reuseable Temperature Probe, 1/ea
ECG (Optional Accessory)	
VA0734	5 Lead ECG Cable, All In One, Grap, US Type, 1/ea
VA0749	SpO ₂ Veterinary SpO ₂ Sensor, 1/ea

Item Code	Description
EtCO₂ Kingst (Optional Accessory)	
VA0703	EtCO ₂ Sidestream Sampling Line, Kingst, 1/ea
VA0704	EtCO ₂ Sidestream Airway Adapter, Kingst, 1/ea
VA0701	EtCO ₂ Sidestream Water Filter, Kingst, 1/ea
VA0706	EtCO ₂ Adult Mainstream Airway Adapter, Kingst, 1/ea
VA0702	EtCO ₂ Sidestream Module Receptacle, Kingst, 1/ea
VA0707	EtCO ₂ Infant Mainstream Airway Adapter, Kingst, 1/ea
VM0083	EtCO ₂ Mainstream Sensor, Kingst, 1/ea
VM0084	EtCO ₂ Sidestream Sensor, Kingst, 1/ea
EtCO₂ Respirationics (Optional Accessory)	
VA0204	Medium/Large Loflo Cannula, 1/ea
VA0229	Module Mounting Bracket, 1/ea
VA0221	Large Airway Adapter, Single Use, 1/ea
VA0222	Extra Small Airway Adapter, Single Use, 1/ea
VA0219	Nafion Straight Sampling Line, 1/ea
VM0019	EtCO ₂ CAPNOSTAT 5 Mainstream Sensor, Respirationics, 1/ea
VM0018	EtCO ₂ Loflo Side Stream Sensor, Respirationics, 1/ea
Other (Optional Accessory)	
VA0774	ECG Alligator Clip, 1/ea
VA0458	AU Power Cord, 1/ea
VM6083	4hr Premium Rechargeable Lithium Ion Battery, 1/ea
VX6232	47 mm Thermal Printer Paper, 10/cs
VA0935	EU Power Cord, 1/ea
VA0802	SurgiVet™ VitaLink™ / VitaLink™ Mini Monitor Roll Stand, 1/ea
SPO₂ (Optional Accessory)	
VA0771	3 Lead Snap Wire ECG Cable, Lead Snap Wire/ Separated, US Type, 1/ea

Vet technical support hotline 888.745.6562

All trademarks are the property of their respective owners.
 FDA does not require submission of a 510(k) or any premarket approval for devices used in veterinary medicine.
 CAUTION: Federal law restricts this device to use on or on the order of a licensed veterinarian

icumedical
 human connections