



LifeShield™

Infusion Safety Software

Bring Teams Together for Safer,
Streamlined and Synchronized Care



Unite medication management in a single platform, removing barriers and streamlining collaboration to amplify the impact of every role.

For Pharmacists

Comprehensive Medication Management

- › Update medications, rulesets and drug libraries quickly and easily across multiple facilities with streamlined workflows.
- › Access near real-time infusion data to uncover insights, drive safety improvements and ensure protocols are always precise.
- › Simplify and streamline drug library review with a change summary report that clearly shows modifications.
- › The cloud environment allows for rapid, frequent and hands-free drug library updates.

For Clinicians

Simplified Safety at the Bedside

- › Reduce errors and spend less time programming with pre-populated doses, rates and durations.
- › Get DERS guidance with depth, like being alerted when a medication is increased or decreased by a certain percentage—even within soft and hard limits.
- › Enable bedside features that enhance infusion performance including infuse to empty, auto flush, syringe characterization and the ability to enable and deliver bolus doses even throughout a running infusion.



For Biomedical Teams

Seamless Fleet Management

- › See your entire fleet on a single screen with a centralized platform to manage, troubleshoot and optimize your entire device fleet.
- › Push drug library and software updates wirelessly, keeping every device ready without disrupting patient care.
- › Stay up to date on PVT testing results, fleet battery health, connectivity status and more with biomed reports.
- › Keep track of devices with real-time location systems (RTLS) or wireless access point naming.

For IT Teams

Control and Security Made Simple

- › Manage security and user accounts in one place, saving time and reducing the burden on your team.
- › Provide an additional layer of security with authentication verification.
- › Protect your system with industry-standard cybersecurity, including FDA clearance, HITRUST certification, independent third-party testing and ongoing vigilance.

All-inclusive software package, maintenance and support

With LifeShield software, there are no hidden fees or surprises. Everything is included—comprehensive cloud-based software, clinical training and implementation support—without extra hardware costs or implementation fees.