



# Important

Please open this form in Adobe Acrobat.  
It has interactive elements which may not function correctly  
if opened in your web browser.

If you don't have Adobe Acrobat Reader installed, you can download it here:

<https://get.adobe.com/reader/>

# Bivona™ Customized FlexTend™ Tracheostomy Tube Template

## 1. Shaft Length

1a. Shaft Length \_\_\_\_\_ mm (5 mm min, 60 mm max)




1b. Shaft Length \_\_\_\_\_ mm (30 mm min, 140 mm max)

## 2. Shaft

I.D.	O.D.								
<input type="radio"/> 2.5	4.0	<input type="radio"/> 4.0	6.0	<input type="radio"/> 5.5	8.0	<input type="radio"/> 7.0	10.0	<input type="radio"/> 8.5	11.8
<input type="radio"/> 3.0	4.7	<input type="radio"/> 4.5	6.7	<input type="radio"/> 6.0	8.7	<input type="radio"/> 7.5	10.4	<input type="radio"/> 9.0	12.3
<input type="radio"/> 3.5	5.3	<input type="radio"/> 5.0	7.3	<input type="radio"/> 6.5	9.4	<input type="radio"/> 8.0	11.0	<input type="radio"/> 9.5	13.3

## 3. Connector (check one box only)

3a. Swivel 2.5–9.5 mm  3b. Fixed 2.5–6.0 mm



## 4. Neck Flange (check one box only)






4a. 2.5–6.0 mm

4c. 2.5–9.5 mm

4d. 6.5–9.5 mm

4e. 2.5–4.0 mm

4f. 2.5–6.0 mm



Neck Flange Length: 60 mm

Neck Flange Length: 83 mm

Neck Flange Length: 67 mm

Neck Flange Length: 47 mm

Neck Flange Length: 60 mm

## 5. Cuff (check one box only)






5a. TTS Tight to Shaft

5b. Aire-Cuf™ Pedi

5b. Aire-Cuf Adult

5c. Fome-Cuf™

5d. Cuffless



2.5–9.5 mm

2.5–5.5 mm


5.0–9.5 mm

2.5–9.5 mm


2.5–9.5 mm

## 6. Shaft Option

6a. Shaft (Only proximal section is wire reinforced.)



6b. Hyperflex



2.5–9.5 mm

## 7. Optional

7a. Proximal (Top) Cuff

TTS Cuff

Aire-Cuf

Fome-Cuf

Pediatric

Adult

Cuff position from

Neck Flange \_\_\_\_\_ mm

Note: Minimum distance 2 mm

7b. Stoma Seal

Solid Silicone

\_\_\_\_\_ mm O.D.

\_\_\_\_\_ mm Length

Tapered

Non Tapered

7c. Distal (Bottom) Cuff

Cuff distance from tip to shaft \_\_\_\_\_ mm

Note: Cuff distance from tip cannot be less than the minimum standard distance listed below:

a. Pediatric TTS, Aire-Cuf & Fome-Cuf – 3 mm

b. Adult Aire-Cuf & Fome-Cuf – 5 mm

c. Adult TTS – 9 mm

7d. Talk Attachment

5.0 mm I.D. Adult shaft or above only

## 8. Approval Please note: ICU Medical will review the completed template for compliance.

Any required modifications will be made with the approval of a clinician.

Hospital/Clinic \_\_\_\_\_ Prescribing Physician \_\_\_\_\_

Physician Title \_\_\_\_\_ Physician License or NPI (USA Only) \_\_\_\_\_

P.O. Number \_\_\_\_\_ Date \_\_\_\_\_

## 9. Select Type of Service

Extended next day service (non-sterile), same day shipping if order received by 12 noon (CST) USA only

Standard service (sterile)

Standard service (non-sterile)

Express service (non-sterile)

Next day service (non-sterile) USA only

**icumedical**

human connections

Phone: 1-800-424-8662

Fax: 219-989-7435

Email: bco@icumed.com

Web: www.icumed.com