



Pediatrics and NICU

Product Catalog

Specialty infusion technology designed for your most delicate patients

icumedical
human connections

Help Keep Your Tiniest Patients Safe from Infection

When it comes to caring for your most precious patients, every detail matters. That's why we have developed a full line of IV therapy products designed to meet the special needs of your most delicate patients.



Clave™ Connector Technology

Help minimize infection risks with manifolds and stopcocks that incorporate clinically preferred Clave connector technology



Vascular Access

Empower your delivery of care to the littlest of patients with a specialized portfolio of vascular access tools.



Airway Management

Access a broad portfolio of tracheostomy, respiratory, breathing circuits, ventilation, and specialty airway products.

Warning: Clave connectors may be incompatible with some male-luer connectors including prefilled glass syringes. To avoid damage to the Clave or syringes or male luers which may result in delays of medication administration and possible serious adverse events, users should confirm mating luers or syringes have an internal diameter range of 0.062" to 0.110". Check the internal diameter of the male-luer connector of the mating syringe prior to using it to access the Clave. Products outside of these dimensional tolerances should not be used.

Table of Contents

Needlefree Connectors5

Disinfection Caps.....6

Single Lumen Extension Sets7

Single Lumen Extension Sets – Filtered.....7

Single Lumen Extension Sets – T-Connector8

Single Lumen Extension Sets – High Pressure9

Multi-Lumen Extension Sets11

Medication Sets12

Medication Sets – Filtered14

Closed Medication Sets16

Stopcocks and Manifolds18

Stopcock and Manifold Extension Sets20

Administration Sets23

Blood Sets.....24

Blood Sets – Bifuse and Y-Type25

Burette Sets.....25

Peripheral IV Catheters26

Tracheostomy.....28

Flexextend TTS.....28

TTS28

Respiratory29

Manual Resuscitation29

HME.....29

Bronchial Hygiene30

Port Access Needles.....31

Gripper Plus™ Safety Port Access Needle31

Gripper Plus™ Power P.A.C. Safety Port Access Needle31

Gripper Micro™ Safety Port Access Needle32

Gripper™ Port Access Needle32

Blood Collection.....33

Blood Diversion.....33

Blood Draw33

Arterial Blood Sampling35

Hypodermic Injections.....37

Enteral Feeding Accessories38



HOME

Needlefree Connectors

Catalog #	Description	PV (mL)	Case Qty
MC100	MicroClave™ Clear Neutral Connector	0.04	100

Catalog #	Description	PV (mL)	Case Qty
B3300	Clave™ Neutral Connector	0.04	100


Catalog #	Description	PV (mL)	Case Qty
A1000	NanoClave™ Connector	0.02	100

Catalog #	Description	PV (mL)	Case Qty
NC100	Neutron™ Connector	0.06	100


Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.




HOME

Catalog #	Description	PV (mL)	Case Qty
D1005	Tego™ Connector	0.06	100
			

Disinfection Caps




Catalog #	Description	Case Qty
SCXT3-2000 ⁺	SwabCap™, 200-ct. Box	2,000
SCXT3-10-2000 ⁺	SwabPack™, 10-ct. Bag	2,000
SCXT3-2400	SwabPack™, 25-ct. Bag	2,400
SCRC3-10-1600	SwabPack Plus™	1,600
SCXT3-5-2000	SwabSleeve™, 5-ct. Sleeve	2,000


Catalog #	Description	Case Qty
SWTP-5	SwabTip™, 5-ct. Strip	200
		


⁺ Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Single Lumen Extension Sets


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
A1007	7" (18 cm) Smallbore Ext Set w/MicroClave™ Clear, NanoClave™, Clamp, Rotating Luer	7	0.29	50
				

Single Lumen Extension Sets — Filtered


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B1478	7" (18 cm) Smallbore Ext Set w/1.2 Micron Filter, Clamp, Rotating Luer	7	0.46	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B1479	8" (20cm) Smallbore Ext Set w/0.2 Micron Filter, Clamp, Luer Lock	8	0.45	50
				

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33933■	7" (18 cm) Appx 0.39 mL, Smallbore Ext Set w/MicroClave™ Clear, 0.2 Micron Filter, Clamp (Blue), Luer Lock	7	0.39	50
				


Single Lumen Extension Sets T Connector


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
A1001	6" (15 cm) Smallbore Ext Set w/Clamp, NanoClave™ T-Connector, Rotating Luer	6	0.26	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
A1099	6.5" (16.5 cm) Approx 0.31 mL, Smallbore Pressure Infusion (400 psig) Ext Set w/Remv MicroClave™ Clear, NanoClave™ T-Connector, Clamp (Purple), Rotating Luer	6.5	0.31	50
				

■ May require longer leadtime
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Single Lumen Extension Sets High Pressure


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33038	7" (18 cm) Appx 0.30 mL, Smallbore Pressure Infusion (400 psig) Ext Set w/Remv MicroClave Clear, Purple Clamp, Luer Lock, Non-DEHP Tubing	7	0.30	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33150●	7" (18 cm) Approx 0.24 mL, Smallbore Pressure Infusion (400 psig) Ext Set w/MicroClave™ Clear, Clamp (Purple), Rotating Luer	7	0.24	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
A1201	6.5" (17 cm) Approx 0.35 mL, Smallbore Pressure Infusion (400 psig) Ext Set w/Remv MicroClave™ Clear, NanoClave™, Clamp (Purple), Rotating Luer	6.5	0.35	50
				


● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33086	7" Smallbore Pressure Infusion (400 psig) Ext Set w/MicroClave™ Clear, Purple Clamp, Luer Lock	7	0.21	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33213	7" (18 cm) Appx 0.31 mL, Smallbore Pressure Infusion (400 psig) Ext Set w/Remv MicroClave Clear, Purple Clamp, Rotating Luer	7	0.31	50
				

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Multi Lumen Extension Sets


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33119	4" (10 cm) Approx 0.36 mL, Smallbore Bifuse Ext Set w/2 Remv MicroClave™ Clear, 2 Clamps, Rotating Luer, Non-DEHP Tubing	4	0.36	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33315	5" (13 cm) Approx 0.3 mL, Smallbore Bifuse Ext Set w/2 MicroClave™ Clear, 2 Clamps, Rotating Luer	5	0.3	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33098	4" (10 cm) Approx 0.4 mL, Smallbore Trifuse Ext Set w/3 MicroClave™ Clear, 3 Clamps, Rotating Luer	4	0.4	50
				


● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33352	60" (152 cm) Approx 0.36 mL, Smallbore Ext Set w/MicroClave™ Clear, Clamp, Luer Lock	60	0.36	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33383	60" (152 cm) Approx 0.54 mL, Smallbore Ext Set w/MicroClave™ Clear, Clamp, Luer Lock	60	0.54	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330425	60" (152 cm) Approx 1.0 mL, Smallbore Ext Set w/MicroClave™ Clear, Clamp (Blue), Luer Lock	60	1.0	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33932	60" (152 cm) Approx 1.0 mL, Smallbore Ext Set w/MicroClave™ Clear, Clamp, Luer Lock	60	1	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B33711	61" (155 cm) Approx 1.0 mL, Smallbore Ext Set w/MicroClave™, Clamp, Luer Lock	61	1	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33514	60" (152 cm) Approx 1.7 mL, Smallbore Ext Set w/MicroClave™ Clear, Clamp (Blue), Luer Lock	60	1.7	50
				


● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B2154	360" (914 cm) Approx 10.1 mL, Smallbore Ext Set w/Clamp, Luer Lock	360	10.1	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B2201	360" PUR Smallbore Coiled Ext Set w/Luer Lock	360	7.2	25
				


Medication Sets — Filtered


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330472 	60" (152 cm) Approx 1.4 mL, PUR Yellow Smallbore Ext Set w/MicroClave™ Clear, 1.2-Micron Filter, Clamp, Luer Lock	60	1.4	50
				

 Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330713	60" (152 cm) Appx 0.55 mL, Smallbore Ext Set w/Clave™ Clear, 0.2 Micron Filter, Clamp, Luer Lock	60	0.55	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330174	60" (152 cm) Approx 1.9 mL, Smallbore Ext Set w/MicroClave™ Clear, 1.2-Micron Filter, 2 Clamps, Luer Lock	60	1.9	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330077 	60" (152 cm) Approx 1.9 mL, Smallbore Ext Set w/MicroClave™ Clear, 0.2-Micron Filter, Clamp, Luer Lock	60	1.9	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC33425 	60" (152 cm) Smallbore Ext Set w/Remv MicroClave™ Clear, 0.2 Micron Filter, Clamp, Luer Lock	60	0.78	50
				

 Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
 May require longer leadtime
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B33895	107" (272 cm) Transfer Set w/Dual Check Valve, MicroClave™, Luer Lock	107	8.7	50




Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330082	109" (277 cm) Transfer Set w/2 MicroClave™ Clear, Dual Check Valve, 3 Clamps (Red, White, Blue), Rotating Luer	109	7.4	50



Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B33065	104" (264 cm) Approx 7.1 mL, Transfer Set w/Clave™, Dual Check Valve, Smallbore Trifuse Ext Set w/2 Clave™, Rotating Luer	104	7.1	50



Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC330442	111" (282 cm) Transfer Set w/MicroClave™ Clear, Dual Check Valve, Smallbore Quadfuse Ext Set w/3 MicroClave™ Clear, 5 Clamps (2 Red, Blue, 2 White), Rotating Luer, Non-DEHP Tubing	111	7.7	50




● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Stopcocks and Manifolds


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AC100	NanoClave™ Stopcock	1.9	0.28	50




Catalog #	Description	Length (in.)	PV (mL)	Case Qty
SC100	Clave™ Stopcock	1.37	0.27	50




Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AC104■	NanoClave™ Stopcock w/NanoClave™	2.22	0.3	50




Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AC200	2 Gang 4-Way NanoClave™ Stopcock Manifold w/Rotating Luer	3.81	0.55	50



Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AC205	5 Gang 4-Way NanoClave™ Stopcock w/Baseplate, Rotating Luer	7	1.4	50




Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM3016	3-Port NanoClave™ Manifold w/Check Valve, Rotating Luer	4	0.88	50





■ May require longer leadtime
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6099	6-Port NanoClave™ Manifold, Check Valve, Rotating Luer	3.1	0.75	50
				


Stopcock and Manifold Extension Sets


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
SC124	30" (76 cm) Approx 3.6 mL, Ext Set w/Remv 2 Gang Clave™ Stopcock, Remv Clave™ Stopcock, NanoClave™ T-Connector, Rotating Luer, 1 Ext	29.314	3.6	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM3105	47" (119 cm) Approx 6.6 mL, Ext Set w/3-Port NanoCalve™ Manifold, 2 MicroClave™ Clear, Check Valve, 2 Clamps, Rotating Luer	47	6.6	50
				


■ May require longer leadtime
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM3127	11" (28 cm) Approx 1.1 mL, Smallbore Ext Set w/3-Port NanoClave™ Manifold, Check Valve, Clave Clear, Clamp, Luer Lock	11	1.1	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6109	10" (25 cm) Approx 0.95 mL, Smallbore Ext Set w/6-Port NanoClave™ Manifold, Check Valve, Clamp, Rotating Luer	10	0.95	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6100	10" (25 cm) Approx 1 mL, Smallbore Ext Set w/6-Port NanoClave™ Manifold (Orange, Red, Blue, Purple, Yellow Rings), Check Valve, Clamp, Rotating Luer	10	1	50
				


● Included in ICU Medical IV Set Essential portfolio of infusion technologies to help support continuity and standardization
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6127■	11" (28 cm) Smallbore Ext Set w/NanoClave™, 6-Port NanoClave™ Manifold, Check Valve, Clamp, Rotating Luer	11	0.97	50
				


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6140■	8" (20 cm) Approx 1.2, mL Smallbore Ext Set w/NanoClave™, 6-Port NanoClave™ Manifold, NanoClave™ 4-Way Stopcock, Check Valve, Clamp, Rotating Luer	8	1.2	50
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6152	20" (51 cm) Approx 1.9 mL, Ext Set, Smallbore w/Clave™ Clear, 2 6-Port NanoClave™ Manifolds, 2 Check Valves, Spin Luer	20	1.9	25
				


■ May require longer leadtime
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
AM6154	10" (25 cm) Appx 1.0 mL, Smallbore Ext Set w/6-Port NanoClave™ Manifold, Check Valve (Orange, Glow, Blue, Purple, Yellow, Green Rings), NanoClave™, Luer Lock	60	1	50
				


Administration Sets

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
20671-28	Primary Microdrip Set, Piggyback with Backcheck Valve, 3 Clave Y-Sites, Secure Lock, 100 Inch	100	15	48
				

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
SC125	120" (305 cm) Approx 15.4 mL, 10-Drop Primary Set w/2 Clave™, Remv 2 Gang Clave™ Stopcock, Remv Clave™ Stopcock, NanoClave™ T-Connector, Rotating Luer, 2 Check Valves, 2 Ext	120	15.4	25
				


Blood Sets

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B31753	27" (69 cm) 10-Drop Blood Set w/200-Micron Filter, Rotating Luer	27	3.6	25
				


The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Blood Sets — Bifuse and Y Type


Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B30061	82" (208 cm) Y-Type 10-Drop Blood Set w/200-Micron Filter, Rotating Luer	82	16.4	25
				

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
MC9011	99" (251 cm) 20-Drop Y-Type Blood Set w/200-Micron Filter, 2 Clave™ Clear, Rotating Luer	99	20.32	25
				

Burette Sets

Catalog #	Description	Length (in.)	PV (mL)	Case Qty
B9409	119" (302 cm) 60 Drop 150 mL Burette Set (w/Shut-Off), 4 Clave™, Rotating Luer	119	14.7	25
				

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.



Item	Description	Gauge	Length (in.)	Material	Winged
SP125-26-14T	SuperCath® 26G X 0.56 in Active Straight	26	0.56	PUR	N
SP125-26-14WT	SuperCath® 26G X 0.56 in Active Straight Winged	26	0.56	PUR	Y
SP125-26-19T	SuperCath® 26G X 0.75 in Active Straight	26	0.75	PUR	N
SP125-24-14T	SuperCath® 24G X 0.56 in Active Straight	24	0.56	PUR	N
SP125-24-14-WT	SuperCath® 24G X 0.56 in Active Straight Winged	24	0.56	PUR	Y
SP125-24-19T	SuperCath® 24G X 0.75 in Active Straight	24	0.75	PUR	N
SP125-24-19-WT	SuperCath® 24G X 0.75 in Active Straight Winged	24	0.75	PUR	Y
SP120-24-19T	SuperCath® 24G X 0.75 in Active Straight Blood Containment	24	0.75	PUR	N
SP120-24-19-WT	SuperCath® 24G X 0.75 in Active Straight Blood Containment Winged	24	0.75	PUR	Y
SP120-22-25T	SuperCath® 22G X 1.00 in Active Straight Blood Containment	22	1	PUR	N
SP120-22-25-WT	SuperCath® 22G X 1.00 in Active Straight Blood Containment Winged	22	1	PUR	Y
SP120-22-31T	SuperCath® 22G X 1.25 in Active Straight Blood Containment	22	1.25	PUR	N
SP120-22-31-WT	SuperCath® 22G X 1.25 in Active Straight Blood Containment Winged	22	1.25	PUR	Y
SP120-20-25T	SuperCath® 20G X 1.00 in Active Straight Blood Containment	20	1	PUR	N
SP120-20-25-WT	SuperCath® 20G X 1.00 in Active Straight Blood Containment Winged	20	1	PUR	Y
SP120-20-31T	SuperCath® 20G X 1.25 in Active Straight Blood Containment	20	1.25	PUR	N

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Item	Description	Gauge	Length (in.)	Material	Winged
SP120-20-31-WT	SuperCath® 20G X 1.25 in Active Straight Blood Containment Winged	20	1.25	PUR	Y
SP120-20-50T	SuperCath® 20G X 2.00 in Active Straight Blood Containment	20	2	PUR	N
SP120-20-50-WT	SuperCath® 20G X 2.00 in Active Straight Blood Containment Winged	20	2	PUR	Y
SP120-18-31T	SuperCath® 18G X 1.25 in Active Straight Blood Containment	18	1.25	PUR	N
SP120-18-31-WT	SuperCath® 18G X 1.25 in Active Straight Blood Containment Winged	18	1.25	PUR	Y
SP120-18-50T	SuperCath® 18G X 2.00 in Active Straight Blood Containment	18	2	PUR	N
SP120-18-50-WT	SuperCath® 18G X 2.00 in Active Straight Blood Containment Winged	18	2	PUR	Y
SP120-16-31T	SuperCath® 16G X 1.25 in Active Straight Blood Containment	16	1.25	PUR	N
SP120-16-31-WT	SuperCath® 16G X 1.25 in Active Straight Blood Containment Winged	16	1.25	PUR	Y
SP125-16-31T	SuperCath® 16G X 1.25 in Active Straight	16	1.25	PUR	N
SP125-16-31-WT	SuperCath® 16G X 1.25 in Active Straight Winged	16	1.25	PUR	Y
SP125-14-31T	SuperCath® 14G X 1.25 in Active Straight	14	1.25	PUR	N
SP125-14-31-WT	SuperCath® 14G X 1.25 in Active Straight Winged	14	1.25	PUR	Y

SuperCath is a registered trademark of MediKit Co., Ltd.

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Tracheostomy

Flextend™ TTS™

	Flextend TTS		
	Product Code	Description	Units Per Case
	67NFPS25	NEONATAL TRACHEOSTOMY TUBE 2.5 TTS FLEXTEND STRAIGHT FLANGE 1/EA	1
	67NFPS30	NEONATAL TRACHEOSTOMY TUBE 3.0 TTS FLEXTEND STRAIGHT FLANGE 1/EA	1
	67NFPS35	NEONATAL TRACHEOSTOMY TUBE 3.5 TTS FLEXTEND STRAIGHT FLANGE 1/EA	1
	67NFPS40	NEONATAL TRACHEOSTOMY TUBE 4.0 TTS FLEXTEND STRAIGHT FLANGE 1/EA	1

	Flextend TTS		
	Product Code	Description	Units Per Case
	67PFSS25	PEDIATRIC TRACHEOSTOMY TUBE 2.5 TTS FLEXTEND STRAIGHT NECK 1/EA	1
	67PFSS30	PEDIATRIC TRACHEOSTOMY TUBE 3.0 TTS FLEXTEND STRAIGHT NECK 1/EA	1
	67PFSS35	PEDIATRIC TRACHEOSTOMY TUBE 3.5 TTS FLEXTEND STRAIGHT NECK 1/EA	1
	67PFSS40	PEDIATRIC TRACHEOSTOMY TUBE 4.0 TTS FLEXTEND STRAIGHT NECK 1/EA	1
	67PFSS45	PEDIATRIC TRACHEOSTOMY TUBE 4.5 TTS FLEXTEND STRAIGHT NECK 1/EA	1
	67PFSS50	PEDIATRIC TRACHEOSTOMY TUBE 5.0 TTS FLEXTEND STRAIGHT NECK 1/EA	1
	67PFSS55	PEDIATRIC TRACHEOSTOMY TUBE 5.5 TTS FLEXTEND STRAIGHT NECK 1/EA	1

TTS

	TTS		
	Product Code	Description	Units Per Case
	67SN025	2.5 TTS CUFFED NEONATAL TRACHEOSTOMY TUBE 1/EA	1
	67SN030	3.0 TTS CUFFED NEONATAL TRACHEOSTOMY TUBE 1/EA	1
	67SN035	3.5 TTS CUFFED NEONATAL TRACHEOSTOMY TUBE 1/EA	1
	67SN040	4.0 TTS CUFFED NEONATAL TRACHEOSTOMY TUBE 1/EA	1

	TTS		
	Product Code	Description	Units Per Case
	67P025	2.5 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1
	67P030	3.0 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1
	67P035	3.5 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1
	67P040	4.0 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1
	67P045	4.5 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1
	67P050	5.0 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1
	67P055	5.5 TTS CUFFED PEDIATRIC TRACHEOSTOMY TUBE 1/EA	1

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Respiratory

Manual Resuscitation

	Manual Resuscitation		
	Product Code	Description	Units Per Case
	8528BM	INFANT 1st Response™ MANUAL RESUSCITATOR WITH OXYGEN RESERVOIR BAG, MANOMETER AND NEONATE MASK	9
	8528M	INFANT REUSCITATOR WITH INFANT MASK, 20" RESERVOIR BAG, 9/CA	9
	8529BM	INFANT 1st Response™ MANUAL RESUSCITATOR WITH OXYGEN RESERVOIR BAG, MANOMETER AND INFANT MASK	9
	8529M	RESUSCITATOR, RESERVOIR BAG INFANT, FLEX 20IN, INFANT MASK 9/CA	9

	Manual Resuscitation		
	Product Code	Description	Units Per Case
	8520	RESUSCITATOR, PAED, P&I ELBOW CHILD MASK, TAIL RESERVOIR 9/CA	9
	8520B	PAED. MANUAL RESUSCITATOR 9/CA	9

HME

Catalog #	Description	Case Qty
570016	THERMOVENT™ T + 50/CA	50
		

Catalog #	Description	Case Qty
100/570/022	THERMOVENT T2 STERILE + 50/CA	50
		

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Catalog #	Description	Case Qty
580011	THERMOVENT 600 + 50/CA	50
		

Catalog #	Description	Case Qty
002815	10ML NEONATAL HME 50/CA	50
		

Bronchial Hygiene

Catalog #	Description	Case Qty
21-1015	ACAPELLA™ DM BLUE MOUTHPIECE NS + 10/CA	10
		

Catalog #	Description	Case Qty
22-2000	COACH 2 FOR KIDS™ NS + 12/CA	12
		

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Port Access Needles

Gripper Plus™ Safety Port Access Needle

Catalog #	Description	Case Qty
21-2767-24	NDL, GRIPPER PLUS, 20G X .75" (19 mm)	12
21-2965-24	NDL, GRIPPER PLUS, 20G X .75" (19 mm), NEEDLELESS Y-SITE	12
21-2865-24	NDL, GRIPPER PLUS, 20G X .75" (19 mm), Y-SITE	12
21-2762-24	NDL, GRIPPER PLUS, 20G X 1" (25 mm)	12
21-2966-24	NDL, GRIPPER PLUS, 20G X 1" (25 mm), NEEDLELESS Y-SITE	12
21-2866-24	NDL, GRIPPER PLUS, 20G X 1" (25 mm), Y-SITE	12
		

Gripper Plus™ Power P.A.C. Safety Port Access Needle

Catalog #	Description	Case Qty
21-3367-24	NDL, GRIPPER PLUS, PWR INJECTABLE, 20G X .75"	12
21-3465-24	NDL, GRIPPER PLUS, PWR INJECTABLE, 20G X .75", NEEDLELESS Y-SITE	12
21-3362-24	NDL, GRIPPER PLUS, PWR INJECTABLE, 20G X 1"	12
21-3466-24	NDL, GRIPPER PLUS, PWR INJECTABLE, 20G X 1", NEEDLELESS Y-SITE	12
		

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Gripper Micro™ Safety Port Access Needle

Catalog #	Description	Case Qty
21-3256-24	NDL, GRIPPER MICRO, 20G X .75"	12
21-3276-24	NDL, GRIPPER MICRO, 20G X .75", LUER ACTIVATED NEEDLELESS Y-SITE	12
21-3257-24	NDL, GRIPPER MICRO, 20G X 1"	12
21-3277-24	NDL, GRIPPER MICRO, 20G X 1", LUER ACTIVATED NEEDLELESS Y-SITE	12
		

Gripper™ Port Access Needle

Catalog #	Description	Case Qty
21-2734-24	NDL, GRIPPER, 20G X .75" (19 mm)	12
21-2947-24	NDL, GRIPPER, 20G X .75" (19 mm), Y-SITE	12
21-2717-24	NDL, GRIPPER, 20G X 1" (25 mm)	12
21-2948-24	NDL, GRIPPER, 20G X 1" (25 mm), Y-SITE	12
		

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Blood Collection

Blood Diversion

VI ByPass Syringe® Blue		
Catalog #	Description	Case Qty
VI-103-LC	Syringe, By-pass, Luer Connection	200
		

Blood Draw

neoheel™ Safety Lancets					
Catalog #	Description	Incision depth and length	Blood Volume	Color	Case Qty
1060M	Micro-Preemie — Low birth weight (<1kg/2.2lbs)	0.65 mm X 1.40 mm	Low Flow	●	4 boxes of 50 each
1061P	Preemie — Low birth weight (1-2.5kg/2.2-5.5lbs)	0.85 mm X 1.75 mm	Medium Flow	●	4 boxes of 50 each
1062N	Newborn — Birth to 6 months (2.5-9kg /5.5-19.8lbs)	1.0 mm X 2.5 mm	High Flow	●	4 boxes of 50 each
1063T	Toddler — 6 to 24 months (>9kg - >19.8lbs)	2.0 mm X 3.0 mm	Maximum Flow	●	4 boxes of 50 each
					

VI ByPass Syringe is a registered trademark of Vascular Integrity.
The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Safe-T Lance™ Plus Safety Lancets				
Catalog #	Description	Blade Width	Hub Color	Case Qty
9662	Safe-T-Lance Pediatric Blade	1.5 mm	●	10 boxes of 200 each
				

Saf-T Wing™ Blood Collection and Infusion Set			
Catalog #	Description	Color	Case Qty
972306	23G x ¾" with 6" Tubing	Blue	4 boxes of 50
972312	23G x ¾" with 12" Tubing	Blue	4 boxes of 50
972506	25G x ¾" with 6" Tubing	Orange	4 boxes of 50
972512	25G x ¾" with 12" Tubing	Orange	4 boxes of 50
982306	23G x ¾" with 6" Tubing and Saf-T Holder Needle Protection Device	Blue	4 boxes of 50
982312	23G x ¾" with 12" Tubing and Saf-T Holder Needle Protection Device	Blue	4 boxes of 50
982506	25G x ¾" with 6" Tubing and Saf-T Holder Needle Protection Device	Orange	4 boxes of 50
982512	25G x ¾" with 12" Tubing and Saf-T Holder Needle Protection Device	Orange	4 boxes of 50
			

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Saf T Holder™ Devices		
Catalog #	Description	Case Qty
96000	Saf-T Holder with Multi-Sample Male Luer Adapter	4 boxes of 50
96000S	Short Saf-T Holder with Multi-Sample Male Luer Adapter	4 boxes of 50
96002	Saf-T Holder with Female Luer Adapter	4 boxes of 50
96002S	Short Saf-T Holder with Female Luer Adapter	4 boxes of 50
96004	Saf-T Holder Blood Culture Device with Male Luer Adapter	4 boxes of 50
96005	Saf-T Holder Blood Culture Device with Female Luer Adapter	4 boxes of 50
96007	Pediatric Tube Insert	4 boxes of 50
		

Arterial Blood Sampling

	Portex™ LineDraw™ Plus Arterial Blood Sampling Kits Dry Heparin			
	200 Units Per Case		Product Codes	
	Components	4040-2 Umbilical	4041LH	4042LH
	1 mL Luer Slip Syringe		•	
	3 mL Luer Lock Syringe	•		•
	22G x 1 ½"	•		
	Filter-Pro™ Air Bubble Removal Device	•	•	•
	Point-Lok™ Needle Protection Device	•		
	Tip Cap		•	
	Sterilized			
	Fully Sterile	•	•	•

*LineDraw Plus syringes contain approximately 23.5 I.U. per mL of advanced formula dry lithium heparin.
*LineDraw Plus Low Heparin syringes (LH Codes) contain approximately 7 I.U. per mL of dry lithium heparin.

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

	Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits–Full*		
	100 Units Per Case Components	4611LH	4640LH
	1 mL Luer Slip Syringe	•	•
	Needle		
	25G x 5/8" Pre-Attached to Syringe	•	
	25G x 5/8" Pre-Attached with Needle-Pro™ to Syringe		•
	23G x 1" Secondary Needle with Needle-Pro™		•
	Filter-Pro™ Air Bubble Removal Device	•	•
	Point-Lok™ Needle Protection Device	•	
	Alcohol Prep Pad	•	•
	Iodophor Prep Pad	•	•
	Gauze Pad (2"x 2")	•	•
	Patient Label	•	•
	Biohazard Labeled Ice Bag	•	•
	Adhesive Bandage (1" x 3")	•	•
	Sterilized		
	Fully Sterile	•	•

	Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits–Strip*		
	200 Units Per Case Components	Product Codes	
		4610LH	4619LH
	1 mL Luer Slip Syringe	•	•
	Needle		
	25G x 5/8" Pre-Attached to Syringe	•	
	25G x 5/8" Pre-Attached with Needle-Pro™ to Syringe		•
	Filter-Pro™ Air Bubble Removal Device	•	•
	Point-Lok™ Needle Protection Device	•	
	Sterilized		
	Fully Sterile	•	•

*Pro-Vent™ Low Heparin (LH Codes) 1 mL syringes contain approximately 7 I.U. per mL of dry lithium heparin.

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Hypodermic Injections

Hypodermic Needle Pro EDGE™ Safety Device

Catalog #	Description	Needle Gauge	Needle Length	Hub Color	Case Qty
401810	HYPODERMIC NEEDLE-PRO EDGE 18G x 1"	18	1"	Pink	10 boxes of 100
402510	HYPODERMIC NEEDLE-PRO EDGE 25G x 1"	25	1"	Orange	10 boxes of 100
402558	HYPODERMIC NEEDLE-PRO EDGE 25G x 5/8"	25	5/8"	Orange	10 boxes of 100
402605	HYPODERMIC NEEDLE-PRO EDGE 26G x 0.5"	26	0.5"	Brown	10 boxes of 100
402705	HYPODERMIC NEEDLE-PRO EDGE 27G x 0.5"	27	0.5"	Gray	10 boxes of 100
4027125	HYPODERMIC NEEDLE-PRO EDGE 27G x 1.25"	27	1.25"	Gray	10 boxes of 100
403005	HYPODERMIC NEEDLE-PRO EDGE 30G x 0.5"	30	0.5"	Yellow	10 boxes of 100




TB Syringes with Needle Pro EDGE Hypodermic Needles


Catalog #	Description	Needle Gauge	Needle Length	Hub Color	Case Qty
412558	NEEDLE-PRO EDGE TB 25G x 5/8" with 1 mL Syringe	25	5/8"	Orange	8 boxes of 50
412705	NEEDLE-PRO EDGE TB 27G x 1/2" with 1 mL Syringe	27	5/8"	Grey	8 boxes of 50
432510	NEEDLE-PRO EDGE TB 25G x 1" with 3 mL Syringe	25	1"	Orange	8 boxes of 50
432558	NEEDLE-PRO EDGE TB 25G x 5/8" with 3 mL Syringe	25	5/8"	Orange	8 boxes of 50




The images shown are for illustration purposes only and may not be an exact representation of the product sizing.


Blunt Fill Needle and Blunt Fill Needle with Filter			
Catalog #	Description	Hub Color	Case Qty
BN1815	18G x 1 1/2" Luer Lock	Red	10 boxes of 100
BN1815N	18G x 1 1/2" Luer Lock Filtered Needle	Pink	10 boxes of 100
			


Enteral Feeding Accessories

Catalog #	Description	Case Qty
M9000	Lopez Valve™, Non-Sterile	50
		

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

Catalog #	Description	Case Qty
M9000-T	Lopez Valve, Tethered Cap, Sterile	50
		

Catalog #	Description	Case Qty
M9000-EFD	ENFit™ Lopez Valve, Non-Sterile	50
		

Catalog #	Description	Case Qty
M9000-SEFD	ENFit Lopez Valve, Sterile	50
		

The images shown are for illustration purposes only and may not be an exact representation of the product sizing.

A Specialized Portfolio for Your Smallest Patients

The smallest patients need specialized care, and clinicians like you need the right infusion systems, consumables, and solutions to help deliver it. ICU Medical offers a broad portfolio specifically designed to help support your NICU/PICU workflows and enhance the delivery of patient care.



Infusion pumps and software

Infusion systems including syringe pumps, large-volume pumps, pain and ambulatory pumps, and infusion safety software




Infusion consumables

IV sets featuring the infection control properties of advanced Clave™ needlefree connectors



Vascular access

Vascular access tools including blood collection technology, small-gauge catheters, hypodermics, and implantable ports



IV solutions

IV solutions for irrigation, infusion, and parenteral nutrition, with non-DEHP/ non-PVC options and supported by robust manufacturing capacity, product quality, warehousing, and distribution capabilities




Renal care

Innovative technologies designed to help protect hemodialysis catheters from contamination and support decreased risk of CRBSI.



Airway Management

A broad portfolio of tracheostomy, respiratory, breathing circuits, ventilation, and specialty airway products.



icumedical
human connections