




Oncology Quick Reference Catalogue

Table of Contents

ChemoClave™ Closed System Transfer Device

Closed Male Luer	3
Vial Spikes.....	4
Bag Spikes	7
Add-On & Transfer Sets	8
Secondary Sets	9
Adaptors	10

MicroClave™ Needlefree Connector..... 11

ChemoLock™ Closed System Transfer Device


Injector.....	11
Port.....	12
Vial Spikes.....	13
Bag Spikes	15
Secondary Sets	17
Add-On Set.....	19
Adaptor	19




HOME

ChemoClave Closed System Transfer Device


Closed Male Luer


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011- CH2000S-C	Spinning Spiros™	FSB1277	0.1	50
				


Item Number	Description	PV (mL)	Case Qty
011- CH2000S-PC	Spinning Spiros Closed Male Luer, Purple Cap	0.1	50
			


Item Number	Description	PV (mL)	Case Qty
011-CH-CAP	Cap, Sterile	0	50
			


Vial Spikes


Item Number	Description	PV (mL)	Case Qty
011-CS-50	Clave™ Connector Multidose Vial Access Spike	0.19	50
			


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CS-51	Universal Vented Vial Spike w/ChemoClave	FSW1491	0.11	50
				


Item Number	Description	PV (mL)	Case Qty
CS-52	Chemo Dispensing Pin w/5 Micron Filter, Clave, Appx 0,38 mL	0.38	50
			


Item Number	Description	PV (mL)	Case Qty
011-CH-53	ChemoClave Vented Vial Spike — No Skirt	0.14	50
			


Item Number	Description	PV (mL)	Case Qty
011- CH-62	ChemoClave Vial Spike, 13 mm	0.07	50
			


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CH-70	ChemoClave Universal Vented Vial Spike	FSB1283	0.14	50
				


Item Number	Description	PV (mL)	Case Qty
011-CH-70S	ChemoClave Vented Vial Spike, 20 mm	0.14	50
			

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CH-72	ChemoClave Vented Vial Spike, 13 mm	FSB1412	0.08	50
				


Item Number	Description	PV (mL)	Case Qty
011-CH-74	ChemoClave Universal Vented Vial Spike	0.14	50
			


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CH-80	ChemoClave Closed Vial Spike	FSB2437	0.14	50
				


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CH-80S	ChemoClave Vial Spike, 20 mm	FSB2052	0.14	50
				


Item Number	Description	PV (mL)	Case Qty
011-CH-82	ChemoClave 13 mm Closed Vial Spike	0.11	50
			


Bag Spikes

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CS-10	Clave Connector IV Bag Access Device	FSB1230	0.39	50
				


Item Number	Description	PV (mL)	Case Qty
011-CS-11	Clave, IV Bag Access Device	0.52	50
			


Item Number	Description	PV (mL)	Case Qty
011-CH-12	ChemoClave Bag Spike w/Port, Dry Spike	0.35	50
			


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CH-14	ChemoClave Vented Bag Spike	FSB1285	0.34	50
				


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CS-15	Clave Conn, IV Bag Access Device, BCV	FSB991	0.53	50
				

Add-On & Transfer Sets



Item Number	Description	NHSSC Code	Length (cm)	PV (mL)	Case Qty
011-CH3034	Transfer Set w/Spiros	FSB1278	11	1.3	50
					



Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-H3161	30 cm (12") Add-On Set w/Spiros, Red Cap, 0.2 Micron Filter, Vented Cap	30	5.1	50
				


Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-H3830	53 cm (21") Transfer Set w/2 Bag Spikes	53	2.6	50
				


Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-C3002	33 cm (13") Transfer Set w/Clamp, Rotating Luer	33	2.0	50
				

Secondary Sets


Item Number	Description	NHSSC Code	Length (cm)	PV (mL)	Case Qty
011-CH3261	76 cm (30") Appx 3.3 mL, Admin Set, 2 Spiros w/Red Cap, 20 Drop In-Line Drip Chamber w/15 Micron Filter, Bag Hanger	FSB2319	76	3.3	50
  					


Item Number	Description	NHSSC Code	Length (cm)	PV (mL)	Case Qty
CH3235	76 cm (30") Appx 4.1 mL, Yellow 20 Drop Admin Set w/15 Micron Filter, Spiros	FSB1240	76	4.1	50
 					


Item Number	Description	NHSSC Code	Length (cm)	PV (mL)	Case Qty
011-CH3304	76 cm (30") 20 Drop Admin Set w/Integrated ChemoClave Drip Chamber, Bag Hanger, Rotating Luer	FSB1973	76	3.3	50
 					

Item Number	Description	NHSSC Code	Length (cm)	PV (mL)	Case Qty
011-CH3751	Admin Set, Amber, Spiros	FSB2490	76	6.9	50
					

Adaptors

Item Number	Description	PV (mL)	Case Qty
CH-33	Double Clave Syringe Transfer Set	0.15	50
			

Item Number	Description	PV (mL)	Case Qty
CH-33C	Double MicroClave Clear Syringe Transfer Set	0.15	50
			


Item Number	Description	PV (mL)	Case Qty
011-CH3900	Graduated Adapter with Clave	0.91	50
			


MicroClave Needlefree Connector

Item Number	Description	PV (mL)	Case Qty
011-MC100	MicroClave™ Clear Connector	0.05	100
			

ChemoLock Closed System Transfer Device

Injector

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL2000S	ChemoLock	FSW1565	0.35	50
				


Item Number	Description	PV (mL)	Case Qty
011-CL2000S-RC	ChemoLock w/Red Cap	0.35	50
			




HOME


Item Number	Description	Case Qty
011-CL-CAP	Cap, Sterile, Red	100
		


Port


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL2100	ChemoLock Port	FSW1492	0.1	50
				

Vial Spikes

Item Number	Description	PV (mL)	Case Qty
011-CL-51	ChemoLock Vented Vial Spike	0.15	50
			


Item Number	Description	PV (mL)	Case Qty
011-CL-53	Spike with Bonded ChemoLock Port	0.2	50
			


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-62	ChemoLock Vial Spike, 13 mm	FSW85025	0.13	50
				


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-70	ChemoLock Universal Vented Vial Spike	FSW85023	0.18	50
				





HOME


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-70S	ChemoLock Vented Vial Spike, 20 mm	FSW85026	0.18	50
				

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-72	ChemoLock Vented Vial Spike, 13 mm	FSW85027	0.13	50
				


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-80	ChemoLock Closed Vial Spike	FSW85028	0.18	50
				


Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-80S	ChemoLock Vial Spike, 20 mm	FSW85029	0.18	50
				


Item Number	Description	PV (mL)	Case Qty
011-CL-80SL	ChemoLock Vial Spike, 28 mm	0.18	50
			

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-82	ChemoLock 13 mm Closed Vial Spike	FSW85030	0.16	50
				


Bag Spikes

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-10	ChemoLock Bag Spike	FSW85024	0.43	50
				

Item Number	Description	PV (mL)	Case Qty
011-CL-12	ChemoLock Bag Spike w/Port, Dry Spike	0.39	50
			


Item Number	Description	PV (mL)	Case Qty
011-CL-13	ChemoLock Bag Spike with Additive Port	1.54	50
			

Item Number	Description	NHSSC Code	PV (mL)	Case Qty
011-CL-14	ChemoLock Vented Bag Spike	FSW1568	0.38	50
				


Item Number	Description	PV (mL)	Case Qty
011-CL3367	ChemoLock Bag Spike w/Port, ChemoClave	0.45	50
			

Secondary Sets


Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-CL3261	76 cm (30") Appx 3.9 mL Admin Set w/2 ChemoLock, 20 Drop In-Line Drip Chamber w/15 Micron Filter, Bag Hanger, Drop-In ChemoLock Port	76	3.9	50




Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-CL3511	76 cm (30") Appx 3.7 mL, 20 Drop Admin Set w/Integrated ChemoLock Port Drip Chamber, ChemoLock, Bag Hanger, Drop-In ChemoLock Port	76	3.7	50




Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-CL3520	102 cm (40") Appx 5.3 mL Admin Set w/2 ChemoLock, 20 Drop In-Line Drip Chamber, Bag Hanger, Drop-in ChemoLock Port	102	5.3	50



Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-CL3028	Admin Set w/ChemoLock Port, Spiros, Hanger	76	6.1	50




Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-CL4100	152 cm (60") Appx 7.4 mL, Admin Set w/ChemoLock, 20 Drop In-Line Chamber, Spiros w/Red Cap, 2 Bag Hangers	152	7.4	50




Add-On Set

Item Number	Description	NHSSC Code	Length (cm)	PV (mL)	Case Qty
011-CL3534	13 cm (5") Appx 1.8 mL, Bag Spike Adapter w/ChemoLock, Vented Cap	FSW1566	13	1.8	50




Item Number	Description	Length (cm)	PV (mL)	Case Qty
011-CL3535	7" (18 cm) Appx 0.93 mL, Ext Set w/ChemoLock, Y-Connector, Rotating Luer	18	0.93	50



Adaptor

Item Number	Description	PV (mL)	Case Qty
011-CL3900	Graduated Adapter w/ChemoLock Port	0.95	50



The product complies with current legislation and has the corresponding CE marking.
For additional information, warnings and /or safety precautions, refer to the manufacturer's Instructions for Use.

icumedical
human connections