

0.9% Sodium Chloride Irrigation USP

STERICARE SOLUTIONS

Item 6380

500 mL (16.9 fl oz)

Indications:
For moistening of wound dressings, wound debridement, and device irrigation.

Not for injection.
No antimicrobial or other substance added. Sterile, isotonic, single use. Discard unused portion.

Sterility Integrity Test:
Steris Ltd. (Steris Ltd.)
Nurse Assist, LLC
4409 Holcom Road,
Haltom City, TX 76117
800-649-6800
www.stericaresolutions.com

AP04131

Isotonic Solution for Irrigation Only.

pH: 4.5 - 7.0

PVC-free and DEHP-free

Osmolarity Calc.: 308 mOsmo/liter

Caution: Warm in oven to not more than 50°C (122°F) for a maximum of 60 days. Discard after 60 days of warming.

Do not use unless solution is clear and seal is intact. Recommended Storage: Room temperature. Avoid excessive heat. Protect from freezing.

Rx Only

Made in USA



2.00 in

3.00 in