

Lifecare PCA™ Set Conversions to Filtered Sets

June 2019




The three primary sets for LifeCare PCA will now include a 0.2 micron filter. Details about the set changes are outlined below. None of the current list numbers for the sets will change. These are permanent changes to the sets.

List Number 14277-28 (list number will not change)

Y-Type PCA Set, Minibore, PAV, 0.2 Micron Filter, Backcheck Valve, 67 Inch

	Normal Length	Filter	Priming Capacity (Approximate)	Tubing I.D.	Case Size	Latex Information*	Non-DEHP
Previous	67" / 170.18 cm	No	2.3 mL	0.033"	1 x 50	Yes	Yes
Conversion	67" / 170.18 cm	Yes	2.8 mL	0.033"	1 x 50	Yes	Yes



* Natural rubber latex has not been used in the manufacture of this device.

List Number 14279-28 (list number will not change)


Y-Type PCA Set, Minibore, PAV, 0.2 Micron Filter, Backcheck Valve, 86 Inch

	Normal Length	Filter	Priming Capacity (Approximate)	Tubing I.D.	Case Size	Latex Information*	Non-DEHP
Previous	86" / 218.44 cm	No	2.6 mL	0.054"	1 x 50	Yes	Yes
Conversion	86" / 218.44 cm	Yes	3.1 mL	0.033"	1 x 50	Yes	Yes

Previous (unfiltered)



Conversion (filtered)




List Number 14280-28 (list number will not change)

Straight-Line PCA Set, Minibore, PAV, 0.2 Micron Filter, 80 Inch

	Normal Length	Filter	Priming Capacity (Approximate)	Tubing I.D.	Case Size	Latex Information*	Non-DEHP
Previous	80" / 203.20 cm	No	1.5 mL	0.033"	1 x 50	Yes	Yes
Conversion	80" / 203.20 cm	Yes	2.0 mL	0.033"	1 x 50	Yes	Yes

Previous (unfiltered)



Conversion (filtered)

