




Transpac™ IV

Disposable Pressure Transducer System

Provides real-time invasive blood pressure measurements needed to manage critically ill patients.

icumedical
human connections

The Transpac IV pressure monitoring system provides you with over 30 years of technology in critical care settings and is available in multiple configurations accommodating all patient populations.



Benefits of Transpac IV

- › Available in multiple standard and custom configurations to meet your clinical needs
 - Includes single, double, triple and quad line options
- › Kit design for adult, pediatric and neonatal patient populations
- › Available with SafeSet™ blood conservation system
- › Clear fluid pathway facilitates priming during setup
- › Optional color coded lines for easy identification
- › High durometer pressure tubing ensures optimal dynamic response
- › Two Flush devices to choose from
 - 3mL squeeze flush
 - 30mL squeeze flush
- › Meets requirements for accuracy established by the Association for the Advancement of Medical Instrumentation (AAMI)

Product Specifications

Operating Pressure	-50 to 300 mm Hg
Overpressure Limits	-400 to 5000 mm Hg
Sensitivity	5 microvolt/volt/mm Hg $\pm 1\%$ @ 6 VDC and 22° C
Zero Offset	<25 mm Hg
Zero Drift	<2 mm Hg in 8 hours
Input Impedance	300-350 ohms
Output Impedance	300 ohms Nominal



To learn more about Transpac IV pressure monitoring system, please visit www.icumed.com

icumedical
human connections