



Vascular Access

Product Catalogue

icumedical
human connections

Table of Contents

Arterial Blood Sampling	4
Pro-Vent Plus™ 3 mL Arterial Blood Sampling Kits	4
Pro-Vent Plus™ 1 mL Arterial Blood Sampling Kits	4
Pulsator™ Arterial Blood Sampling Full and Strip Kits	5
Portex™ LineDraw™ Plus Arterial Blood Sampling Kits Dry Heparin	5
Blood Collection	6
Saf-T Holder™ Devices and Blood Transfer Devices	6
Sol-Care™ Safety Blood Collection Needle with Pre-Attached Holder	6
Sol-Care™ Safety Blood Collection Needle with Multi-Sample Luer	6
SOL-M™ Blood Collection Tube Holder	6
Sol-Care™ Safety Blood Collection Needle and Infusion Set	6
SOL-M™ Blood Collection Tube Holder with Luer Lock Adapter	6
Blood Draw Hypodermic Needle-Pro™ Devices	7
Safe-T-Lance™ Plus Safety Lancets	7
neoheel™ Safety Lancets	7
Hypodermics	8
Needle and Syringe Combinations	8
Needle-Pro Edge™ Safety Hypodermic Needles	8
Needle-Pro Edge™ Safety Hypodermic Needle/Syringe Combinations	9
Needle-Pro™ Safety Hypodermic Needles	9
Needle-Pro™ Safety Hypodermic Needle/Syringe Combinations	10
Needle-Pro Edge™ Safety Insulin Syringes	10
Insulin Syringes with Needle-Pro Edge™ Safety Hypodermic Needles	10
Needle-Pro™ Safety TB Syringes	10
TB Syringes with Needle-Pro Edge™ Safety Hypodermic Needles	10
Blunt Fill Needle	10
Peripheral IV Catheters	11
Deltaven® Single Port	11
Deltaven® Single Port with Stopcock and Needleless Connector	11
Deltaven® Dual Port with End Caps	11
Deltaven® Dual Port with Dual Needleless Connectors	12
ProtectIV Safety™ Safety IV Catheters	12
ProtectIV Plus™ Safety IV Catheters	12
ProtectIV™ Winged Safety IV Catheters	12
ProtectIV Plus™ Winged Safety IV Catheters	13
SuperCath® 5 Safety IV Catheter — Non-Winged	13
SuperCath® 5 Safety IV Catheter — Winged	14
Jelco™ IV Catheter — Straight	14
Jelco™ IV Catheter Clear — Straight	14
Jelco™ 2 IV Catheter — Injection Valve	14
Optiva™ IV Catheter — Straight	15
Optiva™ W IV Catheter — with Fixation Wings	15
Optiva™ 2 IV Catheter — with Injection Valve	15

Ports	16
Port-A-Cath II Power P.A.C.™ Polysulfone/Titanium	16
Port-A-Cath II Power P.A.C.™ Systems with Valved Peel-Away Introducers	17
Port-A-Cath Power P.A.C.™ Titanium	17
P.A.S. Port™ T2 P.A.C. Titanium	18
Port-A-Cath II Systems with Valved Peel-Away Introducers	18
Port-A-Cath II™ Polysulfone/Titanium	19
Port-A-Cath™ Titanium Implantable Venous Access Systems	20
P.A.S. Port™ T2 Titanium Implantable Peripheral Venous Access Systems	20
P.A.S. Port Elite — Plastic/Titanium Implantable Peripheral Venous Access System	20
Non-Valved Peel-Away Introducer Sets	21
Valved Peel-Away Introducer Sets	21
Subcutaneous Tunneling Device	21
Port Access Needles	22
Gripper Plus™ Power P.A.C. Safety Huber Needles	22
Gripper Micro™ Blunt Cannula, Non-Coring Safety Needles	23
Gripper Plus™ Safety Huber Needles	24
Gripper™ Huber Needles	25
Port-A-Cath™ Plastic Hub Huber Needles	26
Subcutaneous Infusion Set	27
Cleo™ 90 Subcutaneous Infusion set	27

Arterial Blood Sampling

Pro Vent Plus 3 mL Arterial Blood Sampling Kits

Product Code	Description	Needle Gauge	Needle Length	Safety Device	Filter-Pro™	Heparin	Units/Box	Units/Case
4645PG	3 mL Luer Slip Pro-Vent Plus	22G	25 mm (1")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4698PG	3 mL Luer Slip Pro-Vent Plus	23G	25 mm (1")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4590PG	3 mL Luer Lock Pro-Vent Plus	22G	25 mm (1")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4558PG	3 mL Luer Lock Pro-Vent Plus	23G	25 mm (1")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4645PE	3 mL Luer Slip Pro-Vent Plus	22G	25 mm (1")	Needle-Pro	•	23.5 IU/ml dry lithium-balanced heparin	50	200
4698LH	3 mL Luer Slip Pro-Vent Plus	23G	25 mm (1")	Needle-Pro	•	7 IU/mL dry lithium heparin	50	200
4698PE	3 mL Luer Slip Pro-Vent Plus	23G	25 mm (1")	Needle-Pro	•	23.5 IU/ml dry lithium-balanced heparin	50	200
G1689	3 mL Luer Slip Pro-Vent Plus	25G	16 mm (5/8")	Needle-Pro	•	23.5 IU/ml dry lithium-balanced heparin	50	200
4646E	3 mL Luer Slip Pro-Vent Plus	—	—	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200
4648PE	3 mL Luer Slip Pro-Vent Plus	22G	25 mm (1")	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200
4658PE	3 mL Luer Slip Pro-Vent Plus	23G	25 mm (1")	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200
4558PE	3 mL Luer Lock Pro-Vent Plus	23G	25 mm (1")	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200

Pro Vent Plus 1 mL Arterial Blood Sampling Kits

Product Code	Description	Needle Gauge	Needle Length	Safety Device	Filter-Pro	Heparin	Units/Box	Units/Case
4629PG	1 mL Luer Slip Pro-Vent Plus	23G	25 mm (1")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4619PG	1 mL Luer Slip Pro-Vent Plus	25G	16 mm (5/8")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4629PE	1ml Luer Slip Pro-Vent Plus	23G	25 mm (1")	Needle-Pro	•	23.5 IU/ml dry lithium-balanced heparin	50	200
4619PE	1ml Luer Slip Pro-Vent Plus	25G	16 mm (5/8")	Needle-Pro	•	23.5 IU/ml dry lithium-balanced heparin	50	200
4616E	1 mL Luer Slip Pro-Vent Plus	—	—	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200
4628PE	1 mL Luer Slip Pro-Vent Plus	23G	25 mm (1")	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200
4618PE	1 mL Luer Slip Pro-Vent Plus	25G	16 mm (5/8")	Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200

Pulsator Plus Arterial Blood Sampling Kits								
Product Code	Description	Needle Gauge	Needle Length	Safety Device	Filter-Pro	Heparin	Units/Box	Units/Case
Dry Pulsator Plus 3 mL Arterial Blood Sampling Kits								
4080PG	3 mL Luer Lock Pulsator Plus	23G	25mm (1")	Needle-Pro Edge	•	23.5 IU/mL dry lithium balanced heparin	50	200
4080PE	3 mL Luer Lock Pulsator Plus	23G	25mm (1")	Needle-Pro	•	23.5 IU/mL dry lithium balanced heparin	50	200
4046E	3 mL Luer Lock Pulsator Plus	—		Disk	•	23.5 IU/mL dry lithium balanced heparin	50	200

Portex LineDraw Plus Arterial Blood Sampling Kits								
Product Code	Description	Needle	Safety Device	Filter-Pro	Heparin		Units/Box	Units/Case
4041E	1 mL Luer Slip Aspirating Syringe	—	—	•	23.5 IU/mL dry lithium balanced heparin		50	200
4042E	3 mL Luer Lock Aspirating Syringe	—	—	•	23.5 IU/mL dry lithium balanced heparin		50	200
4043E	3 mL Luer Slip Aspirating Syringe	—	—	•	23.5 IU/mL dry lithium balanced heparin		50	200
4041LH	1 mL Luer Slip Aspirating Syringe	—	—	•	7 I.U. per mL of dry lithium heparin		50	200



Pro-Vent Plus Arterial Blood Sampling Kits



Portex LineDraw Plus Arterial Blood Sampling Kits



Pulsator Arterial Blood Sampling Kits



Filter-Pro



Needle-Pro



Needle-Pro Edge



Disk

Blood Collection

	Saf-T Holder Devices and Blood Transfer Devices			
	Product Code	Description	Units/Box	Units/Case
	96000	Saf-T Holder with Multi-Sample Male Luer Adapter	50	200
	96000S	Short Saf-T Holder with Multi-Sample Male Luer Adapter	50	200
	96002	Saf-T Holder with Female Luer Adapter	50	200
	96002S	Short Saf-T Holder with Female Luer Adapter	50	200
	96004	Saf-T Holder Blood Culture Device with Male Luer Adapter	–	50
	96005	Saf-T Holder Blood Culture Device with Female Luer Adapter	–	50
	96007	Paediatric Tube Insert	10	100

All Saf-T Holder devices incorporate a multi-sample luer adapter.

	Sol Care Safety Blood Collection Needle with Pre Attached Holder						
	Product Code	Needle Gauge	Needle Length	Tubing Length	Colour Coded	Units/Box	Units/Case
	110201030002	21G	19 mm (¾")	30 cm (12")	●	25	200
	110201030006	23G	19 mm (¾")	30 cm (12")	●	25	200

	Sol Care Safety Blood Collection Needle with Multi Sample Luer						
	Product Code	Needle Gauge	Needle Length	Tubing Length	Colour Coded	Units/Box	Units/Case
	110201030030	21G	19 mm (¾")	18 cm (7")	●	50	200
	110201030031	21G	19 mm (¾")	30 cm (12")	●	50	200
	110201030032	23G	19 mm (¾")	18 cm (7")	●	50	200
	110201030033	23G	19 mm (¾")	30 cm (12")	●	50	200
	110201030034	25G	19 mm (¾")	18 cm (7")	●	50	200
	110201030035	25G	19 mm (¾")	30 cm (12")	●	50	200

	Sol M Blood Collection Tube Holder			
	Product Code	Description	Units/Box	Units/Case
	110201050001	Single-Use Blood Collection Tube Holder	250	1000

	Sol Care Safety Blood Collection Needle and Infusion Set						
	Product Code	Needle Gauge	Needle Length	Tubing Length	Colour Coded	Units/Box	Units/Case
	213412SI	21G	19mm (¾")	30 cm (12")	●	50	200
	233412SI	23G	19mm (¾")	30 cm (12")	●	50	200
	253412SI	25G	19mm (¾")	30 cm (12")	●	50	200

	Sol M Blood Collection Tube Holder with Luer Lock Adapter			
	Product Code	Description	Units/Box	Units/Case
	110201050003	Multi-Sample Tube Holder with Back-End Needle and Luer Connector	50	800

	Blood Draw Hypodermic Needle Pro Devices					
	Product Code	Needle Gauge	Needle Length	Colour Code	Units Per Box	Units Per Case
	4271	20G	25 mm (1")	●	100	800
	4270	20G	38 mm (1 ½")	●	100	800
	4273	21G	25 mm (1")	●	100	800
	4272	21G	38 mm (1 ½")	●	100	800
	4275	22G	25 mm (1")	●	100	800
	4274	22G	38 mm (1 ½")	●	100	800
	4276	23G	25 mm (1")	●	100	800
	4277	25G	16 mm (5/8")	●	100	800
	4278	25G	25 mm (1")	●	100	800

	Safe T Lance Plus Safety Lancets						
	Product Code	Description	Blade Width / Needle Gauge	Depth	Colour Code	Units/Box	Units/Case
	1028-1	Paediatric	1.5 mm Blade	1.2 mm	●	200	2000
	1005-1	Max Flow	1.5 mm Blade	1.6 mm	●	200	2000
	1008-1	High Flow	18G Needle	1.8 mm	●	200	2000
	1014-1	Normal Flow	21G Needle	1.8 mm	●	200	2000
	1020-1	Low Flow	25G Needle	1.4 mm	●	200	2000
	1025-1	Micro Flow	28G Needle	1.6 mm	●	200	2000

	neoheel Safety lancets							
	Product Code	Description	Incision depth	Incision length	Blood volume	Colour Code	Units/Box	Units/Case
	1050M	Micro-Premie Low birth weight <1kg	0.65 mm	1.40 mm	Medium	●	50	200
	1051P	Premie Low birth weight 1-2.5 kg	0.85 mm	1.75 mm	Medium	●	50	200
	1052N	Newborn Birth-6 months 2.5-9 kg	1.00 mm	2.50 mm	Medium	●	50	200
	1053T	Toddler 6-24 months >9 kg	2.00 mm	3.00 mm	High	●	50	200

Hypodermics

Needles and Needle/Syringe Combinations

	Needle Pro Edge Safety Hypodermic Needles					
	Product Code	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	401810	18G	25 mm (1")	●	100	1000
	401815	18G	38 mm (1 1/2")	●	100	1000
	401910	19G	25 mm (1")	●	100	1000
	401915	19G	38 mm (1 1/2")	●	100	1000
	402010	20G	25 mm (1")	●	100	1000
	402015	20G	38 mm (1 1/2")	●	100	1000
	402110	21G	25 mm (1")	●	100	1000
	402115	21G	38 mm (1 1/2")	●	100	1000
	402210	22G	25 mm (1")	●	100	1000
	402215	22G	38 mm (1 1/2")	●	100	1000
	402310	23G	25 mm (1")	●	100	1000
	402315	23G	38 mm (1 1/2")	●	100	1000
	402558	25G	16 mm (5/8")	●	100	1000
	402510	25G	25 mm (1")	●	100	1000
	402515	25G	38 mm (1 1/2")	●	100	1000
	402605	26G	13 mm (1/2")	●	100	1000
	402705	27G	13 mm (1/2")	●	100	1000
	4027125	27G	32 mm (1 1/4")	●	100	1000
	403005	30G	13 mm (1/2")	●	100	1000

	Needle Pro Edge Safety Hypodermic Needle/Syringe Combinations						
	Product Code	Description	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	431810	3 mL Luer Lock	18G	25 mm (1")	●	50	400
	431815	3 mL Luer Lock	18G	38 mm (1 ½")	●	50	400
	431910	3 mL Luer Lock	19G	25 mm (1")	●	50	400
	431915	3 mL Luer Lock	19G	38 mm (1 ½")	●	50	400
	432010	3 mL Luer Lock	20G	25 mm (1")	●	50	400
	432015	3 mL Luer Lock	20G	38 mm (1 ½")	●	50	400
	432110	3 mL Luer Lock	21G	25 mm (1")	●	50	400
	432115	3 mL Luer Lock	21G	38 mm (1 ½")	●	50	400
	432210	3 mL Luer Lock	22G	25 mm (1")	●	50	400
	432215	3 mL Luer Lock	22G	38 mm (1 ½")	●	50	400
	432310	3 mL Luer Lock	23G	25 mm (1")	●	50	400
	432510	3 mL Luer Lock	25G	25 mm (1")	●	50	400
	432558	3 mL Luer Lock	25G	16 mm (⅝")	●	50	400
	452010	5 mL Luer Lock	20G	25 mm (1")	●	50	400
	452015	5 mL Luer Lock	20G	38 mm (1 ½")	●	50	400
	452110	5 mL Luer Lock	21G	25 mm (1")	●	50	400
	452115	5 mL Luer Lock	21G	38 mm (1 ½")	●	50	400
	452215	5 mL Luer Lock	22G	38 mm (1 ½")	●	50	400
	4102010	10 mL Luer Lock	20G	25 mm (1")	●	50	400
	4102015	10 mL Luer Lock	20G	38 mm (1 ½")	●	50	400
	4102115	10 mL Luer Lock	21G	38 mm (1 ½")	●	50	400

	Needle Pro Safety Hypodermic Needles					
	Product Code	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	4281	18G	25 mm (1")	●	100	800
	4280	18G	38 mm (1 ½")	●	100	800
	4283	19G	25 mm (1")	●	100	800
	4282	19G	38 mm (1 ½")	●	100	800
	4285	20G	25 mm (1")	●	100	800
	4284	20G	38 mm (1 ½")	●	100	800
	4287	21G	25 mm (1")	●	100	800
	4286	21G	38 mm (1 ½")	●	100	800
	4289	22G	25 mm (1")	●	100	800
	4288	22G	38 mm (1 ½")	●	100	800
	4290	23G	25 mm (1")	●	100	800
	4291	25G	16 mm (⅝")	●	100	800
	4292	25G	25 mm (1")	●	100	800
	4293	25G	38 mm (1 ½")	●	100	800

	Needle Pro Safety Hypodermic Needle/Syringe Combination (3 mL Luer Lock)					
	Product Code	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	4231	20G	25 mm (1")	●	50	400
	4230	20G	38 mm (1 1/2")	●	50	400
	4233	21G	25 mm (1")	●	50	400
	4232	21G	38 mm (1 1/2")	●	50	400
	4235	22G	25 mm (1")	●	50	400
	4234	22G	38 mm (1 1/2")	●	50	400
	4236	23G	25 mm (1")	●	50	400
	4237	25G	16 mm (5/8")	●	50	400
	4238	25G	25 mm (1")	●	50	400

	Needle Pro Safety Insulin Syringes (permanently attached)						
	Product Code	Description	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	43125-03	0.3mL U-100 insulin syringe	31G	6 mm (1/4")	●	100	600
	42950-03	0.3mL U-100 insulin syringe	29G	12.5 mm (1/2")	●	100	600
	43125-05	0.5mL U-100 insulin syringe	31G	6mm (1/4")	●	100	600
	42950-05	0.5mL U-100 insulin syringe	29G	12.5 mm (1/2")	●	100	600
	43125-1	1mL U-100 insulin syringe	31G	6mm (1/4")	●	100	600
	42950-1	1mL U-100 insulin syringe	29G	12.5 mm (1/2")	●	100	600

	Insulin Syringes with Needle Pro Edge Safety Hypodermic Needles						
	Product Code	Description	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	4317	1 mL Luer Slip U-100 Insulin	26G	12.5 mm (1/2")	●	50	400
	4319	1 mL Luer Slip U-100 Insulin	27G	12.5 mm (1/2")	●	50	400

	Needle Pro Safety TB Syringes (permanently attached)						
	Product Code	Description	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	42638-TB	1mL TB syringe	26G	10 mm (3/8")	●	100	600
	42750-TB	1mL TB syringe	27G	12.5 mm (1/2")	●	100	600

	TB Syringes with Needle Pro Edge Safety Hypodermic Needles						
	Product Code	Description	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	412558	1 mL Luer Lock	25G	16 mm (5/8")	●	50	400
	412705	1 mL Luer Lock	27G	12.5 mm (1/2")	●	50	400

	Blunt Fill Needle						
	Product Code	Description	Needle Gauge	Needle Length	Colour Code	Units/Box	Units/Case
	BN1815	Luer Lock	18G	38 mm (1 1/2")	●	100	1000
	BN1815F	Luer Lock Filtered Needle	18G	38 mm (1 1/2")	●	100	1000

Peripheral IV Catheters

	Deltaven Single Port							
	Product Code	Needle Gauge	Catheter Length	Flow (mL/min)	PSI Rating	Colour Code	Units/Box	Units/Case
	3856222	16	32 mm (1 ¼")	125	330	●	25	100
	3845222	18	25 mm (1")	65	330	●	25	100
	3846222	18	32 mm (1 ¼")	63	330	●	25	100
	3847222	18	45 mm (1 ¾")	58	330	●	25	100
	3835222	20	25 mm (1")	47	330	●	25	100
	3836222	20	32 mm (1 ¼")	42	330	●	25	100
	3837222	20	45 mm (1 ¾")	38	330	●	25	100
	3824222	22	19 mm (¾")	37	330	●	25	100
	3825222	22	25 mm (1")	31	330	●	25	100
	3814222	24	19 mm (¾")	24	330	●	25	100
	3804222	26	19 mm (¾")	15	N/A	●	25	100

	Deltaven Single Port with Stopcock and Needleless Connector							
	Product Code	Gauge	Catheter Length	Flow (mL/min)	PSI Rating	Colour Code	Units/Box	Units/Case
	3856422	16	32 mm (1 ¼")	125	330	●	20	80
	3845422	18	25 mm (1")	65	330	●	20	80
	3846422	18	32 mm (1 ¼")	63	330	●	20	80
	3847422	18	45 mm (1 ¾")	58	330	●	20	80
	3835422	20	25 mm (1")	47	330	●	20	80
	3836422	20	32 mm (1 ¼")	42	330	●	20	80
	3837422	20	45 mm (1 ¾")	38	330	●	20	80
	3824422	22	19 mm (¾")	37	330	●	20	80
	3825422	22	25 mm (1")	31	330	●	20	80
	3814422	24	19 mm (¾")	24	330	●	20	80
	3804422	26	19 mm (¾")	15	N/A	●	20	80

	Deltaven Dual Port with End Caps							
	Product Code	Gauge	Catheter Length	Flow (mL/min)	PSI Rating	Colour Code	Units/Box	Units/Case
	3852022	16	32 mm (1 ¼")	125	330	●	25	100
	3841022	18	25 mm (1")	65	330	●	25	100
	3842022	18	32 mm (1 ¼")	63	330	●	25	100
	3843022	18	45 mm (1 ¾")	58	330	●	25	100
	3831022	20	25 mm (1")	47	330	●	25	100
	3832022	20	32 mm (1 ¼")	42	330	●	25	100
	3833022	20	45 mm (1 ¾")	38	330	●	25	100
	3820022	22	19 mm (¾")	37	330	●	25	100
	3821022	22	25 mm (1")	31	330	●	25	100
	3810022	24	19 mm (¾")	24	330	●	25	100
	3800022	26	19 mm (¾")	15	N/A	●	25	100

* These SKUs do not have CE marking.



	Deltaven Dual Port with Dual Needleless Connectors							
	Product Code	Gauge	Catheter Length	Flow (mL/min)	PSI Rating	Colour Code	Units/Box	Units/Case
	3852522	16	32 mm (1 1/4")	125	330	●	20	80
	3841522	18	25 mm (1")	65	330	●	20	80
	3842522	18	32 mm (1 1/4")	63	330	●	20	80
	3843522	18	45 mm (1 3/4")	58	330	●	20	80
	3831522	20	25 mm (1")	47	330	●	20	80
	3832522	20	32 mm (1 1/4")	42	330	●	20	80
	3833522	20	45 mm (1 3/4")	38	330	●	20	80
	3820522	22	19 mm (3/4")	37	330	●	20	80
	3821522	22	25 mm (1")	31	330	●	20	80
	3810522	24	19 mm (3/4")	24	330	●	20	80
	3800522	26	19 mm (3/4")	15	N/A	●	20	80

	ProtectIV Safety IV Catheters (Radiopaque/Fluorinated Ethylene Propylene Polymer)						
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case
	3048	14G	32 mm (1 ¼")	325	●	50	200
	3042	16G	32 mm (1 ¼")	215	●	50	200
	3055	18G	32 mm (1 ¼")	110	●	50	200
	3057	20G	25 mm (1")	65	●	50	200
	3056	20G	32 mm (1 ¼")	63	●	50	200
	3050	22G	25 mm (1")	38	●	50	200
	3053	24G	19 mm (¾")	24	●	50	200

	ProtectIV Plus Safety IV Catheters (Radiopaque/Polyurethane)						
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case
	3068	14G	32 mm (1 ¼")	325	●	50	200
	3062	16G	32 mm (1 ¼")	215	●	50	200
	3065	18G	32 mm (1 ¼")	110	●	50	200
	3067	20G	25 mm (1")	65	●	50	200
	3066	20G	32 mm (1 ¼")	63	●	50	200
	3060	22G	25 mm (1")	38	●	50	200
	3063	24G	19 mm (¾")	24	●	50	200

	ProtectIV Winged Safety IV Catheters (Radiopaque/Fluorinated Ethylene Propylene Polymer)						
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case
	3078	14G	32 mm (1 ¼")	325	●	50	200
	3072	16G	32 mm (1 ¼")	215	●	50	200
	3075	18G	32 mm (1 ¼")	110	●	50	200
	3077	20G	25 mm (1")	65	●	50	200
	3076	20G	32 mm (1 ¼")	63	●	50	200
	3070	22G	25 mm (1")	38	●	50	200
	3073	24G	15 mm (⅝")	24	●	50	200

	ProtectIV Plus Winged Safety IV Catheters (Radiopaque/Polyurethane)						
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case
	3088	14G	32 mm (1 ¼")	325	●	50	200
	3082	16G	32 mm (1 ¼")	215	●	50	200
	3085	18G	32 mm (1 ¼")	110	●	50	200
	3087	20G	25 mm (1")	65	●	50	200
	3086	20G	32 mm (1 ¼")	63	●	50	200
	3080	22G	25 mm (1")	38	●	50	200
	3083	24G	15 mm (⅝")	24	●	50	200

	SuperCath 5 Safety IV Catheter — Non Winged						
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case
	SP125-26-14T*	26G	14 mm (½")	14	●	50	200
	SP125-26-19T*	26G	19 mm (¾")	12	●	50	200
	SP125-24-14T*	24G	14 mm (½")	21	●	50	200
	SP125-24-19T*	24G	19 mm (¾")	16	●	50	200
	SP120-24-19T						
	SP120-22-25T	22G	25 mm (1")	31	●	50	200
	SP120-22-31T	22G	31 mm (1 ¼")	28	●	50	200
	SP120-20-25T	20G	25 mm (1")	62	●	50	200
	SP120-20-31T	20G	31 mm (1 ¼")	55	●	50	200
	SP120-20-50T	20G	50 mm (2")	52	●	50	200
	SP120-18-31T	18G	31 mm (1 ¼")	100	●	50	200
	SP120-18-50T	18G	50 mm (2")	90	●	50	200
	SP120-16-31T	16G	31 mm (1 ¼")	183	●	50	200
	SP125-16-31T*						
	SP125-14-31T*	14G	31 mm (1 ¼")	272	●	50	200

*No blood-containment valve

	SuperCath 5 Safety IV Catheter — Winged						
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case
	SP125-26-14-WT*	26G	14 mm (½")	14	●	50	200
	SP-125-26-19WT*	26G	19 mm (¾")	12	●	50	200
	SP125-24-14-WT*	24G	14 mm (½")	21	●	50	200
	SP125-24-19-WT*	24G	19 mm (¾")	16	●	50	200
	SP120-24-19-WT						
	SP120-22-25-WT	22G	25 mm (1")	31	●	50	200
	SP120-22-31-WT	22G	31 mm (1 ¼")	28	●	50	200
	SP120-20-25-WT	20G	25 mm (1")	62	●	50	200
	SP120-20-31-WT	20G	31 mm (1 ¼")	55	●	50	200
	SP120-20-50-WT	20G	50 mm (2")	52	●	50	200
	SP120-18-31-WT	18G	31 mm (1 ¼")	100	●	50	200
	SP120-18-50-WT	18G	50 mm (2")	90	●	50	200
	SP120-16-31-WT	16G	31 mm (1 ¼")	183	●	50	200
	SP125-16-31-WT*						
	SP125-14-31-WT*	14G	31 mm (1 ¼")	272	●	50	200

*No blood-containment valve

	Jelco IV Catheter (FEP radiopaque) – Straight							
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case	Obturator Codes
	4038	14G	50 mm (2")	315	●	50	600	4215
	4032	16G	50 mm (2")	210	●	50	600	4217
	4035	18G	32 mm (1 ¼")	110	●	50	600	4219
	4034	18G	45 mm (1 ¾")	100	●	50	600	4218
	4036	20G	32 mm (1 ¼")	65	●	50	600	4220
	4039	20G	45 mm (1 ¾")	60	●	50	600	–
	4030	22G	25 mm (1")	38	●	50	600	4222
	4033	24G	19 mm (¾")	24	●	50	600	4224

	Jelco IV Catheter Clear (FEP non radiopaque) – Straight							
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case	Obturator Codes
	4018	14G	50 mm (2")	315	●	50	600	4215
	4012	16G	50 mm (2")	210	●	50	600	4217
	4014	18G	45 mm (1 ¾")	100	●	50	600	4218
	4016	20G	32 mm (1 ¼")	65	●	50	600	4220
	4010	22G	25 mm (1")	38	●	50	600	4222
	4013	24G	19 mm (¾")	24	●	50	600	4224

	Jelco 2 IV Catheter (FEP radiopaque) – Injection Valve							
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case	Obturator Codes
	MR14	14G	45 mm (1 ¾")	315	●	50	600	4514
	MR16	16G	45 mm (1 ¾")	210	●	50	600	4516
	MR18	18G	45 mm (1 ¾")	100	●	50	600	4518
	MR19	18G	32 mm (1 ¼")	110	●	50	600	4519
	MR20	20G	32 mm (1 ¼")	65	●	50	600	4520
	MR22	22G	25 mm (1")	38	●	50	600	4522

	Optiva IV Catheter (Polyurethane) – Straight							
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case	Obturator Codes
	5068	14G	50 mm (2")	315	●	50	600	4215
	5061	16G	32 mm (1 ¼")	225	●	50	600	–
	5062	16G	50 mm (2")	210	●	50	600	4217
	5065	18G	32 mm (1 ¼")	110	●	50	600	4219
	5064	18G	45 mm (1 ¾")	100	●	50	600	4218
	5066	20G	32 mm (1 ¼")	65	●	50	600	4220
	5069	20G	45 mm (1 ¾")	60	●	50	600	–
	5060	22G	25 mm (1")	38	●	50	600	4222
	5063	24G	19 mm (¾")	24	●	50	600	4224

	Optiva W IV Catheter (Polyurethane) – with Fixation Wings							
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case	Obturator Codes
	1114	14G	45 mm (1 ¾")	310	●	50	600	4514
	1116	16G	45 mm (1 ¾")	200	●	50	600	4516
	1119	18G	32 mm (1 ¼")	95	●	50	600	4519
	1118	18G	45 mm (1 ¾")	95	●	50	600	4518
	1120	20G	32 mm (1 ¼")	60	●	50	600	4520
	1122	22G	25 mm (1")	35	●	50	600	4522
	1124	24G	19 mm (¾")	20	●	50	600	4514

	Optiva 2 IV Catheter (Polyurethane) – with Injection Valve							
	Product Code	Needle Gauge	Needle Length	Flow (mL/min)	Colour Code	Units/Box	Units/Case	Obturator Codes
	1014	14G	45 mm (1 ¾")	315	●	50	600	4514
	1016	16G	45 mm (1 ¾")	210	●	50	600	4516
	1019	18G	32 mm (1 ¼")	110	●	50	600	4519
	1018	18G	45 mm (1 ¾")	100	●	50	600	4518
	1020	20G	32 mm (1 ¼")	65	●	50	600	4520
	1022	22G	25 mm (1")	38	●	50	600	4522

Port-A-Cath II Power P.A.C. Polysulfone/Titanium Implantable Venous Access Systems

MR Conditional up to 3 Tesla, power injection up to 5 mL/sec and 300 psi

	Polysulfone/Titanium Portals with PolyFlow Polyurethane Catheters								
		Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
	Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4452-24 Kit 21-4453-24 Tray	1 0610586 03267 7 1 0610586 03268 4	each each
	Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4454-24 Kit 21-4455-24 Tray	1 0610586 03269 1 1 0610586 03270 7	each each
	Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4482-24 Kit 21-4483-24 Tray	1 0610586 03276 9 1 0610586 03277 6	each each
	Single-lumen Low Profile Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4470-24 Kit 21-4471-24 Tray	1 0610586 03273 8 1 0610586 03274 5	each each
	Dual-lumen Low Profile Portal	10 Fr	Polyurethane	3.2	1.4/1.4	76	21-8468-24 Tray	1 0610586 03285 1	each
	Dual-lumen Low Profile Portal	7 Fr	Polyurethane	2.2	1.0/1.0	76	21-8466-24 Tray	1 0610586 03283 7	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick and Point-Lok™ device.

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus Power P.A.C. safety needle, 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

Port-A-Cath II Power P.A.C. Systems with Valved Peel-Away Introducers

Introducer sheath length – 13 cm Dilator sheath length – 18 cm

Port A Cath Power P.A.C. Systems with Valved Peel Away Introducers									
	Introducer Size	Catheter Material	OD (mm)	Length (cm)	Connector	Product Code	GTIN	Unit	
	Single-lumen Portal	8.5 Fr	Polyurethane	2.6	76	Ultra-Lock™ connector	21-4474-24 Tray	1 0610586 03996 6	each
	Single-lumen Low Profile Portal	8.5 Fr	Polyurethane	2.6	76	Ultra-Lock connector	21-4475-24 Tray	1 0610586 03997 3	each
	Single-lumen Portal	6 Fr	Polyurethane	1.9	76	Ultra-Lock connector	21-4476-24 Tray	1 0610586 03998 0	each
	Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	76	Ultra-Lock connector	21-4477-24 Tray	1 0610586 03999 7	each
	Dual-lumen Low Profile Portal	7 Fr	Polyurethane	2.2	76	Wing-Lock connector	21-8469-24 Tray	1 0610586 04000 9	each
	Dual-lumen Low Profile Portal	10 Fr	Polyurethane	3.2	76	Wing-Lock connector	21-8470-24 Tray	1 0610586 04001 6	each

Introducers are available in ICU Medical's Port-A-Cath system trays and also as separate sets.

Port-A-Cath Power P.A.C. Titanium Implantable Venous Access Systems

MR Conditional up to 3 Tesla, power injection up to 5 mL/sec and 300 psi

Titanium Contoured Portals with PolyFlow Polyurethane Catheters									
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit	
	Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4422-24 Kit	1 0610586 03261 5	each
							21-4423-24 Tray	1 0610586 03262 2	each
	Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4424-24 Kit	1 0610586 03263 9	each
							21-4425-24 Tray	1 0610586 03264 6	each
	Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4436-24 Kit	1 0610586 03265 3	each
							21-4437-24 Tray	1 0610586 03266 0	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick and Point-Lok device.

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus Power P.A.C. safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

P.A.S. Port T2 Power P.A.C. Titanium Implantable Peripheral Venous Access Systems

MR Conditional up to 3 Tesla, power injection up to 5 mL/sec and 300 psi

	Titanium Contoured Portals with PolyFlow Polyurethane Catheters								
		Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
	Peripheral Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4872-24 Kit	1 0610586 03280 6	each
21-4873-24 Tray							1 0610586 03281 3	each	

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, and Point-Lok device.

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus Power P.A.C. safety needle(s), 2 syringes, tunneling tool, and Point-Lok device(s).

All systems are unassembled.

Port-A-Cath II Systems with Valved Peel-Away Introducers

	Port A Cath II Systems with Valved Peel Away Introducers Introducer sheath length – 13 cm Dilator sheath length – 18 cm								
		Introducer Size	Catheter Material	OD (mm)	Length (cm)	Connector	Product Code	GTIN	Unit
	Single-lumen Portal	8.5 Fr	Polyurethane	2.6	76	Ultra-Lock connector	21-4066-24 Tray	1 0610586 03990 4	each
	Single-lumen Portal	6 Fr	Polyurethane	1.9	76	Ultra-Lock connector	21-4067-24 Tray	1 0610586 03991 1	each
	Single-lumen Low Profile Portal	8.5 Fr	Polyurethane	2.6	76	Ultra-Lock connector	21-4068-24 Tray	1 0610586 03992 8	each
	Single-lumen Low Profile Portal	6 Fr	Polyurethane	1.9	76	Ultra-Lock connector	21-4069-24 Tray	1 0610586 03993 5	each

Introducers are available in ICU Medical's Port-A-Cath system trays and also as separate sets.

Port-A-Cath II Polysulfone/Titanium Implantable Venous Access Systems

MR Conditional up to 3 Tesla

	Polysulfone/Titanium Portals with PolyFlow Polyurethane Catheters								
		Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
	Single-lumen Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4052-24 Kit	1 0610586 02352 1	each
							21-4053-24 Tray	1 0610586 02353 8	each
	Single-lumen Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4054-24 Kit	1 0610586 02354 5	each
							21-4055-24 Tray	1 0610586 02355 2	each
	Single-lumen Low Profile™ Portal	6 Fr	Polyurethane	1.9	1.0	76	21-4082-24 Kit	1 0610586 02385 9	each
							21-4083-24 Tray	1 0610586 02361 3	each
	Single-lumen Low Profile™ Portal	8.5 Fr	Polyurethane	2.6	1.6	76	21-4070-24 Kit	1 0610586 01782 7	each
							21-4071-24 Tray	1 0610586 01784 1	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick and Point-Lok device(s).

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, “J” guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle, 2 syringes, tunneling tool and Point-Lok device.

All systems are unassembled.

Polysulfone/Titanium Portals with Silicone Catheters									
		Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
Single-lumen Portal		9 Fr	Silicone	2.8	1.0	76	21-4050-24 Kit	1 0610586 02350 7	each
							21-4051-24 Tray	1 0610586 02351 4	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick and Point-Lok device(s).

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, “J” guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle, 2 syringes, tunneling tool and Point-Lok device.

All systems are unassembled.

Port-A-Cath Titanium Implantable Venous Access Systems

MR Conditional up to 3 Tesla

	Titanium Contoured Portals with PolyFlow Polyurethane Catheters							
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
	Single-lumen Portal 6 Fr	Polyurethane	1.9	1.0	76	21-4024-24 Kit	1 0610586 02379 8	each
						21-4025-24 Tray	1 0610586 02380 4	each
	Single-lumen Portal 8.5 Fr	Polyurethane	2.6	1.6	76	21-4022-24 Kit	1 0610586 02377 4	each
						21-4023-24 Tray	1 0610586 02378 1	each
	Single-lumen Low Profile Portal 6 Fr	Polyurethane	1.9	1.0	76	21-4036-24 Kit	1 0610586 02383 5	each
						21-4037-24 Tray	1 0610586 02384 2	each

Titanium Contoured Portals with Silicone Catheters								
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
Single-lumen Portal 9 Fr	9 Fr	Silicone	2.8	1.0	76	21-4008-24 Kit	1 0610586 02374 3	each
						21-4009-24 Tray	1 0610586 02375 0	each
Single-lumen Low Profile Portal 9 Fr	9 Fr	Silicone	2.8	1.0	76	21-4034-24 Kit	1 0610586 02381 1	each
						21-4035-24 Tray	1 0610586 02382 8	each

Kits include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick and Point-Lok device(s).

Trays include: Portal, catheter, Port-A-Cath straight needle, blunt needle, vein pick, 18G extra thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

All systems are unassembled.

P.A.S. Port T2 Titanium Implantable Peripheral Venous Access Systems

MR Conditional up to 3 Tesla

	Titanium Portal with Polyurethane Catheter							
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Product Code	GTIN	Unit
	P.A.S. Port T2 Titanium Portal 6 Fr	Polyurethane	1.9	1.0	76	21-4572-24 Kit	1 0610586 02516 7	each
						21-4573-24 Tray	1 0610586 02517 4	each

P.A.S. Port Elite – Plastic/Titanium Implantable Peripheral Venous Access System

	Plastic/ Titanium Portal with Polyurethane Catheter							
	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (cm)	Reorder No.	GTIN	Unit
	Single-lumen Portal 6 Fr	Polyurethane	1.9	1.0	76	21-4590-24 Kit	1 0610586 02154 1	each
						21-4591-24 Tray	1 0610586 02156 5	each

Port-A-Cath System Accessories

	Non Valved Peel Away Introducer Sets					
		Sheath Length (cm)	Dilator Length (cm)	Product Code	GTIN	Unit
	6 French Includes two guidewires	16	19	21-2320-24	3 0610586 02511 6	5/box
	7 French	16	19	21-2471-24*	3 0610586 04015 7	5/box
	8.5 French	16	19	21-2470-24	3 0610586 02514 7	5/box
	9 French	16	19	21-2340-24	3 0610586 02512 3	5/box
	10 French	16	19	21-2472-24*	3 0610586 04016 4	5/box
	11 French	16	19	21-2386-24	3 0610586 02513 0	5/box

Non-valved peel-away introducer sets are available in ICU Medical's Port-A-Cath system trays and also as separate sets.

*Subject to availability

	Valved Peel Away Introducer Sets					
		Sheath Length (cm)	Dilator Length (cm)	Product Code	GTIN	Unit
	6 French	13	18	21-2476-24*	3 0610586 04020 1	5/box
	7 French	13	18	21-2477-24*	3 0610586 04021 8	5/box
	8.5 French	13	18	21-2478-24*	3 0610586 04022 5	5/box
	9 French	13	18	21-2479-24*	3 0610586 04023 2	5/box
	10 French	13	18	21-2480-24*	3 0610586 04024 9	5/box
	11 French	13	18	21-2481-24*	3 0610586 04025 6	5/box

*Subject to availability

	Subcutaneous Tunneling Device		
	Product Code	GTIN	Unit
3.2 mm OD x 24 cm length, solid malleable stainless steel, disposable	21-0125-24	3 0610586 01319 9	5/box

Introducer sets include: 18G extra thin wall introducer needle, "J" guidewire, vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle (except 21-2320), 2 syringes and Point-Lok device.

Gripper Plus Power P.A.C. Safety Huber Needles

MR Conditional up to 3 Tesla, power injection up to 5 mL/sec and 300 psi

	Gripper Plus Power P.A.C. Safety Huber Needles Not made with natural rubber latex, non coring, non siliconised, removable contoured grip, cushioned needle platform, 20 cm extension tubing* with or without Y site				
	Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	GTIN	Unit
	● 20 gauge	16 mm	21-3464-24	3 0610586 03363 6	12/box
		19 mm	21-3465-24	3 0610586 03304 3	12/box
		25 mm	21-3466-24	3 0610586 03305 0	12/box
		32 mm	21-3467-24	3 0610586 03306 7	12/box
	● 19 gauge	19 mm	21-3468-24	3 0610586 03307 4	12/box
		25 mm	21-3469-24	3 0610586 03308 1	12/box
		32 mm	21-3470-24	3 0610586 03309 8	12/box
	Without Y-site	Length	Reorder No.	GTIN	Unit
	● 20 gauge	16 mm	21-3369-24	3 0610586 03302 9	12/box
		19 mm	21-3367-24	3 0610586 03300 5	12/box
		25 mm	21-3362-24	3 0610586 03296 1	12/box
		32 mm	21-3363-24	3 0610586 03297 8	12/box
	● 19 gauge	19 mm	21-3368-24	3 0610586 03301 2	12/box
		25 mm	21-3364-24	3 0610586 03298 5	12/box
		32 mm	21-3365-24	3 0610586 03299 2	12/box

* Manufactured with non-DEHP plasticised tubing

Gripper Micro Blunt Cannula, Non-Coring Safety Needles

MR Conditional up to 3 Tesla

Gripper Micro Blunt Cannula, Non-Coring Safety Needles Not made with natural rubber latex, non coring, non siliconised, cushioned infusion site, 20 cm extension tubing* with or without Y site					
	Needleless Y-site Luer-activated needleless injection site	Length	Reorder No.	GTIN	Unit
	● 22 gauge	19 mm	21-3271-24	3 0610586 02988 6	12/box
		25 mm	21-3272-24	3 0610586 02989 3	12/box
		32 mm	21-3273-24	3 0610586 02990 9	12/box
	● 20 gauge	19 mm	21-3276-24	3 0610586 02993 0	12/box
		25 mm	21-3277-24	3 0610586 02994 7	12/box
		32 mm	21-3278-24	3 0610586 02995 4	12/box
	● 19 gauge	19 mm	21-3280-24	3 0610586 02997 8	12/box
		25 mm	21-3281-24	3 0610586 02998 5	12/box
		32 mm	21-3282-24	3 0610586 02999 2	12/box
	Without Y-site	Length	Reorder No.	GTIN	Unit
	● 22 gauge	19 mm	21-3251-24	3 0610586 02974 9	12/box
		25 mm	21-3252-24	3 0610586 02975 6	12/box
		32 mm	21-3253-24	3 0610586 02976 3	12/box
	● 20 gauge	19 mm	21-3256-24	3 0610586 02979 4	12/box
		25 mm	21-3257-24	3 0610586 02980 0	12/box
		32 mm	21-3258-24	3 0610586 02981 7	12/box
	● 19 gauge	19 mm	21-3260-24	3 0610586 02983 1	12/box
		25 mm	21-3261-24	3 0610586 02984 8	12/box
		32 mm	21-3262-24	3 0610586 02985 5	12/box

* Manufactured with non-DEHP plasticised tubing

Gripper Plus Safety Huber Needles

MR Conditional up to 3 Tesla

Gripper Plus Safety Huber Needles Not made with natural rubber latex, non coring, non siliconised, removable contoured grip, cushioned needle platform, 20 cm extension tubing* with or without Y site				
Needleless Y-site luer-activated needleless injection site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2960-24	3 0610586 02768 4	12/box
	19 mm	21-2961-24	3 0610586 02769 1	12/box
	25 mm	21-2962-24	3 0610586 02770 7	12/box
	32 mm	21-2963-24	3 0610586 02771 4	12/box
● 20 gauge	16 mm	21-2964-24	3 0610586 02772 1	12/box
	19 mm	21-2965-24	3 0610586 02773 8	12/box
	25 mm	21-2966-24	3 0610586 02774 5	12/box
	32 mm	21-2967-24	3 0610586 02775 2	12/box
● 19 gauge	19mm	21-2968-24	3 0610586 02776 9	12/box
	25 mm	21-2969-24	3 0610586 02777 6	12/box
	32 mm	21-2970-24	3 0610586 02778 3	12/box
Y-site removable split septum injection cap on luer-lock Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2860-24	3 0610586 02564 2	12/box
	19 mm	21-2861-24	3 0610586 02565 9	12/box
	25 mm	21-2862-24	3 0610586 02566 6	12/box
	32 mm	21-2863-24	3 0610586 02567 3	12/box
● 20 gauge	16 mm	21-2864-24	3 0610586 02568 0	12/box
	19 mm	21-2865-24	3 0610586 02569 7	12/box
	25 mm	21-2866-24	3 0610586 02570 3	12/box
	32 mm	21-2867-24	3 0610586 02571 0	12/box
● 19 gauge	19 mm	21-2868-24	3 0610586 02572 7	12/box
	25 mm	21-2869-24	3 0610586 02573 4	12/box
	32 mm	21-2870-24	3 0610586 02574 1	12/box
Without Y-site	Length	Reorder No.	GTIN	Unit
● 22 gauge	16 mm	21-2770-24	3 0610586 02563 5	12/box
	19 mm	21-2766-24	3 0610586 02559 8	12/box
	25 mm	21-2760-24	3 0610586 02553 6	12/box
	32 mm	21-2761-24	3 0610586 02554 3	12/box
● 20 gauge	16 mm	21-2769-24	3 0610586 02562 8	12/box
	19 mm	21-2767-24	3 0610586 02560 4	12/box
	25 mm	21-2762-24	3 0610586 02555 0	12/box
	32 mm	21-2763-24	3 0610586 02556 7	12/box
● 19 gauge	19 mm	21-2768-24	3 0610586 02561 1	12/box
	25 mm	21-2764-24	3 0610586 02557 4	12/box
	32 mm	21-2765-24	3 0610586 02558 1	12/box

* Manufactured with non-DEHP plasticised tubing



Gripper Huber Needles

MR Conditional up to 3 Tesla

Gripper Huber Needles Not made with natural rubber latex, non coring, non siliconised, removable contoured grip, cushioned needle platform, 20 cm extension tubing* with or without Y site					
Y-site Removable injection cap on luer-lock Y-site	Length	Reorder No.	GTIN	Unit	
● 22 gauge	16 mm	21-2938-24	3 0610586 02321 1	12/box	
	19 mm	21-2939-24	3 0610586 01674 9	12/box	
	25 mm	21-2940-24	3 0610586 01675 6	12/box	
	32 mm	21-2941-24	3 0610586 01676 3	12/box	
● 20 gauge	16 mm	21-2946-24	3 0610586 02323 5	12/box	
	19 mm	21-2947-24	3 0610586 01677 0	12/box	
	25 mm	21-2948-24	3 0610586 01678 7	12/box	
	32 mm	21-2949-24	3 0610586 01679 4	12/box	
● 19 gauge	19 mm	21-2955-24	3 0610586 01680 0	12/box	
	25 mm	21-2956-24	3 0610586 01681 7	12/box	
	32 mm	21-2957-24	3 0610586 01682 4	12/box	
Without Y-site	Length	Reorder No.	GTIN	Unit	
● 22 gauge	16 mm	21-2737-24	3 0610586 02320 4	12/box	
	19 mm	21-2733-24	3 0610586 01584 1	12/box	
	25 mm	21-2714-24	3 0610586 01578 0	12/box	
	32 mm	21-2715-24	3 0610586 01579 7	12/box	
● 20 gauge	16 mm	21-2736-24	3 0610586 02322 8	12/box	
	19 mm	21-2734-24	3 0610586 01585 8	12/box	
	25 mm	21-2717-24	3 0610586 01580 3	12/box	
	32 mm	21-2718-24	3 0610586 01581 0	12/box	
● 19 gauge	19 mm	21-2735-24	3 0610586 01586 5	12/box	
	25 mm	21-2720-24	3 0610586 01582 7	12/box	
	32 mm	21-2721-24	3 0610586 01583 4	12/box	

* Manufactured with non-DEHP plasticised tubing



Port-A-Cath Plastic Hub Huber Needles

MR Conditional up to 3 Tesla



Port A Cath Plastic Hub Huber Needles Non coring, non siliconised				
Straight	Length	Reorder No.	GTIN	Unit
● 22 gauge	25 mm	21-2005-24	3 0610586 02399 0	12/box
	38 mm	21-2007-24	3 0610586 02401 0	12/box
● 20 gauge	25 mm	21-2013-24	3 0610586 02403 4	12/box
	38 mm	21-2015-24	3 0610586 02405 8	12/box
● 19 gauge	25 mm	21-2021-24	3 0610586 02407 2	12/box
	51 mm	21-2027-24	3 0610586 02846 9	12/box
90° Bent	Length	Reorder No.	GTIN	Unit
● 22 gauge	13 mm	21-2316-24	3 0610586 02411 9	12/box
	25 mm	21-2006-24	3 0610586 02400 3	12/box
	38 mm	21-2008-24	3 0610586 02402 7	12/box
● 20 gauge	13 mm	21-2317-24	3 0610586 02412 6	12/box
	19 mm	21-2287-24	3 0610586 02408 9	12/box
	25 mm	21-2014-24	3 0610586 02404 1	12/box
	38 mm	21-2016-24	3 0610586 02406 5	12/box
● 19 gauge	25 mm	21-2022-24	3 0610586 02805 6	12/box
	38 mm	21-2026-24	3 0610586 02806 3	12/box

Subcutaneous Infusion Set



Cleo 90 All In One Subcutaneous Infusion set		
Product Code	Cannula Length (mm)	Tubing Length
21-7220-24	6	61 cm (24")
21-7221-24	6	79 cm (31")
21-7222-24	6	107 cm (42")
21-7230-24	9	61 cm (24")
21-7231-24	9	79 cm (31")
21-7232-24	9	107 cm (42")

Supplied sterile in cases of 10 units

The product complies with current legislation and has the corresponding CE marking.
For additional information, warnings and/or safety precautions, refer to the manufacturer's Instructions for Use.

icumedical
human connections