




P.A.S. Port™

Implantable Venous Access System

icumedical
human connections

P.A.S. Port

Implantable Venous Access System



P.A.S. Port T2 Power P.A.C.™—Titanium Portal						
Product Code	Kit/Tray	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (in)

21-4872-24	Kit	6 Fr	Polyurethane	1.9	1.0	30
21-4873-24	Tray					

P.A.S. Port Elite - Plastic/Titanium Portal						
Product Code	Kit/Tray	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (in)

21-4590-24	Kit	6 Fr	Polyurethane	1.9	1.0	30
21-4591-24	Tray					

Kits include portal, catheter, Port-A-Cath™ straight needle, blunt needle, vein pick and Point-Lok™ device.

Trays include portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra-thin wall introducer needle, "I" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

icumedical
human connections