



PharmAssist™ Pro

Peristaltic Pump & Gravimetric Scale

Convenient and efficient fluid compounding technology designed for the fast transfer of fluids and the preparation of ambulatory cassettes, stock bags, and syringes

icumedical
human connections

Fast and Convenient Fluid Compounding Technology

The PharmAssist Pro peristaltic pump is designed to help you accurately and efficiently transfer large volumes of diluents, prepare admixtures, and reconstitute drug vials. Its user-friendly, easy-to-read colour touch screen allows you to enter preparation data more efficiently, while a powerful motor helps streamline pharmacy operations and helps maximise the speed and accuracy of transfers.

Colour Touch Screen

Enter information more efficiently and accurately with an easy-to-read touch screen.

Automated Air Removal

Fills ambulatory cassettes, IV bags, and elastomeric balls while mitigating the presence of air bubbles.

Integrated Scale

Receive a real-time indication of accuracy with the gravimetric scale accessory.*

*sold separately

Manage Compounding Logs

Document all your preparations and device calibrations and easily export log files in the event of an audit.



PharmAssist Pro Consumables

Product	Description
BX01	Single Lead Heavy Duty Tubing
BX01LL	Single Lead Heavy Duty Tubing w/ Double Luer
BX02	Dual Lead Heavy Duty Tubing
BX03	Triple Lead Heavy Duty Tubing
LPA01	Single Lead Tubing for Use w/ Lipids
PA01	Single Lead Tubing
PA01-XS	Single Lead Short Tubing
PA01LL	Single Lead Tubing w/ Double Luer
PA02	Double Lead Tubing
PA03	Triple Lead Tubing
SA250	Syringe Adaptor

PharmAssist Pro Accessories

Product	Description
DS1000-04-01	PharmAssist Pro Pump Module
DS1030-24-01	Gravimetric Scale Kit
DS1020-04-01	PharmAssist Pro Printer Kit
DS1010-04-01	PharmAssist Pro Scanner Kit
SH01	Syringe Holder
DS1040	PharmAssist Pro Foot Pedal

Reference Specifications

Description	Range
Operating Temperature	+10°C to +40°C
Operating Humidity	15% to 60% (noncondensing)
Operating Altitude	0 to 2,000 meters above sea level
Storage Temperature	0°C to +60°C
Storage Humidity	15% to 85% (noncondensing)
Transportation Temperature	0°C to +60°C
Transportation Humidity	15% to 85% (noncondensing)
Physical Dimensions	145 mm x 249 mm x 218 mm
Weight	5,09 kg
Electrical Operating Power Range	100 – 240 V~ 1.2 – 0.6 A 50 – 60 Hz
Electrical Fuse	 250V F 2.5A 5x20 mm
Accuracy Specifications	+/- 1% >10 mL; +/- 2% <10 mL >= 5 mL; +/- 4% <5 mL >= 1 mL; +/- 10% <1 mL >= 0.5 mL;

The PharmAssist Pro system is not intended to screen, monitor, treat, diagnose or prevent any specific condition or disease.

The product complies with current legislation and has the corresponding CE marking. For additional information, warnings and / or safety precautions, refer to the Instructions for Use.

icumedical
human connections