

Get More From Our New Line of Procedure Trays with Hypodermic Needle-Pro Edge™ Safety Device

ICU Medical offers a new line of procedure tray configurations with the Needle-Pro Edge safety device to meet your needs. Based on expert input to help make every procedure safer and easier, every new kit includes added technology and supplies to support safety, sterility, and convenience while delivering patient care.



Convenient for patients. Convenient for clinicians.

- › Help support safety initiatives with needlestick prevention devices provided for every needle
- › Help enhance patient safety during procedures with glass particle filtration
- › Get added convenience with tools that may include luer-lock syringes, surgical skin markers, and an adequate supply of larger gauze



Built-in Needlestick Prevention

Both the anesthetic and procedure needles include safety devices to help protect against needlestick injury.

Our procedure trays are available in a wide range of configurations to help meet clinician and patient needs.

Lumbar Puncture Standard Tray

Item Number	Description
4823-20	Portex™ Lumbar Puncture Tray, Adult (Lancet Point Spinal Needle 18G x 3 ½")
4824-20	Portex Lumbar Puncture Tray, Adult (Lancet Point Spinal Needle 20G x 3 ½")
4825-20	Portex Lumbar Puncture Tray, Adult (Lancet Point Spinal Needle 22G x 3 ½")
4826-20	Portex Lumbar Puncture Tray, Child (Lancet Point Spinal Needle 22G x 2 ½")
4827-20	Portex Lumbar Puncture Tray, Child (Lancet Point Spinal Needle 22G x 1 ½")
Lumbar Puncture Standard Tray Components	
(1) 25G x 5/8" Needle (1) 22G x 1 ½" Needle (1) 18G x 1 ½" Blunt Filter Needle (1) Plastic Syringe (3 mL, Luer Lock) (1) Plastic Syringe (5 mL, Luer Lock) (1) Manometer (1) 4-Way Stopcock	(4) Numbered Specimen Vials with Screw Caps (1) 6" Flexible Extension Unit (1) Lidocaine HCl (1%), 5 mL (2) Towels (5) Gauze Sponges (3" x 3") (3) Sponge Applicators (1) Povidone-Iodine Solution, ¾ oz
	(1) Needle Stick Pad (1) Medication Label (1) Adhesive Bandage (4) Pressure Sensitive Labels (1) Fenestrated Drape

Lumbar Puncture Safety Tray

Item Number	Description
4824PL-20	Portex Lumbar Puncture Tray, Adult (Lancet Point Spinal Needle 20G x 3 ½")
4825PL-20	Portex Lumbar Puncture Tray, Adult (Lancet Point Spinal Needle 22G x 3 ½")
4826PL-20	Portex Lumbar Puncture Tray, Child (Lancet Point Spinal Needle 22G x 2 ½")
4827PL-20	Portex Lumbar Puncture Tray, Child (Lancet Point Spinal Needle 22G x 1 ½")
Lumbar Puncture Safety Tray Components	
(1) 25G x 1" Needle-Pro Edge™ Safety Device (1) 22G x 1 ½" Needle-Pro Edge Safety Device (1) Plastic Syringe (3 mL, Luer Lock) (1) Plastic Syringe (5 mL, Luer Lock) (1) Manometer (1) 4-Way Stopcock (4) Numbered Specimen Vials with Screw Caps	(1) 6" Flexible Extension Unit (1) Filter Straw (1) Lidocaine HCl (1%), 5 mL (1) Point-Lok™ Needle Protection Device (2) Towels (10) Gauze Sponges (3" x 3") (3) Sponge Applicators
	(1) Povidone-Iodine Solution, ¾ oz (1) Needle Stick Pad (1) Medication Label (1) Adhesive Bandage (4) Pressure Sensitive Labels (1) Fenestrated Drape (1) Surgical Skin Marker

Lumbar Puncture Standard Tray (Infant)

Item Number	Description
4963-20	22G x 1 ½" Lancet Point
Lumbar Puncture Standard Tray (Infant) Components	
(1) 25G x 5/8" Needle (1) 22G x 1 ½" Needle (1) 18G x 1½" Blunt Filter Needle (1) Plastic Syringe (3 mL, Luer Lock) (1) Plastic Syringe (5 mL, Luer Lock) (3) Numbered Specimen Vials with Screw Caps	(1) Lidocaine HCl (1%) 2 mL (2) Towels (5) Gauze Sponges (3" x 3") (3) Sponge Applicators (1) Povidone-Iodine Solution, ¾ oz (1) Needle Stick Pad
	(1) Medication Label (1) Adhesive Bandage (3) Pressure Sensitive Labels (1) Fenestrated Drape

Lumbar Puncture Safety Tray (Infant)

Item Number	Description
4963PL-20	Portex Lumbar Puncture Tray, Infant (Lancet Point Spinal Needle 22G x 1 ½")
Lumbar Puncture Safety Tray (Infant) Components	
(1) 25G x 1" Needle-Pro Edge Safety Device (1) 22G x 1 ½" Needle-Pro Edge Safety Device (1) Plastic Syringe (3 mL, Luer Lock) (1) Plastic Syringe (5 mL, Luer Lock) (1) Filter Straw (3) Numbered Specimen Vials with Screw Caps	(1) Lidocaine HCl (1%) 2 mL (1) Point-Lok Needle Protection Device (2) Towels (10) Gauze Sponges (3" x 3") (3) Sponge Applicators (1) Povidone-Iodine Solution, ¾ oz
	(1) Needle Stick Pad (1) Medication Label (1) Adhesive Bandage (3) Pressure Sensitive Labels (1) Fenestrated Drape (1) Surgical Skin Marker

Amniocentesis Standard Tray

Item Number	Description
4885-20	Portex Amniocentesis Tray (Lancet Point Needle 20G x 3 1/2")
4885-22	Portex Amniocentesis Tray (Lancet Point Needle 22G x 3 1/2")
Amniocentesis Standard Tray Components	
(1) 25G x 5/8" Needle (1) 22G x 1 1/2" Needle (1) 18G x 1 1/2" Blunt Filter Needle (1) Plastic Syringe (5 mL, Luer Lock) (1) Plastic Syringe (20 mL, Luer Lock) (3) Specimen Tubes (15 mL) with Caps (1) Lidocaine HCl (1%), 5 mL	(1) Towel (5) Gauze Sponges (3" x 3") (3) Sponge Applicators (1) Povidone-Iodine Solution, 3/4 oz (1) Needle Stick Pad (1) Medication Label (1) Adhesive Bandage

- (1) Fenestrated Drape
- (1) Sheet of 3 Adhesive Labels
- (1) Alcohol Prep Pad
- (3) UV Protector Sleeves
- (1) 6" Ruler

Amniocentesis Safety Tray (Adult)

Item Number	Description
4379-20	Portex Amniocentesis Tray (Lancet Point Needle 20G x 3 1/2")
Amniocentesis Safety Tray Components	
(1) 25G x 1" Needle-Pro™ EDGE™ Safety Device (1) 22G x 1 1/2" Needle-Pro EDGE Safety Device (1) Plastic Syringe (10 mL, Luer Slip) (1) Plastic Syringe (5 mL, Luer Slip) with 22G 1 1/2" Needle (attached) (1) Filter Straw (2) Glass Specimen Vials, Amber with Screw Caps (1) Lidocaine HCl (1%) 5 mL	(1) Point-Lok™ Needle Protection Device (1) Towel (10) Gauze Sponges (3" x 3") (3) Sponge Applicators (1) Povidone-Iodine Solution, 3/4 oz (1) Needle Stick Pad (1) Medication Label

- (1) Adhesive Bandage
- (3) Pressure Sensitive Labels
- (1) Light-Protective Pouch for Syring
- (1) Fenestrated Drape
- (1) 6" Ruler
- (1) Surgical Skin Marker

Myelogram Safety Tray

Item Number	Description
4968-20	Portex Myelogram Tray (Lancet Point Spinal Needle 22G x 3 1/2")
Myelogram Safety Tray Components	
(1) 25G x 1" Needle-Pro Edge Safety Device (1) 22G x 1 1/2" Needle-Pro Edge Safety Device (1) 18G x 1 1/2" Blunt Filter Needle (1) Plastic Syringe (5 mL, Luer Lock) (1) Plastic Syringe (20 mL, Luer Lock) (3) Numbered Specimen Vials with Screw Caps (1) 6" Flexible Extension Unit	(1) 20" Extension Set (1) Filter Straw (1) Lidocaine HCl (1%), 5 mL (2) Point-Lok Needle Protection Device (1) Towel (10) Gauze Sponges (3" x 3") (3) Sponge Applicators

- (1) Povidone-Iodine Solution, 3/4 oz
- (1) Needle Stick Pad
- (1) Medication Label
- (1) Adhesive Bandage
- (3) Pressure Sensitive Labels
- (1) Fenestrated Drape
- (1) Surgical Skin Marker

To learn more about safety procedure trays from
ICU Medical, visit www.icumed.com/pain-management

