

# Gripper Plus™

## Port Access Needle



### Safety Huber Needle

- › A safety port access needle to administer or withdraw fluids from implanted ports
- › Available in different sizes and lengths
- › Safety mechanism to be activated when removing the needle — audible click and visible engagement confirm that the safety arm is activated into locked position
- › Thick cushioned needle platform sits against the skin when in-situ
- › Removable contoured grip for easy insertion
- › Non-siliconized, non-coring needle
- › The 20 cm-long extension set is available with or without Y-site with removable injection cap on luer lock or needleless Y-site luer-activated needleless injection site
- › Extension set with color-coded clamp(s)
- › Priming volume 0.3 mL, 0.5 mL, or 0.8 mL depending on the extension set with or without the Y-site versions
- › Pressure-rated up to 45 psi
- › Positive ANTT evaluation

### Ordering Information

Needleless Y-site — luer-activated needleless injection site				Y-site removable split septum injection cap on luer-lock Y-site				Without Y-site			
Product Code	Gauge	Length	Units	Product Code	Gauge	Length	Units	Product Code	Gauge	Length	Units
21-2960-24		5/8 in		21-2860-24		5/8 in		21-2770-24		5/8 in	
21-2961-24	● 22	3/4 in	12/box	21-2861-24	● 22	3/4 in	12/box	21-2766-24	● 22	3/4 in	12/box
21-2962-24		1 in		21-2862-24		1 in		21-2760-24		1 in	
21-2963-24		1 1/4 in		21-2863-24		1 1/4 in		21-2761-24		1 1/4 in	
21-2964-24		5/8 in		21-2864-24		5/8 in		21-2769-24		5/8 in	
21-2965-24	● 20	3/4 in	12/box	21-2865-24	● 20	3/4 in	12/box	21-2767-24	● 20	3/4 in	12/box
21-2966-24		1 in		21-2866-24		1 in		21-2762-24		1 in	
21-2967-24		1 1/4 in		21-2867-24		1 1/4 in		21-2763-24		1 1/4 in	
21-2968-24		3/4 in		21-2868-24		3/4 in		21-2768-24		3/4 in	
21-2969-24	● 19	1 in	12/box	21-2869-24	● 19	1 in	12/box	21-2764-24	● 19	1 in	12/box
21-2970-24		1 1/4 in		21-2870-24		1 1/4 in		21-2765-24		1 1/4 in	

