




es. 2



SuperCath™ 5


Safety IV Catheter

Put precision into practice with peripheral IV catheters utilising integrated blood-containment technology and fully encapsulated needle protection.

icumedical
human connections

Help Minimise Blood Exposure While Reducing the Risk of Accidental Needlestick Injury

The SuperCath 5 safety IV catheter combines blood-containment technology with fully encapsulated needle protection to help your efforts to reduce the risk of needlestick injury and comply with the EU Sharps Directive.¹



Help minimise the risk of blood exposure

Avoid blood leakage during catheter insertion and setup with an integrated blood-containment valve.

Reduce the risk of needlestick injury

Avoid accidental needlesticks and needle exposure with push-button, spring-loaded needle encapsulation.

Avoid multiple insertion attempts

Help ensure the catheter is properly in the vessel on the first stick with dual flash technology.

Product Information

Product Code Non-Winged	Product Code Winged	Gauge	Length		O.D. (mm)	I.D. (mm)	Flow Rate (mL /min)
			(in)	(mm)			
SP125-26-14T*	SP125-26-14-WT*	●26	0.56	14	0.6	0.4	14
SP125-26-19T*	SP125-26-19-WT*	●26	0.75	19	0.6	0.4	12
SP125-24-14T*	SP125-24-14-WT*	●24	0.56	14	0.7	0.4	21
SP125-24-19T*	SP125-24-19-WT*	●24	0.75	19	0.7	0.4	16
SP120-24-19T	SP120-24-19-WT	●24	0.75	19	0.7	0.4	16
SP120-22-25T	SP120-22-25-WT	●22	1.00	25	0.9	0.6	31
SP120-22-31T	SP120-22-31-WT	●22	1.25	31	0.9	0.6	28
SP120-20-25T	SP120-20-25-WT	●20	1.00	25	1.1	0.7	62
SP120-20-31T	SP120-20-31-WT	●20	1.25	31	1.1	0.7	55
SP120-20-50T	SP120-20-50-WT	●20	2.00	50	1.1	0.7	52
SP120-18-31T	SP120-18-31-WT	●18	1.25	31	1.3	0.9	100
SP120-18-50T	SP120-18-50-WT	●18	2.00	50	1.3	0.9	90
SP120-16-31T	SP120-16-31-WT	●16	1.25	31	1.7	1.3	183
SP125-16-31T*	SP125-16-31-WT*	●16	1.25	31	1.7	1.3	183
SP125-14-31T*	SP125-14-31-WT*	●14	1.25	31	2.1	1.6	272

*No blood-containment valve.



icumedical
human connections

¹ Council Directive 2010/32/EU - 10 May 2010
SuperCath is a trademark of MediKit Co., Ltd.
© 2026 ICU Medical, Inc. | P26-4549-EMEA

The product complies with current legislation and has the corresponding CE marking. For additional information, warnings and / or safety precautions, refer to the manufacturer's user manual.