




P.A.S. Port™

Implantable Venous Access System

icumedical
human connections

P.A.S. Port

Implantable Venous Access System



A portal with a small footprint
Placed in the arm or chest



P.A.S. Port™ T2 Power P.A.C.



P.A.S. Port™ T2

The titanium portal is available also in a power version (300 psi): "CT" embedded all through the silicone septum.



Gripper Plus™ Safety Huber Needle



Gripper Plus™ Power P.A.C. Safety Huber Needle

Please contact your local sales representative for further information on our wide range of port access needles.

P.A.S. Port T2 Power P.A.C.—Titanium Portal

Product Code	Kit/Tray	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (in)
21-4872-24	Kit	6 Fr	Polyurethane	1.9	1.0	30
21-4873-24	Tray					

P.A.S. Port T2/Titanium Portal

Product Code	Kit/Tray	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (in)
21-4572-24	Kit	6 Fr	Polyurethane	1.9	1.0	30
21-4573-24	Tray					

P.A.S. Port Elite - Plastic/Titanium Portal

Product Code	Kit/Tray	Introducer Size	Catheter Material	OD (mm)	ID (mm)	Length (in)
21-4590-24	Kit	6 Fr	Polyurethane	1.9	1.0	30
21-4591-24	Tray					

Kits include portal, catheter, Port-A-Cath™ straight needle, blunt needle, vein pick and Point-Lok™ device.

Trays include portal, catheter, Port-A-Cath straight needle, blunt needle (except preassembled systems), vein pick, 18G extra-thin wall introducer needle, "J" guidewire(s), vessel dilator/sheath introducer assembly, 90° Port-A-Cath needle, Gripper Plus safety needle(s), 2 syringes, tunneling tool and Point-Lok device(s).

icumedical
human connections