

























CADD™ Ambulatory Tubing Sets with Flow Stop Free Flow Protection

CADD Medication Cassette Reservoirs for CADD-Solis™, CADD-Solis VIP™ and CADD-Legacy™ Pumps - Must be used with a CADD extension set							
	Reference No.	Units	Description	Tubing Length	Medication Bag Plasticiser	Replaces	
	21-7301	12/box	50 mL CADD™ medication cassette reservoir with flow stop, clamp, female luer and non-vented cap	8 in (20 cm)	DEHP	21-7001	
	21-7302	12/box	100 mL CADD™ medication cassette reservoir with flow stop, clamp, female luer and non-vented cap	8 in (20 cm)	DEHP	21-7002	
	21-7300	12/box	100 mL yellow CADD™ medication cassette reservoir with flow stop, clamp, female luer and non-vented cap	8 in (20 cm)	DEHP		
	21-7308	12/box	250 mL clear CADD™ medication cassette reservoir with flow stop, clamp, female luer and non-vented cap	8 in (20 cm)	Non-DEHP		
	21-7309	12/box	250 mL yellow CADD™ medication cassette reservoir with flow stop, clamp, female luer and non-vented cap	8 in (20 cm)	Non-DEHP		
	21-7310	12/box	250 mL blue CADD™ medication cassette reservoir with flow stop, clamp, female luer and non-vented cap	8 in (20 cm)	Non-DEHP		
	21-7600	12/box	100 mL yellow CADD™ medication cassette reservoir with NRFit™ connector, flow stop, female luer, clamp, and non-vented stopper	8 in (20 cm)	DEHP	21-7300	
	21-7609	12/box	250 mL yellow CADD™ medication cassette reservoir with NRFit™ connector, flow stop, female luer, clamp, and non-vented stopper	8 in (20 cm)	Non-DEHP	21-7309	
	21-7664-24	50/box	CADD Filling Adapter with male NRFit™ connector and female luer				
	24-1001-64	10/box	Infusion Adapter to convert female luer to NRFit™ female				
	24-1009-64	50/box	Yellow cap, female				
	24-1004-64	50/box	Yellow cap, male, non-vented				

CADD Administration Sets for CADD-Solis, CADD-Solis VIP and CADD-Legacy Pumps							
	Reference No.	Units	Description	Tubing Length	Priming Vol.	Fluid Path Plasticiser	Replaces
	21-7321	12/box	CADD™ administration set with female luer, flow stop, clamp, and one-way checkvalve with male luer	102 in (259 cm)	2.7 mL	Non-DEHP	21-7021
	21-7322	12/box	CADD™ administration set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	78 in (198 cm)	2.3 mL	Non-DEHP	21-7022
	21-7323	12/box	CADD™ blue-striped administration set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	78 in (198 cm)	2.3 mL	Non-DEHP	21-7023
	21-7324	12/box	CADD™ yellow-striped administration set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	123 in (312 cm)	3.2 mL	Non-DEHP	21-7024
	21-7359	12/box	CADD™ administration set with male luer, flow stop, clamp, and one-way checkvalve with male luer	69 in (175 cm)	2.0 mL	Non-DEHP	21-7059
	21-7390	12/box	CADD™ administration set with female luer, flow stop, 7.6 cm y-extension, one-way checkvalve with female luer, one-way checkvalve with male luer, and clamps	102 in with 3 in y-extension (259 cm with 7.6 cm y-extension)	4.0 mL	Non-DEHP	21-7090
	21-7391	12/box	CADD™ administration set with bag spike, flow stop, 7.6 cm y-extension, one-way check valve with female luer, one-way check valve with male luer, and clamps	108 in with 3 in y-extension (274 cm with 7.6 cm y-extension)	4.2 mL	Non-DEHP	21-7091
	21-7624	12/box	CADD™ yellow-striped administration set with NRFit™ connector, bag spike, flow stop, clamp, and one-way checkvalve with male luer	123 in (312 cm)	3.2 mL	Non-DEHP	21-7024 21-7324

CADD Administration Sets with Air-eliminating Filters for CADD-Solis, CADD-Solis VIP and CADD-Legacy Pumps							
	Reference No.	Units	Description	Tubing Length	Priming Vol.	Fluid Path Plasticiser	Replaces
	21-7343	15/box	CADD® administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, needleless injection site, and one-way check valve with male luer	114 in (290 cm)	14 mL	Non-DEHP	
	21-7346	15/box	CADD® administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male luer	94 in (239 cm)	7 mL	Non-DEHP	
	21-7363	15/box	CADD® administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, spiral tubing, clamp, and one-way checkvalve with male luer	94 in (239 cm)	7 mL	Non-DEHP	
	21-7349	12/box	CADD® yellow-striped administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	130 in (330 cm)	6 mL	Non-DEHP	
	21-7394	12/box	CADD® administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	108 in (274 cm)	3.3 mL	Non-DEHP	21-7094
	21-7395	12/box	CADD® administration set with female Luer, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	102 in (259 cm)	3.1 mL	Non-DEHP	21-7095
	21-7649	12/box	CADD® yellow-striped administration set with NRFit™ connector, bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	130 in (330 cm)	6 mL	Non-DEHP	21-7349

CADD High Volume Administration Sets for use with CADD-Solis, CADD-Solis VIP and CADD-Legacy Pumps							
	Reference No.	Units	Description	Tubing Length	Priming Vol.	Fluid Path Plasticiser	Replaces
	21-7345	15/box	CADD® high volume administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, injection site and one-way checkvalve with male luer	118 in (300 cm)	17 mL	Non-DEHP	
	21-7357	15/box	CADD® high volume administration set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	92 in (234 cm)	11 mL	Non-DEHP	21-7057
	21-7364	15/box	CADD® high volume administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, luer-activated needleless injection site, clamp, and one-way checkvalve with male luer	117 in (297 cm)	17 mL	Non-DEHP	
	21-7361	15/box	CADD® high volume administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, injection site, and one-way checkvalve with male luer	117 in (297 cm)	17 mL	Non-DEHP	21-7381
	21-7386	15/box	CADD® high volume administration set with bag spike, flow stop, 1.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	104 in (264 cm)	15 mL	Non-DEHP	
	21-7384	15/box	CADD® yellow-striped high volume administration set with bag spike, flow stop, clamp, and one-way checkvalve with male luer	123 in (312 cm)	16 mL	Non-DEHP	
	21-7365	15/box	CADD® yellow-striped high volume administration set with bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	135 in (343 cm)	20 mL	Non-DEHP	
	21-7684	15/box	CADD® yellow-striped high volume administration set with NRFit™ connector, bag spike, flow stop, clamp, and one-way checkvalve with male luer	123 in (312 cm)	16 mL	Non-DEHP	
	21-7665	15/box	CADD® yellow-striped high volume administration set with NRFit™ connector, bag spike, flow stop, 0.2-micron air-eliminating filter, clamp, and one-way checkvalve with male luer	134 in (340 cm)	20 mL	Non-DEHP	

CADD Extension Sets for CADD-Solis, CADD-Solis VIP and CADD-Legacy Pumps - Designed to be used with a CADD medication cassette reservoir							
	Reference No.	Units	Description	Tubing Length	Priming Vol.	Fluid Path Plasticiser	Replaces
	21-7045-24	50/box	CADD® extension set with male/male luer, integral anti-siphon valve, and clamp	30 in (76 cm)	1.1 mL	Non-DEHP	
	21-7060-24	12/box	CADD® extension set with male/male luer, integral anti-siphon valve, and clamp	45 in (114 cm)	1.4 mL	Non-DEHP	
	21-7046-24	50/box	CADD® extension set with male/male luer, integral anti-siphon valve, and clamp	45 in (114 cm)	1.4 mL	Non-DEHP	
	21-7061-24	12/box	CADD® extension set with male/male luer, integral anti-siphon valve, and clamp	45 in (114 cm)	1.4 mL	Non-DEHP	
	21-7047-24	50/box	CADD® extension set with male/male luer, integral anti-siphon valve, and clamp	60 in (152 cm)	1.7 mL	Non-DEHP	
	21-7062-24	12/box	CADD® extension set with male/male luer, integral anti-siphon valve, and clamp	60 in (152 cm)	1.7 mL	Non-DEHP	
	21-7092-24	12/box	CADD® checkvalve extension set with male/male luer, integral anti-siphon valve, y-extension with one-way checkvalve with female luer, and clamp	68 in with 3 in y-extension (173 cm with 7.6 cm y-extension)	2.9 mL	Non-DEHP	
	21-7115-24	12/box	CADD® checkvalve extension set with male/male luer, 0.2-micron air-eliminating filter, integral anti-siphon valve with clamp, and y-extension with one-way checkvalve with female luer	69 in with 3 in y-extension (175 cm with 7.6 cm y-extension)	3.2 mL	Non-DEHP	
	21-7105-24	12/box	CADD® yellow-striped extension set with male/male luer, integral anti-siphon valve, and clamp	90 in (229 cm)	2.3 mL	Non-DEHP	
	21-7044-24	25/box	CADD® extension set with male/male luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	55 in (140 cm)	4 mL	Non-DEHP	
	21-7052-24	50/box	CADD® microbore extension set with male luer, 0.2-micron air-eliminating filter, clamp, and integral anti-siphon valve with male luer	60 in (152 cm)	1.1 mL	Non-DEHP	
	21-7106-24	50/box	CADD® extension set with male/male luer, 0.2-micron air-eliminating filter, integral anti-siphon valve, and clamp	60 in (152 cm)	1.8 mL	Non-DEHP	
	21-7147-24	12/box	CADD® yellow-striped extension set with male/male luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	100 in (254 cm)	5 mL	Non-DEHP	
	21-7108-24	50/box	CADD® yellow-striped extension set with male/male luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	90 in (229 cm)	2.4 mL	Non-DEHP	
	21-7109-24	50/box	CADD® blue-striped extension set with male/male luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	90 in (229 cm)	2.4 mL	Non-DEHP	
	21-1200-24	12/box	Bag spike set with male luer lock and clamp	12 in (30 cm)	0.5 mL	Non-DEHP	
	21-7605-24	12/box	CADD® yellow-striped extension set with NRFit™ connector, male/male luer, integral anti-siphon valve, and clamp	90 in (229 cm)	2.3 mL	Non-DEHP	21-7105
	21-7608-24	50/box	CADD® yellow-striped extension set with NRFit™ connector, male/male luer, 0.2-micron air-eliminating filter, and integral anti-siphon valve with clamp	90 in (229 cm)	2.4 mL	Non-DEHP	21-7108