

Plum 360™ Infusion System Catalogue

icumedical
human connections

Connect Smarter. Care Safer™

Get award-winning clinical performance.

The Best in KLAS Plum 360 infusion system gives you automatic safety with unmatched efficiency and interoperability in one connected IV system.

The Plum 360 features a unique cassette and delivery system which enables:

- A direct connection of the secondary line to the Plum cassette, so the pump knows which container — primary or secondary — it is delivering from
- Concurrent delivery of two compatible medications at independent rates through a single IV line
- Capture of 1 mL of air before alarming, helping to reduce air-in-line alarms
- Integration with your EHR for enhanced safety and efficiency






Item Number
30010
Description
Plum 360 Infuser compatible with ICU Medical MedNet™ software

Dedicated Plum™ Sets


The Plum cassette is at the heart of the unique Plum delivery system. The direct connection of the secondary line to the cassette provides control and management for both primary and secondary infusions. You can feel confident that both your piggyback and concurrent infusions (of compatible medications) are delivering as programmed.

The Plum cassette includes an air trap that captures up to 1 mL of air before alarming. This has been shown to reduce air-in-line alarms. If an air-in-line alarm does occur, air can be easily removed without disconnecting from the patient. This closed system also helps to reduce the risk of contamination and limit exposure to hazardous medications.




Plum 360 Primary Sets


Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14001	Primary Plum Set, 15 Micron Filter in Sight Chamber, Clave Secondary Port, Clave Y-Site, Secure Lock, 272 cm/19 mL	272	19.0	PVC non-DEHP	50




Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14015	Primary Plum Set, 15 Micron Filter in Sight Chamber, Clave Secondary Port, 0.2 Micron Filter, Clave Y-Site, Polyethylene Lined Tubing, Secure Lock, 272 cm/23 mL	272	23.0	PVC non-DEHP/PE	50




Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14006	Primary Plum Set, 15 Micron Filter in Sight Chamber, Polyethylene Lined Light Resistant Tubing, Distal Microbore Tubing, Secure Lock, 272 cm/12 mL	272	12.0	PVC non-DEHP/PE	50




Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14007	Primary Plum Set, 15 Micron Filter in Sight Chamber, Clave Secondary Port, Clave Y-Site, Polyethylene Lined Light Resistant Tubing, Distal Microbore Tubing, Secure Lock, 272 cm/12 mL	272	12.0	PVC non-DEHP/PE	50




Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14012	Primary Plum Set, 15 Micron Filter in Sight Chamber, Clave Secondary Port, Secure Lock, 225 cm/17 mL	225	17.0	PVC non-DEHP	50




Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14017	Primary Plum Microdrip Set, 15 Micron Filter in Sight Chamber, Clave Secondary Port, Clave Y-Site, Polyethylene Lined Tubing, Secure Lock, 272 cm/20 mL	272	20.0	PVC non-DEHP/PE	50




Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14018	Single Channel, Primary Plum Set, 15 Micron Filter in Sight Chamber, Secure Lock, 272 cm	272	19 mL	PVC non-DEHP	50



Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14024	Primary Plum Microdrip™ Set, 15µm Filter in Sight Chamber, Clave Secondary Port, Clave Y-Site, Distal Microbore Tubing, Secure-Lok, 262 cm/12.7 mL	262	12.7 mL	PVC non-DEHP	50



Item No	Description	Length (cm)	PV (mL)	Tubing Material	Case Qty
14211	Primary Plum™ Blood Set, 200µm Filter, PP Secondary Port, Non-Vented, Secure-Lok, 279 cm/52 mL	279	52.0 mL	PVC non-DEHP	50



Pump Accessories

ICU Medical provides a range of IV poles and mounting accessories for our Infusion pumps to help you to efficiently organise your space around the patient.



Item No	Case Qty
391192502	1
Description	
Provita heavy duty pole 323840. Pole can accommodate up to 5 x Plum 360 pumps. A Triple pump holder and a dual pump holder are required for this configuration	



Item No	Case Qty
391052502	1
Description	
Provita standard pole 175191. Pole can accommodate up to 2 x Plum 360 pumps (attached directly to the pole)	

Provita IV Pole Mount Accessories



Item No	Case Qty
391233602	1
Description	
Provita dual pump holder. Use to attach two Plum 360 pumps to a premium or heavy duty IV pole	



Item No	Case Qty
391202502	1
Description	
Provita triple pump holder. Use to attach three Plum 360 pumps to a premium or heavy duty IV pole	



Item No	Case Qty
391213602	1
Description	
Provita bag holder. Keep bags handy and lines organised when multiple pumps are mounted to a pole	



Item No	Case Qty
391223602	25
Description	
Cable Tidy. Keep cables tidy when multiple pumps are mounted to a pole. Can only be used with Heavy duty pole (391192502)	



Item No	Case Qty
391293602	25
Description	
Provita single cable tidy for 25 mm tube Keep your cable tidy on the pole (single cable)	

ALWAYS READ THE LABEL AND FOLLOW THE DIRECTIONS FOR USE

icumedical
human connections

© 2026 ICU Medical, Inc. | P26-5581-ANZ

www.icumed.com